

Eichrom Method and Application Note Updates

Daniel McAlister, Senior Chemist

Radiobioassay & Radiochemical Measurements Conference



Knoxville, Tennessee

eichrom®

 **RRMC**
Radiobioassay & Radiochemical Measurements Conference

RRMC 2014
Knoxville, Tennessee USA
October 27 - October 31, 2014

Outline

- Updates to Eichrom “classic” methods
- Application Notes for new methods
- Highlights of important steps

Eichrom Methods

- No major changes to the chemistry
- Format/Organization
- Separated CeF_3 /Electrodeposition into Separate Methods
- Improve Reference Material
- Add One Page Flowcharts
- Column and Cartridge Based Methods

Method	Analyte	Matrix	Column	Cartridge
ACS07	U	soil	x	x
ACU02	Am-Cm/Pu/U	urine	x	x
ACW01	U/Th	water	x	x
ACW02	U	water	x	x
ACW03	Am-Cm/Pu/U	water	x	x
ACW04	Am	water	x	x
ACW08	Th/Np	water	x	x
ACW10	Th	water	x	x
ACW11	Total Alpha	water	N/A	
ACW13	Th/Pu/U	water	x	x
ACW-16	Am-Cm/Pu/U/Np/Th	water	x	x
ACW-17	Am-Cm/Pu/U/Np/Th/Sr	water	x	x
FEW01	Fe	water	x	x
NIW01	Ni	water	x	x
OTS01	Pb	soil	x	x
OTW01	Pb/Po	water	x	x
OTW02	Tritium	water	x	x
RAW01	Ra	water	x	x
RAW03	Ra	water	x	x
RAW04	Ra	water	x	x
SPA01	CeF3 ppt		N/A	
SPA02	Electrodeposition		N/A	
SPA03	Am-Rare Earth Separation		x	x
SRW01	Sr	water	x	x
TCS01	Tc	soil	x	x
TSU01	Tc	urine	Disk	
TCW01	Tc	water	x	x
TCW02	Tc	water	Disk	
TP01	Set-up of U-232 Self-Cleaning Tracer			
VBS01	Vacuum Box System Guide			

Americium, Plutonium and Uranium in Urine (with Vacuum Box System)

1. Scope

- 1.1. This procedure describes a method for separation of americium, plutonium and uranium in urine. After completing this method, source preparation for measurement of americium, plutonium and uranium by alpha spectrometry is performed by electrolytic deposition onto stainless steel planchets (Eichrom Method SPA02) or by rare earth fluoride microprecipitation onto polypropylene filters (Eichrom Method SPA01).
- 1.2. This method does not address all aspects of safety, quality control, calibration or instrument set-up. However, enough detail is given for a trained radiochemist to achieve accurate and precise results for the analysis of the analyte(s) from the appropriate matrix, when incorporating the appropriate agency or laboratory safety, quality and laboratory control standards.

2. Summary of Method

- 2.1. Americium, plutonium and uranium are separated by Eichrom UTEVA and TRU resins prior to measurement by alpha spectrometry. A calcium phosphate precipitation technique can be used to concentrate and remove actinides from urine samples. Tracers are used to monitor chemical recoveries and correct results to improve precision and accuracy.

3. Significance of Use

- 3.1. This is a rapid, reliable method for measurement of actinides in urine samples that is more cost-effective and efficient than traditional ion exchange, solvent extraction and precipitation techniques. It also takes advantage of vacuum assisted flow and the ability to stack resin cartridges for simultaneous extraction of U, Pu and Am on UTEVA and TRU Resins. This approach can reduce resin separation time by more than half when compared to methods using gravity flow resin columns.

4. Interferences

- 4.1. Actinides with unresolvable alpha energies such as Am-241 and Pu-238 or Np-237 and U-234 must be chemically separated to enable measurement. This method separates these isotopes effectively.

- 4.2. Very high levels of phosphate may reduce the recovery of actinides during the calcium phosphate precipitation and column separation. Adjusting the amount of phosphate added to coprecipitate the actinides may be necessary in these cases.

5. Apparatus

- 5.1. Analytical balance - 0.0001 g sensitivity
- 5.2. Cartridge reservoirs - 10mL Eichrom Part AC-200-RV10 or 20mL Eichrom Part AC-200-RV20
- 5.3. Centrifuge, with rotor and carriers for 50mL and 250mL tubes
- 5.4. Centrifuge tubes, 50mL and 250mL
- 5.5. Eichrom vacuum box system - Eichrom Part No. AC-12-BOX or AC-24-BOX
- 5.6. Vacuum Box White Inner Support Tube - PE - Eichrom Part No. AC-1000-TUBE-PE
- 5.7. Vacuum Box Yellow Outer Tips - Eichrom Part No. AC-1000-OT
- 5.8. Fume hood
- 5.9. Hotplate
- 5.10. Vacuum Pump - 115 V, 60 Hz Fisher part no. 01-092-25 (or equivalent) or house vacuum
- 5.11. Glass stir rods
- 5.12. Vortex mixer

6. Reagents

Note: Analytical grade or ACS grade reagents and trace metal grade (or equivalent) acids are recommended. Evaluation of key reagents, such as aluminum nitrate and ammonium hydrogen phosphate, for contribution to method background levels from naturally occurring radioactive materials is recommended.

Deionized water, All reagents are prepared with deionized water
Appropriate tracers or standards (U-232, Am-243, Pu-242 or Pu-236)
Ammonium hydroxide (57%), concentrated NH_4OH
Barium Chloride, $\text{BaCl}_2 \cdot 2\text{H}_2\text{O}$
Hydrochloric acid (37%), concentrated HCl
Hydrofluoric acid (49%), concentrated HF
Nitric acid (70%), concentrated HNO_3
Aluminum Nitrate Nonahydrate, $\text{Al}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$
Oxalic Acid Dihydrate, $\text{H}_2\text{C}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$
Ammonium Oxalate Monohydrate, $(\text{NH}_4)_2\text{C}_2\text{O}_4 \cdot \text{H}_2\text{O}$
Ammonium Hydrogen Phosphate, $(\text{NH}_4)_2\text{HPO}_4$
Ammonium Thiocyanate, NH_4SCN
Ascorbic Acid Powder, $\text{C}_6\text{H}_8\text{O}_6$

Calcium Nitrate, CaNO_3
Sodium Nitrite, NaNO_2
Femic Nitrate Nonahydrate, $\text{Fe}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$
Sulfamic Acid, H_2NSO_2
Hydrogen Peroxide (30%), concentrated H_2O_2
1-Octanol, $\text{C}_8\text{H}_{17}\text{OH}$
TRU Resin - 2mL prepacked cartridge, 50-100 μm , Eichrom Part TR-R50-S
UTEVA Resin - 2mL prepacked cartridge, 50-100 μm , Eichrom Part UT-R50-S

- 6.1. Ammonium bioxalate (0.1M) - Dissolve 6.3g of oxalic acid and 7.1g of ammonium oxalate in 900mL of water. Dilute to 1L with water.
- 6.2. Ammonium hydrogen phosphate (3.2M) - Dissolve 208g of $(\text{NH}_4)_2\text{HPO}_4$ in 400mL of water. Heat gently to dissolve. Dilute to 500mL with water.
- 6.3. Ammonium thiocyanate indicator (1M) - Dissolve 7.6g of ammonium thiocyanate in 90mL of water. Dilute to 100mL with water.
- 6.4. Ba carrier (5 mg/mL) - Dissolve 4.37g of $\text{BaCl}_2 \cdot 2\text{H}_2\text{O}$ in 250mL of water. Dilute to 500mL with water.
- 6.5. Calcium nitrate (1.25M) - Dissolve 102g of $\text{Ca}(\text{NO}_3)_2$ in 400mL of water. Dilute to 500mL with water.
- 6.6. Hydrochloric acid (1M) - Add 83mL of concentrated HCl to 800mL of water. Dilute to 1L with water.
- 6.7. Hydrochloric acid (4M) - Add 333mL of concentrated HCl to 500mL of water. Dilute to 1L with water.
- 6.8. Hydrochloric acid (4M) - hydrofluoric acid (0.1M) - Add 333mL of concentrated HCl and 3.6mL of concentrated HF to 500mL of water. Dilute to 1L with water.
- 6.9. Hydrochloric acid (5M) - oxalic acid (0.05M) solution - Dissolve 6.3g oxalic acid dihydrate in 400mL of water. Add 417mL concentrated HCl. Dilute to 1L with water.
- 6.10. Hydrochloric acid (9M) - Add 750mL of concentrated HCl to 100mL of water. Dilute to 1L with water.
- 6.11. Nitric acid (2M) - sodium nitrite (0.1M) solution - Add 31mL of concentrated HNO_3 to 200mL of water, dissolve 1.72g of sodium nitrite in the solution. Dilute to 250mL with of water. **Prepare fresh daily.**
- 6.12. Nitric acid solution (3M) - Add 188mL of concentrated HNO_3 to 700mL of water. Dilute to 1L with water.
- 6.13. Nitric acid (8M) - Add 500mL of concentrated HNO_3 to 45 mL of water. Dilute to 1L with water.
- 6.14. Nitric acid (3M) - Aluminum nitrate (1.0M) solution - Dissolve 375g of $\text{Al}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$ in 500mL of water. Add 188mL of concentrated HNO_3 . Dilute to 1L with of deionized water.
- 6.15. Femic Nitrate Solution (5 mg/mL Fe) in 0.1M HNO_3 - To a 500mL volumetric flask, add 18g $\text{Fe}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$, 400mL water and 3.1mL concentrated HNO_3 . Swirl to dissolve. Dilute to 500mL with water.

REFERENCES


- (1) Horwitz, E.P., et al. "Separation and Preconcentration of Actinides from Acidic Media by Extraction Chromatography," *Analytica Chimica Acta*, 281, 361-372, (1993).
- (2) Horwitz, E.P., et al. "Separation and Preconcentration of Uranium from Acidic Media by Extraction Chromatography," *Analytica Chimica Acta*, 266, 25-37, (1992).
- (3) Maxwell, S.L., et al. "Rapid Analysis of Urine and Water Samples," *Journal of Radioanalytical and Nuclear Chemistry*, 275(3), 447-502 (2008).
- (4) ASTM Method D3648-14, "Standard Practices for the Measurement of Radioactivity."
- (5) ASTM Method D7282-06, "Standard Practice for Set-up, Calibration, and Quality Control of Instruments Used for Radioactivity Measurements."

ASTM International, Standards for Radiochemistry, <http://www.astm.org/>


National Analytical Management Program (NAMP), Webinar Series on Actinide chemistry, http://www.wipp.energy.gov/namp/en_content-30-trainingedu.html

National Academy of Sciences, Monograph Series on Radiochemistry, <http://www.lanl.gov/library/find/ebooks/radiochem.php>


- 1) Aliquot up to 1200mL of urine into glass beaker.
- 2) Add yield tracers.
- 3) Add 2-3 drops 1-octanol, 1mL of 1.25M $\text{Ca}(\text{NO}_3)_2$ and 1mL Ba carrier.
- 4) Add 25mL of conc. HNO_3 .
- 5) Cover with watch glass and heat at medium setting for at least 30-60 min.
- 6) Turn down heat. Add 2.5mL of 3.2M $(\text{NH}_4)_2\text{HPO}_4$.
- 7) While stirring sample, slowly add conc. NH_4OH until reaching pH 9.
- 8) After 30 minutes, turn off heat. Allow sample to cool and precipitate to settle or centrifuge.
- 9) Decant supernate and discard as waste.
- 10) Transfer precipitate to centrifuge tube with DI water.
- 11) Centrifuge ~10minutes at 2000rpm. Decant supernate.
- 12) Add 10mL DI water to precipitate. Mix well. Centrifuge. Decant supernate.
- 13) Rinse precipitation beaker with 5mL conc. HNO_3 . Transfer to precipitate in centrifuge tube. Dissolve precipitate and transfer to 100mL beaker.
- 14) Rinse centrifuge tube with 2x 5mL conc. HNO_3 . Transfer to same 100mL beaker.
- 15) Digest precipitate with $\text{HNO}_3/\text{H}_2\text{O}_2$ until white residue obtained.
- 16) Dissolve residue in 20mL 3M HNO_3 -1M $\text{Al}(\text{NO}_3)_3$.
- 17) Add 1mL 1.5M Sulfamic Acid, 0.5 mL Fe, and 1mL 1M Ascorbic Acid.




- 18) Precondition UTEVA-TRU with 5mL 3M HNO_3 .
- 19) Load sample onto UTEVA-TRU. Allow liquid to drain. UTEVA retains U. TRU retains Am and Pu.
- 20) Rinse sample tube with 5mL 3M HNO_3 . Add rinse to UTEVA-TRU. Allow liquid to drain.
- 21) Rinse UTEVA-TRU with 5mL 3M HNO_3 . Allow liquid to drain.
- 22) Separate UTEVA and TRU cartridges.



- 23) Rinse UTEVA cartridge with 15mL 8M HNO_3 .
- 24) Rinse UTEVA cartridge with 5mL 9M HCl .
- 25) Rinse UTEVA cartridge with 20mL 5M HCl -0.05M oxalic acid. Discard to Waste.
- 26) Place clean centrifuge tube below each cartridge. Strip Uranium with 15mL 1M HCl .



- 27) Rinse TRU with 5mL 2M HNO_3 -0.1M NaNO_2 .
- 28) Place clean centrifuge tube below each cartridge. Strip Am with 15mL 4M HCl .
- 29) Place a clean centrifuge tube below each cartridge. Rinse TRU with 25mL 4M HCl -0.1M HF .
- 30) Place a clean centrifuge tube below each cartridge. Strip Pu with 10mL 0.1M ammonium bioxalate.



Update Progress


25%

50%

75%

100%

Draft Update



Formatting and Editing



Final Proofread



Testing



Application Notes

- Significant changes to the chemistry
- Application to new matrix or sample size
- New analyte or combination of analytes
- Make vast library of applications of Eichrom Products more accessible.
- Two pages highlight:
 - Summary of method
 - Reagent/Equipment lists
 - Flowsheet(s)
 - Testing data
 - Reference to literature citation

Rapid Determination of Sr in 1-2 Liter Seawater Samples

AN-1404

Summary of Method Strontium is separated and concentrated from 1-2L samples of seawater with calcium phosphate precipitation, enhanced with 200mg of iron and 1mg of yttrium. Strontium is separated from matrix impurities and potentially interfering radionuclides in the sample using two stacked 2mL cartridges of Eichrom Sr Resin. Radiostromium is measured on a low background gas flow proportional counter. Chemical yield of strontium is determined by gravimetric recovery of native stable strontium in the seawater. Average chemical recovery of strontium is $89 \pm 5\%$ for 1L samples and $82 \pm 4\%$ for 2L samples. Measured values of ^{80}Sr agreed to within 1% and 4% of reference values, for 1L and 2L, respectively, with two hour count times. The minimum detectable activity for ^{80}Sr for 2L samples with a two hour count time is 9.1Bq/L. A single operator can prepare batches of 12-24 samples for measurement of radiostromium in less than 8 hours.

Reagents

Sr Resin, 2mL Cartridges (Eichrom SR-R50-S)
 Nitric Acid (70%)
 Ammonium Hydroxide (57%)
 Deionized Water
 Iron Carrier (50mg/mL Fe, as ferric nitrate)
 Yttrium Carrier (1mg/mL)
 3.2M $(\text{NH}_4)_2\text{HPO}_4$
 Aluminum Nitrate, Nonahydrate
 ^{80}Sr standard
 Oxalic acid
 Boric acid

Equipment

Vacuum Box (Eichrom AC-24-BOX or AC-12-BOX)
 Inner Support Tubes-PE (Eichrom AC-1000-TUBE-PE)
 Yellow Outer Tips (Eichrom AC-1000-OT)
 50mL Centrifuge Tubes
 250-500mL Centrifuge Tubes
 Centrifuge
 Cupped Stainless Steel Planchets (~5mL volume)
 Gas Flow Proportional Counter
 HotPlate
 Analytical Balance
 600mL Glass Beakers
 Vacuum Pump

Figure 1. Sample Preparation





Figure 2. Strontium Resin Separation (Optional ^{80}Y Ingrowth)


- (1) Precondition Sr Resin with 10mL 8M HNO_3
- (2) Load ^{80}Sr sample at 1-2mL/min.
- (3) Rinse sample tube with 5mL 8M HNO_3
- (4) Add tube rinse to Sr Resin. Elute at 1-2mL/min.
- (5) Rinse Sr Resin sequentially with:
 - 15mL 8M HNO_3
 - 10mL 3M HNO_3 - 0.05 oxalic acid
 - 10mL 8M HNO_3
- (6) Dispose of (1) to (5) as waste.
- (7) Strip Sr with 20mL 0.05M HNO_3 at 1mL/min.

- (8) Evaporate samples to dryness on tared cupped stainless steel planchets.
- (9) Rinse Sr sample vials with 2mL 0.05M HNO_3 . Transfer vial rinse to planchets. Evaporate to dryness.



(10) Weigh planchets on an analytical balance to determine gravimetric yield of stable $\text{Sr}(\text{NO}_3)_2$

(11) Measure radiostromium in samples on low background gas flow proportional counter.



(Options for ^{80}Sr Desorbtion)
 When necessary to obtain ^{80}Sr and ^{80}Y data:

(a) Sr fraction from step (7) can be transferred to a liquid scintillation vial. ^{80}Sr can be measured by Cerenkov counting (without LSC cocktail).

(b) Sr fraction from step (10) can be dissolved in 10mL 8M HNO_3 after >7 days of ^{80}Y ingrowth. ^{80}Sr can be removed on Sr Resin. ^{80}Y will elute in Sr Resin load and can be counted by liquid scintillation or gas flow proportional counting.

Performance of ^{80}Sr Method for 1L and 2L Seawater Samples

Sample Replicates	Sample Volume, L	^{80}Sr , Reference Value (mBq/L)	^{80}Sr , Measured Value (mBq/L)	% Bias	Sr carrier % Recovery
11	1	148	150 ± 11	1.2	89 ± 5
4	2	148	154 ± 5	4.2	82 ± 4


2 hour count times
 MDA = 9.1 mBq/L for 2L sample

References


1) Sherrod L. Maxwell, Brian K. Culligan, Robin C. Utsey, "Rapid determination of radiostromium in seawater samples," *J. Radioanal. Nucl. Chem.*, 298(2), 867-875 (2013).

<u>Number</u>	<u>Title</u>	<u>Analytes</u>	<u>Matrix</u>
AN-1401	Rapid Determination of ²²⁶ Ra in Emergency Urine and Water	²²⁶ Ra	Emergency Urine/Water
AN-1402	Rapid Determination of Sr in Emergency Milk Samples	Sr	Emergency Milk
AN-1403	Rapid Determination of Sr in 50g Soil Samples	Sr	Soil (50g)
AN-1404	Rapid Determination of Sr in 1-2 Liter Seawater Samples	Sr	Seawater (1-2L)
AN-1405	Rapid Determination of Sr in Vegetation Samples	Sr	Vegetation
AN-1406	Rapid Determination of Actinides in Vegetation Samples	U, Pu, Am, Cm	Vegetation
AN-1407	Rapid Determination of Sr in Animal Tissue Samples	Sr	Animal Tissue
AN-1408	Rapid Determination of Actinides in Animal Tissue Samples	U, Pu, Am, Cm	Animal Tissue
AN-1409	Rapid Determination of Sr in Building Materials	Sr	Building Materials
AN-1410	Rapid Determination of Sr in Emergency Urine Samples	Sr	Emergency Urine
AN-1411	Rapid Determination of Sr in Emergency Water Samples	Sr	Emergency Water
AN-1412	Rapid Determination of Actinides in Emergency Urine Samples	U, Pu, Am, Cm	Emergency Urine
AN-1413	Rapid Determination of Actinides in Emergency Water Samples	U, Pu, Am, Cm	Emergency Water
AN-1414	Rapid Determination of ⁹⁰ Sr in Up to 40 Liter Seawater Samples	⁹⁰ Sr	Seawater (40L)
AN-1415	Rapid Determination of ²¹⁰ Po in Water Samples	²¹⁰ Po	Water
AN-1416	Rapid Determination of Actinides and ²¹⁰ Po in Water	²¹⁰ Po, U, Pu, Am, Cm	Water
AN-1417	Rapid Determination of ^{226/228} Ra in Water Samples	^{226/228} Ra	Water
AN-1418	Rapid Determination of ²²⁶ Ra in Water Samples	²²⁶ Ra	Water
AN-1419	Rapid Determination of ²²⁶ Ra in Concrete and Brick	²²⁶ Ra	Concrete/Brick
AN-1420	Rapid Determination of ²²⁶ Ra in Glass Fiber Air Filters	²²⁶ Ra	Glass Fiber Air Filters
AN-1421	Rapid Determination of ²²⁶ Ra in 1g Soil Samples	²²⁶ Ra	Soil
AN-1422	Rapid Determination of ²²⁶ Ra in 5g Vegetation Samples	²²⁶ Ra	Vegetation
AN-1423	Rapid Determination of Pu, Np, and U in 1-8L Seawater Samples	Np, Pu, U	Seawater (1-8L)
AN-1424	Rapid Determination of Pu, Am and Cm in 80L Seawater Samples	Pu, Am, Cm	Seawater (80L)
AN-1425	Rapid Determination of Actinides in 10g Emergency Food Samples	U, Pu, Np, Am, Cm	Emergency Food (10g)
AN-1426	Rapid Determination of Actinides in 100g Emergency Food Samples	U, Pu, Np, Am, Cm	Emergency Food (100g)
<u>Coming Soon</u>			
AN-1427	Rapid Method for Plutonium in Rice	Pu	Rice
AN-1428	Rapid Method for Actinides in fecal samples	Am, Pu, U	Fecal Samples
AN-1429	Rapid Method for Actinides in Asphalt samples	Am, Pu, U	Asphalt
AN-1430	Rapid Method for Actinides in Emergency Soil Samples	Am, Cm, Pu, U	Emergency Soil
AN-1431	Rapid Method for Determination of Pu, Am, Cm in large soil samples	Pu, Am, Cm	100-200g Soil
AN-1432	Rapid Method for Actinides in Emergency Concrete and Brick Samples	Am, Cm, Pu, Np, U	Concrete and Brick
AN-1433	Rapid Method for actinides in emergency air filter samples	Pu, Am, U	Air Filters
AN-1434	Rapid Method for strontium in emergency air filter samples	Sr	Air Filters
AN-1435	Rapid Method for Np and Pu in large soil samples	Np, Pu	20-50g Soil
AN-1436	Rapid Method for Np and Pu in large soil samples by ICP-MS	Np, Pu	20-75g Soil
AN-1437	Rapid Method for Actinides in Urine by ICP-MS and alpha spectrometry	Am, Cu, Pu, Np, U	urine
AN-1438	Rapid Method for Np and Pu in water by ICP-MS and alpha specrometry	Np, Pu	water


Sample Preparation



- 1) Acidify with HNO_3 (pH 1-2).
- 2) Add Tracers.
- 3) Heat. Mix.
- 4) Remove from heat.
- 5) Add $\text{Ca}/(\text{NH}_4)_2\text{HPO}_3$.
- 6) Mix. Add NH_4OH to pH 8-9
- 7) Additional heating.
- 8) Cool. Settle. Centrifuge.
- 9) Wash ppt. w/ dilute NH_4OH .



- 1) HNO_3 protonates PO_4^{3-} , dissolving calcium phosphate ppt.
- 2) $\text{Al}(\text{NO}_3)_3 - \text{HNO}_3$ provide NO_3^- to promote extraction of nitrate-complexes of actinides/Sr.





$\text{Al}(\text{NO}_3)_3\text{-HNO}_3$

- 1) HNO_3 protonates PO_4^{3-} , dissolving calcium phosphate ppt.
- 2) $\text{Al}(\text{NO}_3)_3 - \text{HNO}_3$ provide NO_3^- to promote extraction of nitrate-complexes of actinides/Sr.
- 3) Al^{3+} complexes phosphate, minimizing impact on Th^{4+} .

Valence Adjustment


1) Cool to room temperature. Warm solutions can increase unwanted reactions.

2) Add ferric iron and sulfamic acid. Mix.



3) Add ascorbic acid to reduce Fe(II), Pu(III), Np(IV). Uranium remains U(VI). Am/Cm(III). Sr(II).

Ascorbic Acid Oxidation



4) Add NaNO_2 to oxidize Pu(III) to Pu(IV).

I.L. Jenkins, "Factors governing the choice of $^{237}\text{Np}/^{238}\text{Pu}$ Separation Process, *Actinides Reviews*, 1, 187-211 (1969)

Column/Cartridge Separations

- 1) Precondition w/ 2-3 bed volumes of 3M HNO₃.
- 2) Load/Rinse at 2mL/min.
- 3) OK for cartridges to run dry.
- 4) New/clean reservoirs and vacuum box tips.
- 5) Strip analytes at 1mL/min.

Alpha Spectrometry Source Preparation

Microprecipitation

Adequate resolution

Rapid/Many Samples

Direct Prep. from
Column Strip Solution

Most Routine
Analytical Sources

Electrodeposition

Superior resolution

Slower

Requires evaporation
to remove HNO₃/HCl.



Calibration Sources
Nuclides with difficult to
resolve peaks


Need Geometry Match to Calibration Source

M. Schultz, NAMP Actinide Chemistry Webinar Series, "Source Preparation for alpha spectrometry,"
https://www.icln.org/default/assets//File/Source%20Prep%20Alpha%20Spec%20Final_1-21-13%20slide%20deck.pdf

ASTM Method C1284-10 "Standard Practice for Electrodeposition of the Actinides for Alpha Spectrometry"


C. W. Sill, "Precipitation of actinides as fluorides or hydroxides for high-resolution alpha spectrometry," *Nuclear and Chemical Waste Management*, 7, 201-215 (1987).






CeF₃ Microprecipitation of U-238





20mL 1M HCl, 1mL conc. HF, 0.5mL 10% TiCl₃











CeF₃ Microprecipitation of U-238

20mL 1M HCl, 1mL conc. HF, 0.5mL 10% TiCl₃










CeF₃ Microprecipitation of Th-232

30mL 4M HCl, 3mL conc. HF



CeF₃ Microprecipitation of Th-232

20mL 1M HCl, 1mL conc. HF, 0.5mL 10% TiCl₃

