



# Characterisation of a Cu selective resin for use in the production of Cu isotopes



# Outline

- Resin characterisation
  - Selectivity
  - Interferences with Ni or Zn
- Column Experiments
  - Optimisation
  - Simulated targets
  - Decontamination factors
- Conversion of Cu resin eluate via AIX
- Other applications

# D<sub>w</sub> Cu resin in HCl




Figure 1: D<sub>w</sub> of Cu and selected elements on Cu resin in HCl, varying pH values, measurement via ICP-MS, each element 10µg/ml

Loading conditions (> pH 2):

- Good selectivity for Cu over other elements tested,
- Some selectivity for Fe at high pH

Elution:

- Low Cu D<sub>w</sub> at low pH, further tests with 4, 6 and 8 M HCl

# $D_w$ Cu resin in $HNO_3$




Figure 2:  $D_w$  of Cu and selected elements on Cu resin in  $HNO_3$ , varying pH values, measurement via ICP-MS, each element  $10\mu g/ml$

Loading conditions (> pH 2):

- Good selectivity for Cu over all other elements tested

Elution:

- Low Cu  $D_w$  at low pH

# $D_w$ Cu resin in $H_2SO_4$




Figure 3:  $D_w$  of Cu and selected elements on Cu resin in  $H_2SO_4$ , varying pH values, measurement via ICP-MS, each element  $10\mu\text{g/ml}$

Loading conditions ( $> \text{pH } 2$ ):

- Good selectivity for Cu over most other elements tested
- Limited selectivity for Fe

Elution:

- Low Cu  $D_w$  at low pH

# D<sub>w</sub> Cu – interferences - HCl pH 2

## Ni interference




Figure 4: D<sub>w</sub> of Cu on Cu resin in HCl (pH 2) in presence of various amounts of Ni

## Zn interference




Figure 5: D<sub>w</sub> of Cu on Cu resin in HCl (pH 2) in presence of various amounts of Zn

- For high amounts of Ni or Zn D<sub>w</sub>(Cu) remains > 1000 in pH 2 HCl
- For up to 1 g target material per g resin only negligible interference

# Conclusions I

- Stable, high Cu  $D_w$  values at pH > 2 for all acids
- Good selectivity for Cu over tested elements
- No selectivity for Zn or Ni
- No interference on Cu extraction by elevated amounts of Ni or Zn
  - Tested up to 1g Ni or Zn per g of resin
- HCl system chosen for radiopharmaceutical application
  - Loading solution: HCl pH 2 (or higher)
  - Elution with HCl (elevated concentration)

# First approach




Figure 6: Scheme of elution conditions – first approach



# Elution study – 2 mL CU resin column

## Elution with 6 M HCl





Figure 7/8 : Elution study; simulated Ni targets (200 mg Ni, traces of Zn, Cu), elution with 6M HCl and 8M HCl

## Elution with 8 M HCl




- 8M HCl allows for elution in smaller volume
- Clean Cu fraction

# Elution study - irradiated Ni target

- Irradiation of a Ni foil (10 mg) :
  - ( $\varnothing = 13$  mm, 0.025 mm thickness)
    - Cyclotron BC1710 at Forschungs Zentrum Jülich
    - $E_p = 15$  MeV
    - 1 h; 0.5 mA
- Additionally 170 mg non-irradiated Ni-foil added

# Elution study – irradiated Ni target



~100 % Cu in E1  
~100 % Ni in L1 & R1

## Loading:

L1 : 5 mL HCl pH 2

## Rinsing:

R1 : 10 mL HCl pH 2

R2 : 10 mL HCl pH 2


## Elution:

E1 : 5 mL HCl 8 M


E2 : 5 mL HCl 8 M

Figure 9 : Elution study, irradiated Ni<sub>nat</sub> target


- Very pure Cu fraction
- Elution and rinsing volumes to be optimized



**Eluate from loading : L 1  
5 mL pH 2 HCl**



Rinsing : R1  
10 mL HCl pH 2



Elution : E1  
5 mL HCl 8M

# Optimization of Cu elution volume





Figure 10/11: Elution study, simulated Ni target, varying Cu resin amounts, gravity flow

- Use of smaller columns allows reducing Cu elution volume
  - Near quantitative Cu recovery in ~1 mL 8M HCl
- No impact on Cu purity

# Optimization of Cu elution volume – vacuum-assisted flow





Figure 12/13 : Elution study, simulated Ni and Zn target, 350 mg Cu resin

- Flow rate: 1 mL/min (rinsing up to 6 mL/min)
- 0.35g CU resin columns
- > 90% recovery in 1 mL 8M HCl, quantitative recovery in 1.5 mL
- Simulated Ni and Zn targets (200 mg target material)
- No impact on purity



# Optimized method




Figure 14: Scheme of elution conditions – optimized method

# Decontamination factors $D_f$

- Flow rate: 1 mL/min (rinsing 6 mL/min)
- 0.35 g columns
- Loading solution: Ni, Zn, Co, Ga and Au in 5 mL HCl pH 2
- Separation following optimized method
- ICP-MS measurement
- Calculation of decontamination factors  $D_f$  for Cu fractions
  - Fraction E1 (0.5 mL 8M HCl):
    - $D_f$ : Ni, Co & Zn > 20 000
    - $D_f$ : Au & Ga > 10 000
  - Fraction E2 (0.5 mL 8M HCl):
    - $D_f$  : Ni > 20 000, Co > 40 000, Zn > 70 000, Au > 50 000, Ga > 10 000

# Conclusion II

- 350 mg columns allow for Cu elution in small volume
- Vacuum assisted flow
- Rapid separation
- Quantitative recovery of Cu in 1 – 1,5 mL 8M HCl
  - Cu yield > 90% in 1 mL 8M HCl
  - 97,6% ± 2,3% (k = 1, N=25) in 1.5 mL 8M HCl
- Pure Cu fraction
  - D<sub>f</sub> (ICP-MS)
  - $\gamma$ -spectrometry
- Obtained Cu suitable for labelling (ARRONAX)
- Ni recovered in small volume of 8M HCl
  - 10 – 13 mL load and rinse
  - Further purification for reuse e.g. via direct load on AIX

# Conversion of Cu eluate

- Aim: recovery of Cu in small volume of dilute HCl, water or NaCl solution
- Anion exchange resins (AIX) shows necessary selectivity
- Cu eluate (1 – 1,5 mL 8M HCl) from Cu resin column directly loaded onto small AIX column
- Rinse with 8M HCl
- Elution with deion. water

# Elution study - AIX





Figure 15: Elution study, AIX, various elements, 400 mg A8 Resin (200 – 400 mesh, Eichrom Technologies)

- 400 mg AIX (A8)
- Cu elution in 0.6 – 0.8 mL water
- Anion exchange conversion step gives additional decontamination from Ni, Zn, Au

# Cu fraction conversion step



➤ Direct load of Cu eluate onto anion exchange resin / elution with water or saline solution

➤ Further purification from matrix elements and organics

➤ Cu yield in 0.7 mL water:  
 $93,8 \pm 6,4\%$  ( $k = 1, N=25$ )

Fig 16: Conversion step using anion exchange, 400 mg A8 resin (Eichrom Technologies), 200 – 400 mesh

# Full method

- Vacuum-assisted flow
- CU resin (350 mg)
  - Load from 1 - 2 mL HCl pH2
  - Rinse with 5 mL and 3 mL HCl pH 2
  - Load and rinse contain ~100% Ni (→ Ni-64 recovery)
  - Cu elution in 1 – 1.5 mL 8M HCl
- Cu yield > 95%; high decontamination factors
- Conversion on AIX (400 mg)
  - Load from 1 – 1.5 mL 8M HCl
  - Rinse with 0.6 mL 8M HCl and < 0.2 mL water
  - Cu elution in 0.6 – 0.8 mL water (or saline solution)
- Cu yield > 90%; add. decontamination (Ni, Zn, Au, organics)
- Overall separation time: <10 minutes

# Other applications

- Removal of trace Cu-64 before Ni-64 plating
- Cu concentration and purification for analytical purpose (e.g. Cu in sea water)




Fig 17: Elution study, 10 mL spiked sea water, 350 mg CU resin column, vacuum assisted flow

- 10 mL sea water (pH 2.3)
- Cu yield > 95% in 1 mL 8M HCl
- Pure Cu fraction