

AN-1411-10

Rapid Determination of Sr in Emergency Water Samples

Summary of Method Strontium is separated and concentrated from up to 400mL water samples using calcium phosphate precipitation. The precipitate is dissolved in nitric acid and aluminum nitrate. Strontium is then separated from matrix impurities and potentially interfering radionuclides in the sample using a 2mL cartridge of Eichrom Sr Resin. Radiostrontium is measured on a low background gas flow proportional counter or liquid scintillation counter. Chemical yield of strontium is determined by gravimetric recovery of stable strontium or ICP-AES measurement. Typical chemical recovery of strontium is >80%. Measured values of 90Sr agreed to within 14% of reference values for 10 minute count times, although longer count times can be used to improve detection limits and uncertainty. A single operator can complete the separation method for batches of 12-24 samples in as little as 3-4 hours.

Reagents

Sr Resin, 2mL Cartridges (Eichrom SR-R50-S)
Nitric Acid (70%)
Ammonium Hydroxide (listed as 28% NH₃ or 56% NH₄OH)
Deionized Water
1.25M Ca(NO₃)₂
3.2M (NH₄)₂HPO₄
Strontium Carrier (10mg/mL)
2M Al(NO₃)₃
90Sr standard
Oxalic acid

Equipment

Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)
Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)
Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)
Yellow Outer Tips (Eichrom AR-1000-OT)
50mL and 250mL Centrifuge Tubes
Centrifuge
Cupped Stainless Steel Planchets (~5mL volume)
Gas Flow Proportional Counter
Analytical Balance
Vacuum Pump

Figure 1. Sample Preparation

400mL water
+4mg Stable Sr Carrier*
+2mL 1.25M Ca(NO₃)₃
+5mL 3.2 M (NH₄)₂HPO₄.
Mix Well.

Adjust to pH 10 with NH₄OH.

*Amount of Sr carrier may need to be adjusted to account for native Sr content in sample.

Centrifuge 3500 rpm, 10 min.

Mix Well.


Decant Supernate To Waste

Dissolve precipitate/residue in 8mL 6M HNO₃ and 8mL 2M Al(NO₃)₃.




Figure 2. Load Solution Preparation and Strontium Separation

- (1) Precondition Sr Resin with 10mL 8M HNO₃.
- (2) Load sample at 1-2mL/min.
- (3) Rinse sample tube with 3mL 3M HNO₃.
- (4) Add tube rinse to Sr Resin. Elute at 1-2mL/min.
- (5) Rinse Sr Resin sequentially with:
- 10 mL 8M HNO₃
- 5mL 3M HNO₃ 0.05 oxalic acid
- 5mL 8M HNO₃
- (6) Dispose of (1) to (5) as waste.
- (7) Strip Sr with 15mL 0.05M HNO₃ at 1mL/min.



Gas Flow Proportional Counting:*

- (8) Evaporate samples to dryness on tared cupped stainless steel planchets.
- (9) Rinse Sr sample vials with 2mL 0.05M HNO₃. Transfer vial rinse to planchets. Evaporate to dryness.



- (10) Weigh planchets on an analytical balance to determine gravimetric yield of stable Sr(NO₃)₂.
- (11) Measure radiostrontium in samples on low background gas flow proportional counter.

- *(Options for 89/90Sr Descrimination)
- (a) Sr fraction from step (7) can be transferred to a liquid scintillation vial. ⁸⁹Sr can be measured by Cerenkov counting (no LSC cocktail). ^{89/90}Sr may then be measured after adding liquid scintillation cocktail.
- (b) Sr fraction from step (10) can be dissolved in 10mL 8M HNO₃ after >7 days of ⁹⁰Y ingrowth. ^{89/90}Sr can be removed on Sr Resin. ⁹⁰Y will elute in Sr Resin load and can be counted by liquid scintillation or gas flow proportional counting.

Actinides may also be measured by adding a 2mL TEVA, TRU and DGA cartridges above Sr Resin and following the separation scheme in Eichrom application note AN-1413, "Rapid Determination of Actinides in Emergency Water Samples."

References

1) Sherrod L. Maxwell, Brian K. Culligan, "Rapid separation method for emergency water and urine samples," *J. Radioanal. Nucl. Chem.*, 279(3), 901-907 (2009).