

MnO₂ Resin

A new approach for Radium separation

UGM05 - 09/12/05 - Manchester

Outline

- Radium
- pH study
- Reaction time
- Amount of resin
- Salt concentration effect
- Flow-rate study
- Conclusion

Radium

- ^{226}Ra (1620 y), ^{228}Ra (5,28 y)
- ^{223}Ra (11,4 d), ^{224}Ra (3,6 d)
 - Water analysis for radioprotection and geological tracing
- Same metabolic behaviour as Ca

Ra separations used

- Separation with TRU or LN resin:
 - Isolation of ^{228}Ra daughter ^{228}Ac .
 - ^{228}Ac measured by γ -spectrometry (338 and 911 keV)


Use of MnO₂

- Resin MnO₂ 80-160 μm particle size
- Use of de-ionised , ground and sea waters
- Study with variation of
 - pH
 - Kinetics
 - Ratio mass of resin vs volume of solution
 - Salinity
 - Flow rates
- Use of ¹³³Ba to monitor Ra behaviour

pH variation

- 10 mL de-ionized water spiked with ^{133}Ba
- 25 mg MnO_2 resin
- Magnetic stirrer
- pH adjusted with HCl or NaOH
- 60 minutes contact
- 20°C
- Measure of filtered MnO_2 in NaI well-type gamma counter

pH variation




Kinetics

- 10 mL water samples spiked with ^{133}Ba (de-ionized, 0.2, 3.5 and 35‰ salinities)
- 25 mg MnO_2 resin
- Magnetic stirrer
- pH = 7.0
- 1.0 to 90 minutes contact
- Measure of filtered MnO_2 in NaI well-type gamma counter

Kinetics



Charged sample



Ratio V_s/V_m

- **Batch experiment**
 - 25 mg MnO₂/10 mL solution (25 g/1.0 L)
 - 10 mg MnO₂/10 mL solution (10 g/1.0 L)
- **Column**
 - 1 g MnO₂ resin
 - Geometry: $\varnothing_i = 0.9$ cm, $H = 6,5$ cm
 - Water samples: de-ionised and artificial seawater

Ratio V_s/V_m




Flow rate


- 1.0 L sample water spiked with ^{133}Ba and ^{228}Ra
- 1.0 g MnO_2 resin cartridge
- pH = 7.0
- Peristaltic pump, flow rate 0 -300 mL/min
- Samples: de-ionized and seawater

Flow rates - results


^{133}Ba



^{228}Ra



Flow rates – results: $^{133}\text{Ba}/^{228}\text{Ra}$



Review

- **pH window : 4-8 for 90-100% recovery**
- **Kinetics dependant of salinity**
- **Use of columns**
- **Flow rate dependant of salinity**
 - **20 mL/min OK for all types of water samples**

Results: ^{226}Ra and ^{228}Ra in water

- 1.25 g MnO_2 resin cartridges
- 1.5 L samples spiked with ^{133}Ba
- pH = 6-7
- Flow rate ~20 mL/min

Procedure

- Load of sample on MnO₂
- Strip MnO₂ with 15 mL 4 M HCl-1.5% H₂O₂
- Load 15 mL strip solution from MnO₂ resin on TRU (1 mL)-DGA, Normal (2 mL) stacked
 - TRU (U, Th removal), DGA (²²⁸Ac retention)
 - Collect load + 5 mL 4 M HCl rinse for ²²⁶Ra
 - Discard TRU resin
 - Strip ²²⁸Ac from DGA with 10 mL 0.5 M HCl
 - Micro-precipitate with CeF₃
- α-spectrometry

Summary

- **Rapid and efficient method for Ra separation**
- **Applicable to a variety of water samples**
- **Use of large initial samples**