

URANIUM AND THORIUM IN WATER

1. SCOPE

- 1.1. This is a method for the separation of uranium and thorium from water samples. After completing this method, source preparation for measurement of uranium and thorium by alpha spectrometry is performed by electrolytic deposition onto stainless steel planchets (Eichrom Method SPA02) or by rare earth fluoride microprecipitation onto polypropylene filters (Eichrom Method SPA01).
- 1.2. This method does not address all aspects of safety, quality control, calibration or instrument set-up. However, enough detail is given for a trained radiochemist to achieve accurate and precise results for the analysis of the analyte(s) from the appropriate matrix, when incorporating the appropriate agency or laboratory safety, quality and laboratory control standards.

2. SUMMARY OF METHOD

- 2.1. Uranium and thorium are separated by Eichrom TEVA and UTEVA resins prior to measurement by alpha spectrometry. A calcium phosphate precipitation can be used to concentrate actinides from water samples. Tracers are used to monitor chemical recoveries and correct results to improve precision and accuracy.

3. SIGNIFICANCE OF USE

- 3.1. This is a rapid, reliable method for measurement of actinides in water samples that is more cost-effective and efficient than traditional ion exchange, solvent extraction and precipitation techniques.

4. INTERFERENCES

- 4.1. Radionuclides with unresolvable alpha energies such as ^{241}Am and ^{238}Pu , ^{237}Np and ^{234}U , or ^{210}Po and ^{232}U must be chemically separated to enable measurement. This method separates these isotopes effectively.
- 4.2. Very high levels of phosphate in the sample may cause reduced recovery of actinides during calcium phosphate precipitation and

column separations. Adjusting the amount of phosphate added to co-precipitate the actinides may be necessary in these cases.

- 4.3. The sample preparation procedure outlined in this method will adequately recover actinides from freshly collected, well preserved, homogenous water samples. Older, poorly preserved samples or samples with significant organic or solid matter may require more aggressive treatment to recover actinides which have precipitated or adsorbed to the walls of the storage container or solid matter. Rinsing the empty storage container with warm HNO₃, adjusting the HNO₃ concentration of the sample to 1M HNO₃ and boiling, and/or wet-washing the calcium phosphate precipitate may be required for older, poorly preserved samples.

5. APPARATUS

- Analytical balance, 0.0001 g sensitivity
- Beakers, glass
- Centrifuge tubes, 50mL and 250mL
- Centrifuge, with rotor and carriers for 50mL and 250mL tubes
- Column rack, Eichrom Part: AC-103
- Extension funnels, 25 mL, Eichrom Part: AC-120
- Fume hood
- Hotplate
- Stir rods, glass
- Vortex mixer

6. REAGENTS

Note: Analytical grade or ACS grade reagents are recommended. Evaluation of key reagents, such as aluminum nitrate and ammonium hydrogen phosphate, for contribution to method background levels from naturally occurring radioactive materials is recommended.

Aluminum nitrate nonahydrate, $Al(NO_3)_3 \cdot 9H_2O$
Ammonium hydrogen phosphate, $(NH_4)_2HPO_4$
Ammonium hydroxide(57% NH_4OH or 28% NH_3), concentrated NH_4OH
Ammonium oxalate monohydrate, $(NH_4)_2C_2O_4 \cdot H_2O$
Appropriate tracers or standards (U-232, Th-229)
Ascorbic acid powder, $C_6H_8O_6$
Calcium nitrate, $CaNO_3$

<i>Deionized water, All reagents are prepared with deionized water</i>
<i>Hydrochloric acid (37%), concentrated HCl</i>
<i>Hydrogen peroxide (30%), concentrated H₂O₂</i>
<i>Isopropyl alcohol, C₃H₇OH</i>
<i>Nitric acid (70%), concentrated HNO₃</i>
<i>Oxalic acid dihydrate, H₂C₂O₄·2H₂O</i>
<i>Phenolphthalein pH indicator</i>
<i>TEVA[®] resin, 2mL prepacked column, 100-150μm, Eichrom Part TE-C50-A</i>
<i>UTEVA[®] resin, 2mL prepacked column, 100-150μm, Eichrom Part UT-C50-A</i>

- 6.1. *Ammonium hydrogen phosphate (3.2M)*- Dissolve 104g of (NH₄)₂HPO₄ in 200mL of water, heat gently to dissolve. Dilute to 250mL with water.
- 6.2. *Calcium nitrate (1.25M)* - Dissolve 51g of Ca(NO₃)₂ in 100mL of water and dilute to 250mL with water.
- 6.3. *Hydrochloric acid (9M)* - Add 750mL of concentrated HCl to 10mL of water and dilute to 1L with water.
- 6.4. *Hydrochloric acid (5M) - oxalic acid (0.05M) solution*- Dissolve 6.3g oxalic acid dihydrate in 400mL water. Add 417mL concentrated HCl. Cool to room temperature. Dilute to 1L with water.
- 6.5. *Nitric acid solution (3M)* - Add 188mL of concentrated HNO₃ to 700mL of water. Dilute to 1L with water.
- 6.6. *Nitric acid (3M) - Aluminum nitrate (1M) solution* - Dissolve 375g of Al(NO₃)₃·9H₂O in 500mL of water, add 188mL of concentrated HNO₃ and dilute to 1L with water.
- 6.7. *Phenolphthalein solution* - dissolve 1g phenolphthalein in 100mL 95% isopropyl alcohol and dilute with 100mL of water.

7. PROCEDURE

7.1. Water Sample Preparation:

- 7.1.1. If samples larger than 1L are analyzed, evaporate the sample to approximately 1L.

7.1.2. Aliquot 500 to 1000mL of the filtered sample (or enough to meet required detection limit) into an appropriate size beaker.

7.1.3. Add 5mL concentrated HNO₃.

Note: If using self-cleaning ²³²U tracer (Eichrom Method TP01), mix for 1-2 minutes to suspend the BaSO₄ precipitate to remove ²²⁸Th and daughters, centrifuge and then take aliquot for Uranium tracing.

7.1.4. Add appropriate tracers per lab protocol.

7.1.5. Add 0.5mL of 1.25M Ca(NO₃)₂ to each beaker.

7.1.6. Place each beaker on a hot plate.

7.1.7. Cover each beaker with a watch glass.

7.1.8. Heat samples at medium setting for 30-60 minutes.

7.1.9. Remove watch glass and turn the heat down.

7.1.10. Add 2-3 drops of phenolphthalein indicator and 1mL of 3.2 M (NH₄)₂HPO₄ solution.

7.1.11. While stirring, slowly add enough concentrated NH₄OH with a squeeze bottle to reach the phenolphthalein end point and form a calcium phosphate precipitate. Allow the sample to heat for another 20-30 minutes.

7.1.12. Remove samples from hot plate, cool samples to room temperature, and allow precipitate to settle until solution can be decanted (30 minutes to 2 hours) or centrifuge.

7.1.13. Decant supernate and discard to waste.

7.1.14. Transfer the precipitate to a centrifuge tube and centrifuge the precipitate for approximately 10 minutes at 2000 rpm.

7.1.15. Decant supernate and discard to waste.

7.1.16. Wash the precipitate with an amount of water approximately twice the volume of the precipitate. Mix well on a vortex mixer. Centrifuge for 5-10 minutes. Discard the supernate.

7.1.17. If an ammonia odor persists, repeat water wash of precipitate.

7.1.18. Dissolve precipitate in 5mL of concentrated nitric acid. Transfer solution to a 100mL beaker. Rinse centrifuge tube with 2-3 mL

of concentrated nitric acid and transfer to beaker. Evaporate solution to dryness.

7.2. Actinide Separations using Eichrom Resins:

7.2.1. Dissolve each precipitate with 10mL of 3M HNO₃-1M Al(NO₃)₃.

Note: An additional 5-10mL may be necessary if the volume of precipitate is large.

7.2.2. Th separation using TEVA Resin

7.2.2.1. For each sample dissolved, place a TEVA Resin column in the column rack.

7.2.2.2. Place a beaker below each column, remove the bottom plug from each column, push top frit down to the top of the resin bed, and allow to drain. Attach column reservoirs to each column.

7.2.2.3. Add 5mL of 3M HNO₃ into each column to condition resin. Allow solution to drain.

7.2.2.4. Place a clean, labeled beaker below each column.

7.2.2.5. Transfer each redissolved sample into the appropriate TEVA Resin column by pouring or by using a plastic transfer pipet. Collect eluate.

7.2.2.6. Add 5mL of 3M HNO₃ to rinse to each beaker. Transfer each rinse solution into the appropriate TEVA Resin column reservoir. Collect eluate.

7.2.2.7. Add 10mL of 3M HNO₃ to each TEVA column. Collect eluate. Set beakers aside for uranium separation.

7.2.2.8. Place a new beaker below each TEVA column. Add 15mL of 3M HNO₃ into each column. Discard eluate as waste.

Note: Residual Uranium is removed with the extended 3M HNO₃ rinse.

7.2.2.9. Place a clean, labeled 50mL beaker below each column. Make sure all traces of the previous 3M HNO₃ rinse have been removed from the column reservoir. If necessary, the reservoir (but not the TEVA column itself) can be rinsed with a small volume of water.

7.2.2.10. Add 15mL of 9M HCl into each column and collect eluate.

Note: 9M HCl will strip Th from the column. Pu⁴⁺ and Np⁴⁺ are retained on the column.

7.2.2.11. Set Th samples aside for alpha source preparation.

7.2.3. Uranium separation using UTEVA Resin

- 7.2.3.1. For each sample solution, place a UTEVA Resin column in the column rack.
- 7.2.3.2. Place a beaker below each column, remove the bottom plug and caps from each column, push the top frit down to the top of the resin bed, and allow to drain. Attach column reservoirs to each column.
- 7.2.3.3. Add 5mL of 3M HNO₃ into each column reservoir to condition resin. Allow solution to drain.
- 7.2.3.4. Transfer each solution from step 7.2.2.7 into the appropriate UTEVA Resin column. Allow solution to drain.
- 7.2.3.5. Add 5mL of 3M HNO₃ to rinse to each beaker. Transfer each rinse solution into the appropriate UTEVA Resin column reservoir. Allow solution to drain.
- 7.2.3.6. Add 5mL of 3M HNO₃ into each column. Allow solution to drain.
- 7.2.3.7. Add 5mL of 9M HCl into each column reservoir. Allow solution to drain.

Note: This rinse converts the resin to the chloride system.

- 7.2.3.8. Add 20mL of 5M HCl-0.05M oxalic acid into each column reservoir. Allow solution to drain. Discard the combined eluate to this point as waste.

Note: This rinse removes any traces of plutonium, neptunium and thorium from the column.

- 7.2.3.9. Place a clean, labeled beaker below each column.
- 7.2.3.10. Add 15mL of 1M HCl into each column to strip the uranium. Allow solution to drain.
- 7.2.3.11. Set U samples aside for alpha source preparation.

7.3. Prepare samples for actinide measurement by alpha spectrometry using electrodeposition (Eichrom method SP02) or rare earth fluoride microprecipitation (Eichrom method SP01).

8. CALCULATIONS

Calculate the actinide activity as follows:

Calculate tracer yield:

$$\text{Yield} = \frac{(C_s - B_s)}{E_s \times A_s}$$

where:

C_s = measured actinide tracer, cpm

B_s = background, cpm

E_s = counting efficiency for tracer

A_s = tracer activity, dpm

Note: If any tracer may be present in the sample, a spiked and unspiked sample must be analyzed to determine chemical yield, where:

$$\text{Yield} = \frac{(\text{spiked sample tracer cpm} - \text{unspiked sample tracer cpm})}{E \times \text{actinide spike activity dpm}}$$

$$\text{Percent yield} = \text{Yield} \times 100$$

Calculate actinide isotope activity:

$$\text{Sample dpm / g or dpm / L} = \frac{S - B}{E \times V \times Y}$$

where:

S = sample activity, cpm

B = background, cpm

E = counting efficiency = measured cpm/dpm of isotopic standard
 V = sample weight, g or volume, L

Y = yield

Conversion of dpm/g to pCi/g:

$$\text{pCi/g} = (\text{dpm/g})/2.22$$

9. REFERENCES

- 1) Horwitz, E.P., et al. "Separation and Preconcentration of Actinides from Acidic Media by Extraction Chromatography," *Analytica Chimica Acta*, 281, 361-372 (1993).
- 2) Horwitz, E.P., et al. "Separation and Preconcentration of Uranium from Acidic Media by Extraction Chromatography," *Analytica Chimica Acta*, 266, 25-37 (1992).
- 3) ASTM Method D3972-09, "Standard Test Method for Isotopic Uranium in Water by Radiochemistry."
- 4) ASTM Method D7282-06, "Standard Practice for Set-up, Calibration, and Quality Control of Instruments Used for Radioactivity Measurements."
- 5) ASTM Method D3648-14, "Standard Practices for the Measurement of Radioactivity."

