

AMERICIUM, PLUTONIUM, AND URANIUM IN WATER

1. SCOPE

- 1.1. This is a method for the separation of americium, plutonium and uranium from up to 1L of water. After completing this method, source preparation for measurement of americium, plutonium and uranium by alpha spectrometry is performed by electrolytic deposition onto stainless steel planchets (Eichrom Method SPA02) or by rare earth fluoride micro precipitation onto polypropylene filters (Eichrom Method SPA01).
- 1.2. This method does not address all aspects of safety, quality control, calibration or instrument set-up. However, enough detail is given for a trained radiochemist to achieve accurate and precise results for the analysis of the analyte(s) from the appropriate matrix, when incorporating the appropriate agency or laboratory safety, quality and laboratory control standards.

2. SUMMARY OF METHOD

- 2.1. Americium, plutonium and uranium are separated by Eichrom UTEVA and TRU resins prior to measurement by alpha spectrometry. A calcium phosphate precipitation can be used to concentrate actinides from water samples. Tracers are used to monitor chemical recoveries and correct results to improve precision and accuracy.

3. SIGNIFICANCE OF USE

- 3.1. This is a rapid, reliable method for measurement of actinides in water samples that is more cost-effective and efficient than traditional ion exchange, solvent extraction and precipitation techniques.

4. INTERFERENCES

- 4.1. Nuclides with unresolvable alpha energies such as ^{241}Am and ^{238}Pu , ^{237}Np and ^{234}U , or ^{232}U and ^{210}Po must be chemically separated to enable measurement. This method separates these isotopes effectively.

- 4.2. Very high levels of phosphate in the sample may cause reduced recoveries of actinides during calcium phosphate precipitation and column separations. Adjusting the amount of phosphate added to co-precipitate the actinides may be necessary in these cases.
- 4.3. The sample preparation procedure outlined in this method will adequately recover actinides from freshly collected, well preserved, homogenous water samples. Older, poorly preserved samples or samples with significant organic or solid matter may require more aggressive treatment to recover actinides which have precipitated or adsorbed to the walls of the storage container or solid matter. Rinsing the empty storage container with warm HNO₃, adjusting the HNO₃ concentration of the sample to 1M HNO₃ and boiling, and/or wet-ashing the calcium phosphate precipitate may be required for older, poorly preserved samples.
- 4.4. This method may also be applied using TEVA-TRU separation chemistry as described in Eichrom Application Note AN-1413.

5. APPARATUS

- Analytical balance, 0.0001 g sensitivity
- Beakers, glass
- Centrifuge tubes, 50mL and 250mL
- Centrifuge, with rotor and carriers for 50mL and 250mL tubes
- Column rack, Eichrom Part: AC-103
- Extension funnels, 25 mL, Eichrom Part: AC-120
- Fume hood
- Hotplate
- Stir rods, glass
- Vortex mixer

6. REAGENTS

Note: Analytical grade or ACS grade reagents are recommended. Evaluation of key reagents, such as aluminum nitrate and ammonium hydrogen phosphate, for contribution to method background levels from naturally occurring radioactive materials is recommended.

Aluminum nitrate nonahydrate, Al(NO ₃) ₃ ·9H ₂ O
Ammonium hydrogen phosphate, (NH ₄) ₂ HPO ₄
Ammonium hydroxide(57% NH ₄ OH or 28% NH ₃), concentrated NH ₄ OH
Ammonium oxalate monohydrate, (NH ₄) ₂ C ₂ O ₄ ·H ₂ O

Ammonium thiocyanate, NH ₄ SCN
Appropriate tracers or standards (U-232, Am-243, Pu-242 or Pu-236)
Ascorbic acid powder, C ₆ H ₈ O ₆
Calcium nitrate, CaNO ₃
Deionized water, All reagents are prepared with deionized water
Ferric nitrate nonahydrate, Fe(NO ₃) ₃ ·9H ₂ O
Hydrochloric acid (37%), concentrated HCl
Hydrofluoric acid (49%), concentrated HF
Hydrogen peroxide (30%), concentrated H ₂ O ₂
Isopropyl alcohol, C ₃ H ₇ OH
Nitric acid (70%), concentrated HNO ₃
Oxalic acid dihydrate, H ₂ C ₂ O ₄ ·2H ₂ O
Phenolphthalein pH indicator
Sodium nitrite, NaNO ₂
Sulfamic acid, H ₃ NSO ₃
TRU [®] resin, 2mL prepacked column, 100-150µm, Eichrom Part TR-C50-A
UTEVA [®] resin, 2mL prepacked column, 100-150µm, Eichrom Part UT-C50-A

- 6.1. *Ammonium bioxalate (0.1M)* - Dissolve 6.31g of oxalic acid dihydrate and 7.11g of ammonium oxalate monohydrate in 900mL of water. Dilute to 1L with water.
- 6.2. *Ammonium hydrogen phosphate (3.2 M)* - Dissolve 104g of (NH₄)₂HPO₄ in 200mL of water, heat gently to dissolve. Dilute to 250mL with water.
- 6.3. *Ammonium thiocyanate indicator (1M)* - Dissolve 7.6g of ammonium thiocyanate in 90mL of water. Dilute to 100mL with water.
- 6.4. *Ascorbic acid (1M)* - Dissolve 17.6g of ascorbic acid in 90mL of water. Dilute to 100mL with water. **Prepare fresh weekly.**
- 6.5. *Calcium nitrate (1.25M)* - Dissolve 51g of Ca(NO₃)₂ in 100mL of water. Dilute to 250mL with water.
- 6.6. *Hydrochloric acid (1M)* - Add 83mL of concentrated HCl to 900mL of water. Dilute to 1L with water.

- 6.7. *Hydrochloric acid (4M)* - Add 333mL of concentrated HCl to 500mL of water. Dilute to 1L with water.
- 6.8. *Hydrochloric acid (4M) - hydrofluoric acid (0.1M)*- Add 333mL of concentrated HCl and 3.6mL of concentrated HF to 500mL of water. Dilute to 1L with water.
- 6.9. Hydrochloric acid (5M) - oxalic acid (0.05M) solution - Dissolve 6.3g oxalic acid dihydrate in 400mL water. Add 417 mL concentrated HCl. Dilute to 1L with water.
- 6.10. Hydrochloric acid (9M) - Add 750mL of concentrated HCl to 100mL of water. Dilute to 1L with water.
- 6.11. Nitric acid (2M) - sodium nitrite (0.1M) solution- Add 31mL of concentrated HNO₃ to 200 mL of water. Dissolve 1.72g of sodium nitrite in the solution. Dilute to 250mL with water. Prepare fresh daily.
- 6.12. Nitric acid solution (8M) - Add 500mL of concentrated HNO₃ to 400mL of water. Dilute to 1L with water.
- 6.13. Nitric acid solution (3M) - Add 188mL of concentrated HNO₃ to 700mL of water. Dilute to 1L with water.
- 6.14. Nitric acid (3M) - Aluminum nitrate (1M) solution - Dissolve 375g of Al(NO₃)₃·9H₂O in 500mL of water. Add 188mL of concentrated HNO₃. Dilute to 1L with water.
- 6.15. Ferric Nitrate Solution (5 mg/mL Fe) in 0.1M HNO₃ - To a 500mL volumetric flask, add 18g Fe(NO₃)₃·9H₂O, 400mL of water and 3.12mL of concentrated HNO₃. Swirl to dissolve. Dilute to 500mL with water.
- 6.16. Sulfamic acid (1.5M) - In a 500mL volumetric flask, add 73g of sulfamic acid to 400mL of water. Swirl to dissolve. Dilute to 500mL with water.
- 6.17. Phenolphthalein solution - dissolve 1g phenolphthalein in 100 mL 95% isopropyl alcohol and dilute with 100 mL of water.

7. PROCEDURE

7.1. Water Sample Preparation:

- 7.1.1. Aliquot 500 to 1000mL of the sample (or enough to meet required detection limit) into an appropriate size beaker. If

samples larger than 1L are analyzed, evaporate the sample to approximately 1L.

7.1.2. Add 5mL concentrated HNO₃.

Note: If using self-cleaning ²³²U tracer (Eichrom TP01), vortex mix and centrifuge standard to ensure that ²²⁸Th and its daughters are effectively removed from ²³²U by the BaSO₄ precipitate.

7.1.3. Add appropriate tracers per lab protocol.

7.1.4. Calcium phosphate precipitation:

7.1.4.1. Add 0.5 mL of 1.25M Ca(NO₃)₂ to each beaker.

7.1.4.2. Place each beaker on a hot plate.

7.1.4.3. Cover each beaker with a watch glass.

7.1.4.4. Heat at medium setting for 30-60 minutes.

7.1.4.5. Remove the watch glass and turn the heat down.

7.1.4.6. Add 2-3 drops of phenolphthalein indicator and 1mL of 3.2M (NH₄)₂HPO₄ solution.

7.1.4.7. While stirring, slowly add enough concentrated NH₄OH to reach the phenolphthalein end point (pH~9) and form a calcium phosphate precipitate. Allow the sample to heat for another 20-30 minutes.

7.1.4.8. Remove samples from the hot plate, cool to room temperature, and allow precipitate to settle until solution can be decanted (30 minutes to 2 hours) or centrifuge.

7.1.4.9. Decant supernate and discard to waste.

7.1.4.10. Transfer the precipitate to a centrifuge tube and centrifuge the precipitate for approximately 10 minutes at 2000 rpm.

7.1.4.11. Decant supernate and discard to waste.

7.1.4.12. Wash the precipitate with an amount of water approximately twice the volume of the precipitate. Mix well on a vortex mixer. Centrifuge for 5-10 minutes. Discard the supernate.

7.1.4.13. If an ammonia odor persists, repeat water wash of precipitate.

7.1.4.14. Dissolve precipitate in 5mL concentrated nitric acid. Transfer solution to a 100mL beaker. Rinse centrifuge tube with 2-3 mL of concentrated nitric acid and transfer to beaker. Evaporate the solution to dryness.

7.2. Actinide Separations using Eichrom Resins:

7.2.1. Dissolve each precipitate with 10mL of 3M HNO₃- 1M Al(NO₃)₃.

Note: An additional 5-10mL may be necessary if the volume of precipitate is large.

Note: Make sure that all reagents and the load solution have cooled to room temperature. Warm solutions can cause reactions that will affect oxidation adjustments performed in the following steps.

7.2.2. Add 1mL of 1.5 M sulfamic acid and 0.5mL of 5mgFe/mL to each solution. Swirl to mix.

Note: If the additional 5-10 mL was used to dissolve the sample in step 7.2.1., add a proportionally larger volume of sulfamic acid and Fe solution.

7.2.3. Add 1 drop of 1M ammonium thiocyanate indicator to each sample and mix.

Note: The color of the solution turns deep red, due to the presence of iron in Fe (III).

7.2.4. Add 1mL of 1M ascorbic acid to each solution, swirling to mix. Wait for 2-3 minutes.

Note: The red color should disappear which indicates reduction of Fe(III) to Fe(II). If the red color still persists then additional ascorbic acid solution should be added drop-wise with mixing until the red color disappears. Once the red color disappears, add an additional 0.5mL of 1M ascorbic acid.

Note: If particles are observed suspended in the solution, centrifuge the sample. The supernatant will be transferred to the column in step 7.2.5.5. The precipitate will be discarded.

7.2.5. U separation from Pu, Am using UTEVA Resin

7.2.5.1. For each sample solution, place a UTEVA Resin column in the column rack.

7.2.5.2. Place a beaker below each column, remove the bottom plug and caps from each column, push top frit down to the top of the resin bed, and allow to drain. Attach column reservoirs to each column.

- 7.2.5.3. Add 5mL of 3M HNO₃ into each column to precondition the resin. Allow solution to drain.
- 7.2.5.4. Place a clean, labeled 50 mL beaker below each column.
- 7.2.5.5. Transfer each solution from step 7.2.4. into the appropriate UTEVA Resin column reservoir. Allow solution to drain through column.
- 7.2.5.6. Add 5mL of 3M HNO₃ to rinse to each sample beaker. Transfer each beaker rinse into the appropriate UTEVA Resin column reservoir. Allow solution to drain through column.
- 7.2.5.7. Add 5mL of 3M HNO₃ into each column and collect in the same beaker as in 7.2.5.6.
- 7.2.5.8. Set aside the solutions collected in steps 7.2.5.5, 7.2.5.6, and 7.2.5.7 for Am and Pu separations on TRU.
- 7.2.5.9. Place a new beaker below each UTEVA resin column. Add 15mL of 8M HNO₃ to each UTEVA column. Allow solution to drain.

Note: This rinse will remove Po isotopes, including ²¹⁰Po, which can interfere with the measurement of ²³²U by alpha spectrometry.

- 7.2.5.10. Add 5mL of 9M HCl into each column and allow to drain.

Note: This rinse converts the resin to the chloride system. Some Th may be removed here.

- 7.2.5.11. Add 20mL of 5M HCl- 0.05 M oxalic acid into each column and allow solution to drain.

Note: This rinse removes neptunium and thorium from the column and any residual ferrous ion that might interfere with electrodeposition.

- 7.2.5.12. Place a clean, labeled beaker below each column.

- 7.2.5.13. Add 15mL of 1M HCl into each column to strip the uranium. Allow solution to drain.

- 7.2.5.14. Set U samples aside for alpha source preparation.

7.2.6. Pu, Am Separation Using TRU Resin:

- 7.2.6.1. For each sample dissolved, place a TRU Resin column in the column rack.

- 7.2.6.2. Place a beaker below each TRU column, remove the bottom plug from each column, push the top frit down to the top of the resin bed, and allow each column to drain.
- 7.2.6.3. Add 5mL of 3M HNO₃ into each column to precondition the resin. Allow the solution to drain.
- 7.2.6.4. Transfer solution from step 7.2.5.8. into the appropriate TRU Resin reservoir. Allow the load solution drain through the columns.
- 7.2.6.5. Add 5mL of 3M HNO₃ into the sample beaker and transfer this rinse to the appropriate TRU column reservoir. Allow solution to drain through the columns.
- 7.2.6.6. Add 5mL of 2M HNO₃- 0.1M NaNO₂ into each TRU column. Allow solution to drain through the columns.

Note: Sodium nitrite is used to oxidize Pu(III) to Pu(IV) to enable the Pu/Am separation.

- 7.2.6.7. Discard the load and rinse solutions to waste.
- 7.2.6.8. Ensure that clean, labeled beakers are below each column.
- 7.2.6.9. Make sure that all of the 2M HNO₃-0.1M NaNO₂ is removed from the column reservoir. If necessary, invert each column and rinse the reservoir (but not the column) with a small volume of water.
- 7.2.6.10. Add 15mL of 4M HCl to elute americium. Allow solution to drain through columns.
- 7.2.6.11. Set Am/Cm samples aside for alpha source preparation.
- 7.2.6.12. Place a new beaker below each TRU resin column. Rinse the columns with 25mL of 4M HCl- 0.1 MHF. Discard eluate.

Note: 4M HCl- 0.1 MHF is used to selectively remove any residual Th that may be present on the TRU column.

- 7.2.6.13. Place clean, labeled beakers or centrifuge tubes below each TRU resin column.
- 7.2.6.14. Add 10mL of 0.1M ammonium bioxalate to elute plutonium from each column.
- 7.2.6.15. Set Pu samples aside for alpha source preparation.

- 7.3. Prepare samples for measurement of actinides by alpha spectrometry using electrodeposition (Eichrom SPA02) or rare earth fluoride micro precipitation (Eichrom SPA01).

8. CALCULATIONS

Calculate the actinide activity as follows:

Calculate tracer yield:

$$\text{Yield} = \frac{(C_s - B_s)}{E_s \times A_s}$$

where:

- C_s = measured actinide tracer, cpm
 B_s = background, cpm
 E_s = counting efficiency for tracer
 A_s = tracer activity, dpm

Note: If any tracer may be present in the sample, a spiked and unspiked sample must be analyzed to determine chemical yield, where:

$$\text{Yield} = \frac{(\text{spiked sample tracer cpm} - \text{unspiked sample tracer cpm})}{E \times \text{actinide spike activity dpm}}$$

$$\text{Percent yield} = \text{Yield} \times 100$$

Calculate actinide isotope activity:

$$\text{Sample dpm / g or dpm / L} = \frac{S - B}{E \times V \times Y}$$

where:

- S = sample activity, cpm
 B = background, cpm
 E = counting efficiency = measured cpm/dpm of isotopic standard
 V = sample weight, g or volume, L
 Y = yield

$$\text{Conversion of dpm/g to pCi/g: } \text{pCi/g} = (\text{dpm/g})/2.22$$

9. PRECISION AND BIAS

- 9.1. *Precision*- A relative standard deviation of 4.2% at the 25 dpm level has been reported for uranium. A relative standard deviation of 3.2% at the 60 dpm level has been reported for plutonium.
- 9.2. *Bias*- Mean chemical recoveries of 95% for americium, 93% for plutonium and 86% for uranium have been reported. Since results are corrected based on spike recovery, no significant bias exists for the method.
- 9.3. The data is based on the calcium phosphate precipitation option.

10. REFERENCES

- 1) Horwitz, E.P., et al. "Separation and Preconcentration of Actinides from Acidic Media by Extraction Chromatography," *Analytica Chimica Acta*. 281 (1993), 361-372.
- 2) Horwitz, E.P., et al. "Separation and Preconcentration of Uranium from Acidic Media by Extraction Chromatography," *Analytica Chimica Acta*. 266 (1992), 25-37.
- 3) Maxwell, S.L., "Rapid Column Extraction Methods for Actinides and Sr-89/90 in Water Samples," *Journal of Radioanalytical and Nuclear Chemistry*, 267(3), 537-543 (2006).
- 4) ASTM Method D3972-09, "Standard Test Method for Isotopic Uranium in Water by Radiochemistry."
- 5) ASTM Method D7282-06, "Standard Practice for Set-up, Calibration, and Quality Control of Instruments Used for Radioactivity Measurements."
- 6) ASTM Method D3648-14, "Standard Practices for the Measurement of Radioactivity."
- 7) Eichrom Application Note AN-1413, "Rapid Determination of Actinides in Emergency Water Samples."

