

THORIUM AND NEPTUNIUM IN WATER

1. SCOPE

- 1.1. This is a method for the separation of neptunium and thorium from water samples. After completing this method, source preparation for measurement of thorium and neptunium by alpha spectrometry is performed by electrolytic deposition onto stainless steel planchets (Eichrom Method SPA02) or by rare earth fluoride micro precipitation onto polypropylene filters (Eichrom Method SPA01).
- 1.2. This method does not address all aspects of safety, quality control, calibration or instrument set-up. However, enough detail is given for a trained radiochemist to achieve accurate and precise results for the analysis of the analyte(s) from the appropriate matrix, when incorporating the appropriate agency or laboratory safety, quality and laboratory control standards.

2. SUMMARY OF METHOD

- 2.1. Thorium and neptunium are separated by TEVA Eichrom resin prior to measurement by alpha spectrometry. A calcium phosphate precipitation is used to concentrate actinides from water samples. Tracers are used to monitor chemical recoveries and correct results to improve precision and accuracy.

3. SIGNIFICANCE OF USE

- 3.1. This is a rapid, reliable method for measurement of actinides in water samples that is more cost-effective and efficient than traditional ion exchange, solvent extraction and precipitation techniques.

4. INTERFERENCES

- 4.1. Actinides with unresolvable alpha energies such as ^{241}Am and ^{238}Pu , ^{237}Np and ^{234}U , or ^{232}U and ^{210}Po must be chemically separated to enable measurement. This method separates these isotopes effectively.

- 4.2. Very high levels of phosphate in the sample may lead to reduced actinide recoveries during calcium phosphate precipitation and column separations. Adjusting the amount of phosphate added to co-precipitate the actinides may be necessary in these cases.
- 4.3. The sample preparation procedure outlined in this method will adequately recover actinides from freshly collected, well preserved, homogenous water samples. Older, poorly preserved samples or samples with significant organic or solid matter may require more aggressive treatment to recover actinides which have precipitated or adsorbed to the walls of the storage container or solid matter. Rinsing the empty storage container with warm HNO₃, adjusting the HNO₃ concentration of the sample to 1M HNO₃ and boiling, and/or wet-ashing the calcium phosphate precipitate may be required for older, poorly preserved samples.

5. APPARATUS

- Analytical balance, 0.0001 g sensitivity
- Beakers, glass
- Centrifuge tubes, 50mL and 250mL
- Centrifuge, with rotor and carriers for 50mL and 250mL tubes
- Column rack - Eichrom Part: AC-103
- Extension Funnels, 25 mL, Eichrom Part: AC-120
- Fume hood
- Hotplate
- Stir rods, glass

6. REAGENTS

Note: Analytical grade or ACS grade reagents and are recommended. Evaluation of key reagents, such as aluminum nitrate and ammonium hydrogen phosphate, for contribution to method background levels from naturally occurring radioactive materials is recommended.

Aluminum nitrate nonahydrate, Al(NO ₃) ₃ ·9H ₂ O
Ammonium hydrogen phosphate, (NH ₄) ₂ HPO ₄
Ammonium hydroxide(57% NH ₄ OH or 28% NH ₃), concentrated NH ₄ OH
Ammonium oxalate monohydrate, (NH ₄) ₂ C ₂ O ₄ ·H ₂ O
Ammonium thiocyanate, NH ₄ SCN
Appropriate tracers or standards (Th-229, Am-243/Np-239)

Ascorbic acid powder, $C_6H_8O_6$
Calcium nitrate, $CaNO_3$
Deionized water, All reagents are prepared with deionized water
Ferric nitrate nonahydrate, $Fe(NO_3)_3 \cdot 9H_2O$
Hydrochloric acid (37%), concentrated HCl
Isopropanol, C_3H_7OH
Nitric acid (70%), concentrated HNO_3
Oxalic acid dihydrate, $H_2C_2O_4 \cdot 2H_2O$
Sulfamic acid, H_3NSO_3
TEVA [®] resin, 2mL prepacked column, 100-150 μ m, Eichrom Part TE-C50-A

- 6.1. Ammonium bioxalate (0.1M) - Dissolve 6.31g of oxalic acid dihydrate and 7.11g of ammonium oxalate monohydrate in 900mL of water. Dilute to 1L with water.
- 6.2. Ammonium hydrogen phosphate (3.2M) - Dissolve 104g of $(NH_4)_2HPO_4$ in 200mL of water. Heat gently to dissolve. Dilute to 250mL with water.
- 6.3. Ammonium thiocyanate indicator (1M) - Dissolve 7.6g of ammonium thiocyanate in 90mL of water. Dilute to 100mL with water.
- 6.4. Ascorbic acid (1M) - Dissolve 17.6g of ascorbic acid in 90mL of water. Dilute to 100mL with water. **Prepare fresh weekly.**
- 6.5. Calcium nitrate (1.25M) - Dissolve 51g of $Ca(NO_3)_2$ in 100mL of water. Dilute to 250mL with water.
- 6.6. Hydrochloric acid (9M) - Add 750mL of concentrated HCl to 100mL of water. Dilute to 1L with water.
- 6.7. Nitric acid solution (3M) - Add 188mL of concentrated HNO_3 to 800mL of water. Dilute to 1L with water.
- 6.8. Nitric acid solution (3M) - aluminum nitrate (1M) - Dissolve 375g $Al(NO_3)_3 \cdot 9H_2O$ in 500mL of water. Add 188mL of concentrated HNO_3 . Dilute to 1L with water.
- 6.9. Phenolphthalein indicator - Dissolve 1g of phenolphthalein in 50mL of isopropyl alcohol and add 50mL of water.
- 6.10. Ferric Nitrate Solution (5 mg/mL Fe)in 0.1M HNO_3 - To a 500mL volumetric flask, add 18g $Fe(NO_3)_3 \cdot 9H_2O$, 400mL of water and

- 3.1mL of concentrated HNO₃. Swirl to dissolve. Dilute to 500mL with water.
- 6.11. Sulfamic acid (1.5M) - In a 500mL volumetric flask, add 73g of sulfamic acid to 400mL of water. Swirl to dissolve. Dilute to 500mL with water.
- 6.12. Nitric acid (2.5M) - Sulfamic Acid (0.1M) - Ascorbic Acid (0.1M) - Fe Solution – To 250mL of water, add 78mL of concentrated HNO₃, 33mL 1.5M sulfamic acid and 15mL ferric nitrate solution (5mg/mL Fe). Swirl to mix. Add 50mL 1M ascorbic acid solution. Dilute to 500mL with water. **Prepare fresh daily.**

7. PROCEDURE

7.1. Water Sample Preparation:

- 7.1.1. If samples larger than 1L are analyzed, evaporate the sample to approximately 1L.
- 7.1.2. Aliquot 500 to 1000 mL of the filtered sample (or enough to meet required detection limit) into an appropriate size beaker.
- 7.1.3. Add 5mL concentrated HNO₃ acid.
- 7.1.4. Add appropriate tracers per lab protocol.
- 7.1.5. Calcium phosphate precipitation:
- 7.1.5.1. Add 0.5 mL of 1.25M Ca(NO₃)₂ to each sample.
- 7.1.5.2. Place each beaker on a hotplate. Cover each beaker with a watch glass. Heat samples at medium setting for 30-60 minutes.
- 7.1.5.3. Remove the watch glass and turn the heat down.
- 7.1.5.4. Add 2-3 drops of phenolphthalein indicator and 200 µL of 3.2M (NH₄)₂HPO₄ solution.
- 7.1.5.5. While stirring, slowly add enough concentrated NH₄OH to reach the phenolphthalein end point (pH~9) and form a calcium phosphate precipitate. Allow the sample to heat for another 20-30 minutes.
- 7.1.5.6. Remove samples from hot plate, cool to room temperature, and allow precipitate to settle until solution can be decanted (30 minutes to 2 hours) or centrifuge.

- 7.1.5.7. Decant supernate and discard to waste.
- 7.1.5.8. Transfer the precipitate to a centrifuge tube and centrifuge the precipitate for approximately 10 minutes at 2000 rpm.
- 7.1.5.9. Decant supernate and discard to waste.
- 7.1.5.10. Wash the precipitate with an amount of water approximately twice the volume of the precipitate. Mix well on a vortex mixer. Centrifuge for 5-10 minutes. Discard the supernate.
- 7.1.5.11. If an ammonia odor persists repeat 7.1.5.10.
- 7.1.5.12. Dissolve precipitate in approximately 5mL concentrated nitric acid. Transfer solution to a 100mL beaker. Rinse centrifuge tube with 2-3mL of concentrated nitric acid and transfer to beaker. Evaporate solution to dryness.

7.2. Np, Th separation from using TEVA Resin

- 7.2.1. Dissolve each precipitate with 10 mL of 3M HNO₃-1M Al(NO₃)₃.

Note: An additional 5-10mL may be necessary if the volume of precipitate is large.

Note: Make sure that all reagents and the load solution have cooled to room temperature. Warm solutions can cause reactions that will affect oxidation adjustments performed in the following steps.

- 7.2.2. Add 1mL of 1.5M sulfamic acid and 0.5mL of 5mg/mL Fe (as Fe(NO₃)₃) solution to each sample. Swirl to mix.

Note: If the additional 5-10mL was used to dissolve the sample in step 7.2.1.1, add a proportionately larger amount of sulfamic acid and iron solution.

- 7.2.3. Add 1 drop of 1M ammonium thiocyanate indicator to each sample and mix.

Note: The color of the solution turns deep red, due to the presence of Fe(III).

- 7.2.4. Add 1mL of 1M ascorbic acid to each solution, swirling to mix. Wait for 2-3 minutes.

Note: The red color should disappear which indicates reduction of Fe(III) to Fe(II). If the red color persists then additional ascorbic acid solution should be added drop-wise with mixing

until the red color disappears. Following disappearance of the red color, add an additional 0.5mL of 1M ascorbic acid.

Note: If particles are observed suspended in the solution, centrifuge the sample. The supernate will be transferred to the column in step 7.2.8. The residue will be discarded.

- 7.2.5. For each sample dissolved, place a TEVA Resin column in the column rack.
- 7.2.6. Place a beaker below each column, remove the bottom plug from each column, push top frit down to the top of the resin bed, and allow to drain. Attach column reservoirs to each column.
- 7.2.7. Add 5mL of 3M HNO₃ into each column to condition resin and allow to drain.
- 7.2.8. Transfer each dissolved sample into the appropriate TEVA Resin column reservoir. Allow solution to drain through columns.
- 7.2.9. Add 5mL of 2.5M HNO₃ - 0.1M Sulfamic Acid - 0.1M Ascorbic Acid- Fe Solution to rinse to each sample beaker. Transfer each rinse solution into the appropriate TEVA Resin column reservoir. Allow solution to drain through columns.
- 7.2.10. Add 10mL of 2.5M HNO₃ - 0.1M Sulfamic Acid - 0.1M Ascorbic Acid- Fe Solution to each column reservoir. Allow solution to pass through each column.
- 7.2.11. Add 20mL of 3M HNO₃ to each column reservoir. Discard eluate.

Note: Plutonium, uranium and americium are removed with the load solution and rinses.

- 7.2.12. Place a clean, labeled 50 mL beaker below each column.
- 7.2.13. Add 15mL of 9M HCl into each column reservoir to elute Th. Allow solution to drain through each column.
- 7.2.14. Set Th samples aside for alpha source preparation.
- 7.2.15. Place a clean, labeled 50mL beaker or centrifuge tube below each column.
- 7.2.16. Add 10 mL of 0.1M ammonium bioxalate into each column to strip the neptunium. Allow solution to drain through each column.

- 7.3. Prepare samples for the measurement of actinides by alpha spectrometry using electrodeposition (Eichrom SPA02) or rare earth fluoride micro precipitation (Eichrom SPA01).

8. CALCULATIONS

Calculate the actinide activity as follows:

Calculate tracer yield:

$$\text{Yield} = \frac{(C_s - B_s)}{E_s \times A_s}$$

where:

C_s = measured actinide tracer, cpm

B_s = background, cpm

E_s = counting efficiency for tracer

A_s = tracer activity, dpm

Note: *If any tracer may be present in the sample, a spiked and unspiked sample must be analyzed to determine chemical yield, where:*

$$Y = \frac{(\text{spiked sample tracer cpm} - \text{unspiked sample tracer cpm})}{E \times \text{actinide spike activity, dpm}}$$

Percent yield = Yield x 100

Calculate actinide isotope activity:

$$\text{Sample dpm/L} = \frac{S - B}{E \times V \times Y}$$

where:

S = sample activity, cpm

B = background, cpm

E = counting efficiency = measured cpm/dpm of isotopic standard

V = sample volume, L

Y = yield

Conversion of dpm/g to pCi/g: $\text{pCi/L} = (\text{dpm/L})/2.22$

9. PRECISION AND BIAS

- 9.1. *Precision - A relative standard deviation of 3.5% for Th in the range of 1 pCi/L to 20 pCi/L has been reported. A relative standard deviation of 6.3% for Np in the same range has been reported.*
- 9.2. *Bias - A mean recovery of 102% ± 3.5% for Th and 102% ± 6% for Np have been reported. Since results are corrected based on spike recovery, no significant bias exists for the method.*

10. REFERENCES

- 1) Horwitz, E P, et al. "Separation and Preconcentration of Actinides by Extraction Chromatography Using a Supported Liquid Anion Exchanger: Application to the Characterization of High-Level Nuclear Waste Solutions." *Analytica Chimica Acta*, 310, 63-78 (1995).
- 2) Maxwell, S.L., "Rapid Column Extraction Methods for Actinides and Sr-89/90 in Water Samples," *Journal of Radioanalytical and Nuclear Chemistry*, 267(3), 537-543 (2006).
- 3) ASTM Method C1475-05, "Standard Method for Determination of Neptunium-237 in Soil."
- 4) ASTM Method D7282-06, "Standard Practice for Set-up, Calibration, and Quality Control of Instruments Used for Radioactivity Measurements."
- 5) ASTM Method D3648-14, "Standard Practices for the Measurement of Radioactivity."

