

THORIUM IN WATER

(WITH VACUUM BOX SYSTEM)

1. SCOPE

- 1.1. This is a method for the separation of thorium from water samples. After completing this method, source preparation for measurement of thorium by alpha spectrometry is performed by electrolytic deposition onto stainless steel planchets (Eichrom Method SPA02) or by rare earth fluoride micro precipitation onto polypropylene filters (Eichrom Method SPA01).
- 1.2. This method does not address all aspects of safety, quality control, calibration or instrument set-up. However, enough detail is given for a trained radiochemist to achieve accurate and precise results for the analysis of the analyte(s) from the appropriate matrix, when incorporating the appropriate agency or laboratory safety, quality and laboratory control standards.

2. SUMMARY OF METHOD

- 2.1. Thorium is separated by Eichrom TEVA resin prior to measurement by alpha spectrometry. A calcium phosphate precipitation is used to concentrate actinides from water samples. Tracers are used to monitor chemical recoveries and correct results to improve precision and accuracy.

3. SIGNIFICANCE OF USE

- 3.1. This is a rapid, reliable method for measurement of actinides in water samples that is more cost-effective and efficient than traditional ion exchange, solvent extraction and precipitation techniques.

4. INTERFERENCES

- 4.1. Nuclides with unresolvable alpha energies such as ^{241}Am and ^{238}Pu , ^{237}Np and ^{234}U , or ^{232}U and ^{210}Po must be chemically separated to enable measurement. This method separates these isotopes effectively.
- 4.2. Very high levels of phosphate in the sample may lead to reduced recovery of actinides in calcium phosphate precipitation and column

separations. Adjusting the amount of phosphate added to co-precipitate the actinides may be necessary in these cases.

- 4.3. The sample preparation procedure outlined in this method will adequately recover actinides from freshly collected, well preserved, homogenous water samples. Older, poorly preserved samples or samples with significant organic or solid matter may require more aggressive treatment to recover actinides which have precipitated or adsorbed to the walls of the storage container or solid matter. Rinsing the empty storage container with warm HNO₃, adjusting the HNO₃ concentration of the sample to 1M HNO₃ and boiling, and/or wet-washing the calcium phosphate precipitate may be required for older, poorly preserved samples.

5. APPARATUS

- Analytical balance, 0.0001 g sensitivity
- Beakers, glass
- Cartridge reservoirs, 10mL (Eichrom Part: AR-200-RV10) or 20mL (Eichrom Part: AR-200-RV20)
- Centrifuge tubes, 50mL and 250mL
- Centrifuge, with rotor and carriers for 50mL and 250mL tubes
- Fume hood
- Hotplate
- Stir rods, glass
- Vacuum box system, Eichrom Part: AR-12-BOX or AR-24-BOX
- Vacuum Box white inner support tube-PE, Eichrom Part: AR-1000-TUBE-PE
- Vacuum box yellow outer tips, Eichrom Part: AR-1000-OT
- Vacuum pump - 115 V, 60 Hz Fisher Part: 01-092-25 (or equivalent) or house vacuum
- Optional for collection of load and rinse fractions:
 - Vacuum box Inner liner, Eichrom Part: AR-12-LINER or AR-24-LINER

6. REAGENTS

Note: Analytical grade or ACS grade reagents are recommended. Evaluation of key reagents, such as aluminum nitrate and ammonium hydrogen phosphate, for contribution to method background levels from naturally occurring radioactive materials is recommended.

Aluminum nitrate nonahydrate, $Al(NO_3)_3 \cdot 9H_2O$
Ammonium hydrogen phosphate, $(NH_4)_2HPO_4$

<i>Ammonium hydroxide(57% NH₄OH or 28% NH₃), concentrated NH₄OH</i>
<i>Ammonium oxalate monohydrate, (NH₄)₂C₂O₄·H₂O</i>
<i>Appropriate tracers or standards (Th-229)</i>
<i>Calcium nitrate, CaNO₃</i>
<i>Deionized water, All reagents are prepared with deionized water</i>
<i>Hydrochloric acid (37%), concentrated HCl</i>
<i>Isopropanol, C₃H₇OH</i>
<i>Nitric acid (70%), concentrated HNO₃</i>
<i>Oxalic acid dihydrate, H₂C₂O₄·2H₂O</i>
<i>Phenolphthalein pH Indicator</i>
<i>TEVA® resin, 2mL prepacked cartridge, 50-100µm, Eichrom Part TE-R50-S</i>

- 6.1. *Ammonium hydrogen phosphate (3.2M)* - Dissolve 104g of (NH₄)₂HPO₄ in 200mL of water. Heat gently to dissolve. Dilute to 250mL with water.
- 6.2. *Calcium nitrate (1.25M)* - Dissolve 51g of Ca(NO₃)₂ in 100mL of water and dilute to 250mL with water.
- 6.3. *Hydrochloric acid (9M)* - Add 750mL of concentrated HCl to 100mL of water. Dilute to 1L with water.
- 6.4. *Nitric acid solution (3M)* - Add 188mL of concentrated HNO₃ to 800mL of water. Dilute to 1L with water.
- 6.5. *Nitric acid solution (3M) - aluminum nitrate (1M)* - Add 188mL of concentrated HNO₃ to 500mL of water. Dissolve 375g Al(NO₃)₃·9H₂O. Dilute to 1L with water.
- 6.6. *Phenolphthalein indicator* - Dissolve 1g of phenolphthalein in 50mL of isopropyl alcohol and add 50mL of water.

7. PROCEDURE

7.1. Water Sample Preparation:

- 7.1.1. If samples larger than 1L are analyzed, evaporate the sample to approximately 1L.
- 7.1.2. Aliquot 500 to 1000mL of the filtered sample (or enough to meet required detection limit) into an appropriate size beaker.

- 7.1.3. Add 5mL concentrated HNO₃.
- 7.1.4. Add appropriate tracers and/or reference standards per lab protocol.
- 7.1.5. Calcium phosphate precipitation:
 - 7.1.5.1. Add 0.5mL of 1.25M Ca(NO₃)₂ to each sample.
 - 7.1.5.2. Place each beaker on a hotplate.
 - 7.1.5.3. Cover each beaker with a watch glass.
 - 7.1.5.4. Heat at medium setting for 30-60 minutes.
 - 7.1.5.5. Take the watch glass off the beaker and turn the heat down.
 - 7.1.5.6. Add 2-3 drops of phenolphthalein indicator and 1mL of 3.2M (NH₄)₂HPO₄ solution.
 - 7.1.5.7. While stirring, slowly add enough concentrated NH₄OH to reach the phenolphthalein end point and form a calcium phosphate precipitate. Heat sample for another 20-30 minutes.
 - 7.1.5.8. Remove samples from the hot plate, cool to room temperature, and allow precipitate to settle until solution can be decanted (30 minutes to 2 hours) or centrifuge.
 - 7.1.5.9. Decant supernate and discard to waste.
 - 7.1.5.10. Transfer the precipitate to a centrifuge tube and centrifuge the precipitate for approximately 10 minutes at 2000 rpm.
 - 7.1.5.11. Decant supernate and discard to waste.
 - 7.1.5.12. Wash the precipitate with an amount of water approximately twice the volume of the precipitate. Mix well on a vortex mixer. Centrifuge for 5-10 minutes. Discard the supernate.
 - 7.1.5.13. If an ammonia odor persists repeat 7.1.5.12.
 - 7.1.5.14. Dissolve precipitate in 5mL concentrated HNO₃. Transfer solution to a 100mL beaker. Rinse centrifuge tube with 2-3mL of concentrated HNO₃ and transfer to beaker. Evaporate solution to dryness.

7.2. Th Separation using TEVA Resin

7.2.1. Dissolve each precipitate with 10mL of 3M HNO₃-1M Al(NO₃)₃.

Note: An additional 5-10mL may be necessary if the volume of precipitate is large.

7.2.2. Place the inner tube rack (supplied with vacuum box system) into the vacuum box. Place centrifuge tubes in the rack. Fit the lid to the vacuum system box. Alternatively, a vacuum box inner liner may be used.

7.2.3. Place yellow outer tips into all 12 or 24 openings in the lid of the vacuum box. Fit a white inner support tube into each yellow tip.

7.2.4. For each sample solution, fit a TEVA cartridge on to the inner support tube.

7.2.5. Add syringe barrels (funnels/reservoirs) to the top end of each TEVA cartridge.

Note: The unused openings on the vacuum box should be sealed. Vacuum manifold plugs can be used to plug unused white tips to achieve good seal during the separation. Alternatively, unused vacuum box holes can be sealed with scotch tape.

7.2.6. Connect the vacuum pump to the box. Turn the vacuum pump on and ensure proper fitting of the lid.

7.2.7. Add 5mL of 3M HNO₃ into each TEVA cartridge reservoir to condition the resin. Adjust vacuum to achieve a flow rate of 1-2 mL/min. Allow solution to completely pass through each TEVA cartridge.

7.2.8. Transfer each sample solution into the appropriate TEVA cartridge reservoir. Allow solution to completely pass through each TEVA cartridge at 1-2mL/min. Th is retained on the resin.

7.2.9. Add 5mL of 3M HNO₃ to rinse to each sample tube. Transfer each rinse solution into the appropriate TEVA cartridge reservoir. Allow solution to completely pass through each TEVA cartridge at 1-2mL/min.

Note: Uranium, americium and Np(V) are removed with the load solution and 3M HNO₃ rinses.

7.2.10. Disengage vacuum and empty centrifuge tubes to waste or replace with fresh centrifuge tubes. Place clean reservoirs above each TEVA cartridge.

- 7.2.11. Add 30mL of 3M HNO₃ into each TEVA cartridge reservoir. Engage vacuum. Allow solution to completely pass through each TEVA cartridge at 1-2mL/min. Disengage vacuum. Discard eluate as waste.
- 7.2.12. Place a clean, labeled 50mL centrifuge tube below each TEVA cartridge. Replacing yellow outer tips and inner support tubes at this point can help ensure clean Th fractions in the following steps.
- 7.2.13. Add 15mL of 9M HCl into each TEVA cartridge to strip Th. Engage vacuum. Allow solution to completely pass through each TEVA cartridge at 1-2mL/min.
- 7.3. Prepare samples for the determination of thorium by alpha spectrometry using electrodeposition (Eichrom SPA02) or rare earth fluoride micro precipitation (Eichrom SPA01).

8. CALCULATIONS

Calculate the actinide activity as follows:

Calculate tracer yield:

$$\text{Yield} = \frac{(C_s - B_s)}{E_s \times A_s}$$

where:

- C_s = measured actinide tracer, cpm
 B_s = background, cpm
 E_s = counting efficiency for tracer
 A_s = tracer activity, dpm

Note: *If any tracer may be present in the sample, a spiked and unspiked sample must be analyzed to determine chemical yield, where:*

$$Y = \frac{(\text{spiked sample tracer cpm} - \text{unspiked sample tracer cpm})}{E \times \text{actinide spike activity, dpm}}$$

$$\text{Percent yield} = \text{Yield} \times 100$$

Calculate actinide isotope activity:

$$\text{Sample dpm/L} = \frac{S - B}{E \times V \times Y}$$

where:

- S = sample activity, cpm
B = background, cpm
E = counting efficiency = measured cpm/dpm of isotopic standard
V = sample volume, L
Y = yield

Conversion of dpm/g to pCi/g:

$$\text{pCi/L} = (\text{dpm/L})/2.22$$

9. PRECISION AND BIAS

- 9.1. *Precision* - A relative standard deviation of 3.5% for Th in the range of 1 pCi/L to 20 pCi/L has been reported.
- 9.2. *Bias* - A mean recovery of 102% ± 3.5% for Th have been reported. Since results are corrected based on spike recovery, no significant bias exists for the method.

10. REFERENCES

- 1) Horwitz, E.P., et al. "Separation and Preconcentration of Uranium from Acidic Media by Extraction Chromatography." *Analytica Chimica Acta*, 266, 25-37, (1992).
- 2) Horwitz, E.P., et al. "Separation and Preconcentration of Actinides from Acidic Media by Extraction Chromatography." *Analytica Chimica Acta*, 281, 361-372, (1993).
- 3) Maxwell, S.L., "Rapid Column Extraction Methods for Actinides and Sr-89/90 in Water Samples," *Journal of Radioanalytical and Nuclear Chemistry*, 267(3), 537-543 (2006).
- 4) ASTM Method D7282-06, "Standard Practice for Set-up, Calibration, and Quality Control of Instruments Used for Radioactivity Measurements."
- 5) ASTM Method D3648-14, "Standard Practices for the Measurement of Radioactivity."

