

Rapid Determination of Sr in Emergency Milk Samples

Summary of Method Strontium is separated and concentrated from 100mL milk samples using a calcium phosphate precipitation. The precipitate is dissolved with nitric acid and centrifuged to remove residual protein and fat. The supernate, containing Sr, is wet ashed with $\text{HNO}_3\text{-H}_2\text{O}_2$ and then heated in a muffle furnace at 550°C for 30-60 minutes to destroy any residual organic matter. The muffled residue is wet ashed again with $\text{HNO}_3\text{-H}_2\text{O}_2$ and dissolved in $\text{HNO}_3\text{-Al}(\text{NO}_3)_3$. Strontium is separated from matrix impurities and potentially interfering radionuclides in the sample using stacked 2mL and 1mL cartridges of Eichrom Sr Resin. Radiostrontium is measured on a low background gas flow proportional counter or liquid scintillation counter. Chemical yield of strontium is determined by gravimetric recovery of stable strontium or by ICP-AES measurement. Average chemical recovery of strontium is $75 \pm 17\%$. Measured values of ^{90}Sr agreed to within 3.2% and 0.5% of reference values for 20 minute count times and 60 minute count times, respectively. The lower limit of detection for 100mL samples with 20 minute count times is 0.5Bq/L and with 60 minute count times is 0.16Bq/L. A single operator can prepare batches of 12-24 samples for ^{90}Sr measurement in less than 8 hours.

Reagents

Sr Resin, 2mL Cartridges (Eichrom SR-R50-S)
 Sr Resin, 1mL Cartridges (Eichrom SR1ML-R50-S)
 Ammonium Hydroxide (listed as 28% NH_3 or 56% NH_4OH)
 Nitric Acid (70%)
 Hydrogen Peroxide (30%)
 Deionized Water
 1.25M $\text{Ca}(\text{NO}_3)_2$
 3.2M $(\text{NH}_4)_2\text{HPO}_4$
 Strontium Carrier (10mg/mL)
 2M $\text{Al}(\text{NO}_3)_3$
 ^{90}Sr standard
 Oxalic acid

Equipment

Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)
 Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)
 Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)
 Yellow Outer Tips (Eichrom AR-1000-OT)
 50mL and 250mL Centrifuge Tubes
 Centrifuge
 Cupped Stainless Steel Planchets (~5mL volume)
 Gas Flow Proportional Counter
 Muffle Furnace
 Hot Plate
 Analytical Balance
 250mL Glass Beakers
 Vacuum Pump

Figure 1. Sample Preparation





Figure 2. Load Solution Preparation and Strontium Separation



Performance of Radiostrontium in Milk Method

20 Minute Count Times			
⁹⁰ Sr, reference (Bq/L)	⁹⁰ Sr, measured (Bq/L)	Uncertainty %, k = 2	% Bias
0	0.26	98.9	N/A
0	0.26	81.9	N/A
2.86	2.66	24.1	-7.0
2.86	3.96	24.7	38
2.86	3.31	20.2	15.7
2.86	2.67	18.7	-6.6
5.7	6.11	16.7	7.2
5.7	5.71	13.1	0.2
5.7	5.16	13.9	-9.5
14.3	12.8	9.1	-11
14.3	15.2	8.5	6.3
14.3	14.1	8.6	-1.4

60 Minute Count Times			
⁹⁰ Sr, reference (Bq/L)	⁹⁰ Sr, measured (Bq/L)	Uncertainty %, k = 2	% Bias
0	0.11	130	N/A
0	0.27	59	N/A
2.86	3.09	13.2	8.0
2.86	3.11	16.7	8.7
2.86	2.67	13.6	-6.6
2.86	2.67	11.3	-6.6
5.7	5.85	10.4	2.6
5.7	5.75	8.3	0.9
5.7	6.04	8.2	5.9
14.3	13.6	6.1	-4.9
14.3	14.0	6.1	-2.1
14.3	14.2	6.1	-0.7

References

1) Sherrod L. Maxwell, Brian K. Culligan, "Rapid method for the determination of radiostrontium in emergency milk samples," *J. Radioanal. Nucl. Chem.*, 279(3), 757-760 (2009).