

AN-1422-10

# Rapid Determination of <sup>226</sup>Ra in 5g Vegetation Samples

**Summary of Method** <sup>226</sup>Ra is separated from 5 gram samples of vegetation and measured by alpha spectrometry. Samples are fused with sodium hydroxide at 600°C. The fusion cake is dissolved in water, and radium is precipitated from samples with calcium carbonate. The calcium carbonate precipitate is dissolved in hydrochloric acid, and cation exchange chromatography is used to purify radium and barium from matrix ions. Barium is removed from samples using Eichrom Sr Resin. Eichrom DGA Resin is used to separate other alpha emitting nuclides from radium. Samples are prepared for alpha spectrometry by barium sulfate micro-precipitation onto Eichrom® Resolve Filters. Sample preparation, including alpha spectrometry source preparation, for batches of 12 samples can be completed by a single operator in as little as 6 hours. Yields are traced with <sup>225</sup>Ra(<sup>229</sup>Th) by alpha spectrometry. At least 8 hours of ingrowth time for the alpha emitting <sup>217</sup>At daughter of <sup>225</sup>Ra is required prior to measurement by alpha spectrometry.

#### Reagents

Cation Exchange Resin (Eichrom C8-B500-F-H) Sr Resin, 2mL Cartridges (Eichrom SR-R50-S)

DGA Resin, Normal 2mL Cartridges (Eichrom DN-R50-S)

Nitric Acid (70%) Hydrochloric Acid (37%)

<sup>225</sup>Ra(<sup>229</sup>Th) Tracer 1.25M Ca(NO<sub>3</sub>)<sub>2</sub>

2M Na<sub>2</sub>CO<sub>3</sub> Barium Carrier (1mg/mL) Isopropyl Alcohol Ammonium Sulfate

Denatured Ethanol Ascorbic Acid

Sodium Hydroxide Hydrogen Peroxide (30%)

#### **Equipment**

Plastic Chromatography Column (Eichrom AC-50E-5M)

Column Extension Funnel (Eichrom AC-20X-20M)

Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)

Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)

Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)

Yellow Outer Tips (Eichrom AR-1000-OT)

Resolve Filter in Disposable Funnel (Eichrom RF-DF-25-25PP01)

50mL and 250mL Centrifuge Tubes

Centrifuge

Stainless Steel Planchets with adhesive tape

Hotplate

Alpha Spectrometry System

150mL Glass beakers

Vacuum Pump

250mL Zirconium Crucible w/ lid

Muffle Furnace

**Heat Lamp** 

#### Figure 1. Sample Preparation

5g Vegetation +Tracer <sup>225</sup>Ra(<sup>229</sup>Th) in Zr crucible

Heat at 600°C in muffle furnace for 2 hours.

Wet ash on hotplate with 5mL 70% HNO<sub>3</sub> and 30% H<sub>2</sub>O<sub>2</sub>.

Heat at 600°C in muffle furnace for 5-10 minutes.

Add 10g NaOH.

Fuse at 600°C in muffle furnace for 15 minutes.

Remove from furnace. Cool 10 minutes.

Dissolve fusion cake with 100mL DI water.

Transfer to 250mL centrifuge tube. Add 10mL 37% HCl. Dilute to 150mL.

Add 0.5mL 1.25M  $Ca(NO_3)_2$  and

10mL 2M Na<sub>2</sub>CO<sub>3</sub>. Mix well.

Place in ice bath for 10 minutes.


Centrifuge 3500 rpm, 10 min

Dissolve precipitate in 20mL 1.5M HCl and 3mL 1.5M ascorbic acid.

Proceed to Column Purification

Decant Supernate To Waste

### Figure 2. Column Purification and Alpha Source Preparation



<sup>1</sup>If using <sup>133</sup>Ba tracer, 3.0g of cation exchange resin and proportionally smaller rinse volumes may be used. <sup>2</sup>If tracing with <sup>229</sup>Th, a 20mL 1M HCI-1M H<sub>3</sub>PO<sub>4</sub> rinse following the sample load can improve purity of final <sup>226</sup>Ra fraction.

## Method Performance <sup>226</sup>Ra in 5g Vegetation Samples

	<sup>225</sup> Ra( <sup>217</sup> At)	<sup>226</sup> Ra(mBq/g)	<sup>226</sup> Ra(mBq/g)	
Sample**	% Yield*	Reference	Measured	%Bias
1	91.5	73.8	70.8	-4.1
2	88.3	73.8	73.8	0.0
3	93.1	73.8	69.8	-5.4
4	82.2	73.8	68.5	-7.2
5	80.2	73.8	81.4	10.3
AVG	87 <u>+</u> 6	73.8	73 <u>+</u> 5	-1.1

<sup>\*&</sup>lt;sup>225</sup>Ra tracer is added as <sup>229</sup>Th in equilibrium with its daughters and measured by its alpha emitting <sup>217</sup>At daughter (7.066MeV) after >8 hours ingrowth.

#### References

1) Sherrod L. Maxwell, Brian K. Culligan, "Rapid Determination of <sup>226</sup>Ra in Environmental Samples," *J. Radioanal. Nucl. Chem.*, 293(1), 149-155 (2012).

<sup>\*\*5</sup> grams of blank hay matrix spiked with <sup>226</sup>Ra