

Rapid Determination of Plutonium in Large Rice Samples

Summary of Method Plutonium is separated and measured from up to 1.5kg rice samples. Rice samples are muffled and wet ashed to reduce volume and destroy organic content. The residue is then fused with sodium hydroxide. Precipitation steps remove additional matrix and prepare plutonium for separation on Eichrom TEVA resin. Plutonium is measured by alpha spectrometry following cerium fluoride microprecipitation onto Eichrom Resolve® Filters. Plutonium recovery through the method, determined using ²⁴²Pu tracer, was $87 \pm 4\%$ for 1kg samples. Measured values for ²³⁹Pu and ²³⁸Pu agreed within 6% of reference values, even when refractory ²³⁹Pu was present in the sample. Sample preparation can be completed in less than 48 hours.

Reagents

TEVA Resin, 2mL Cartridges (Eichrom TE-R50-S)
 Iron carrier (50mg/mL Fe, as ferric iron nitrate)
 La carrier (10mg/mL)
 Ce carrier (1mg/mL)
 Deionized Water 1.25M Ca(NO₃)₂
 3.2M (NH₄)₂HPO₄ 2M Al(NO₃)₃
 10% (w:w) TiCl₃ HNO₃ (70%)
 HCl (37%) NaOH
 HF (49%) or NaF Boric acid
 H₂O₂ (30%) NaNO₂
 Denatured ethanol Sulfamic Acid
 Ascorbic Acid ²⁴²Pu tracer

Equipment

Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)
 Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)
 Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)
 Yellow Outer Tips (Eichrom AR-1000-OT)
 Resolve Filters in Funnel (Eichrom RF-DF25-25PP01)
 50mL and 250mL Centrifuge Tubes
 250mL Zirconium crucibles with zirconium lids
 Stainless Steel Planchets with adhesive tape
 Alpha Spectrometry System
 Centrifuge Heat Lamp
 Muffle Furnace Hot Plate
 Analytical Balance Vacuum Pump
 600mL glass beakers

Figure 1. Sample Preparation







Figure 2. Plutonium Separation on TEVA Resin and Source Preparation

<p>(1) Precondition 2mL TEVA, 5mL 3M HNO₃.</p> <p>(2) Load sample solution.</p> <p>(3) Rinse sample tube with 5mL 3M HNO₃. * Add tube rinse to cartridges.</p> <p>(4) Rinse TEVA cartridge with: -15mL 3M HNO₃ -20mL 9M HCl (Th removal) -5mL 3M HNO₃</p> <p>(5) Strip Pu from TEVA cartridge with 20mL 0.1M HCl-0.05MHF-0.01M TiCl₃.</p> <p>-If measuring Pu by ICP-MS, Pu may be stripped from TEVA with 20mL of 0.05M HCl-0.025M HF-0.02M hydroxylamine-HCl.</p> <p>-If preparing Pu sources for alpha spectrometry by electrodeposition, strip Pu with 20mL 0.1M HCl-0.025M HF-0.02M rongalite (sodium-hydroxymethanesulfinate).</p> <p>(6) Add 0.5mL 30% H₂O₂ for Uranium decontamination in rare earth fluoride precipitation alpha source</p>	<p>(7) Add 50ug Ce carrier to all samples. Mix well.</p> <p>(8) Add 1mL 49% HF. Mix well. Wait 15-20 minutes.</p> <p>(9) Set up Resolve® Filter Funnel on vacuum box.</p> <p>(10) Wet filter with 3mL 80% ethanol followed by 3mL DI water.</p> <p>(11) Filter sample.</p> <p>(12) Rinse sample tube with 5mL DI water and add to filter.</p> <p>(13) Rinse filter funnel with 3mL DI water and 2mL 100% ethanol.</p> <p>(14) Draw vacuum until filter is dry.</p>	<p>(15) Remove filter from funnel assembly and mount filter on stainless steel planchet with 2-sided tape.</p> <p>(16) Dry filter under heat lamp for 3-5 minutes.</p> <p>(17) Measure plutonium by alpha spectrometry.</p>
		

*Adding 50uL of 30% H₂O₂ to tube rinse can improve Uranium decontamination.

Method Performance

Sample (kg)	Replicates	²⁴² Pu Tracer % Yield	Reference (mBq/kg)		Measured (mBq/kg)		% Bias	
			²³⁹ Pu	²³⁸ Pu	²³⁹ Pu	²³⁸ Pu	²³⁹ Pu	²³⁸ Pu
1.0	8	87 ± 4	12.5	10.6	11.8 ± 1.0	10.5 ± 0.7	-5.6	-0.7

MDA for 1 kg sample, 30hours count time, 0.37uBq/kg

References

1) Sherrod L. Maxwell, Brian K. Culligan, Jay B. Hutchinson, "Rapid fusion method for determination of plutonium isotopes in large rice samples," *J. Radioanal. Nucl. Chem.*, 298(2), 1367-1374 (2013).