

Rapid Determination of Actinides in Urine by ICP-MS + Alpha Spec.

AN-1437-10

Summary of Method Actinides are separated and concentrated from 100mL urine samples. Actinides are concentrated from urine samples using a calcium phosphate precipitation. Pu, Np, Am-Cm, and U are separated on 2mL cartridges of Eichrom TEVA, TRU and DGA resins. Pu-Np are measured by ICP-MS. Measured values for ²³⁹Pu and ²³⁷Np agreed to within 1-2% of reference values, while ²⁴¹Am and ²⁴⁴Cm agreed to within 2-3%. Decontamination factors of >106 were achieved for Pu over U (²³⁸U-H can interfere with the measurement of ²³⁹Pu by ICP-MS). Sample preparation for batches of 12 samples can be completed by a single operator in <8 hours.

Reagents

TEVA Resin, 2mL Cartridges (Eichrom TE-R50-S)

TRU Resin, 2mL Cartridges (Eichrom TR-R50-S)

DGA Resin, Normal, 2mL Cartridges (Eichrom DN-R50-S)

Iron carrier (50mg/mL Fe, as ferric iron nitrate)

 $^{\rm 242}Pu$ (ICP-MS) or $^{\rm 236}Pu$ (alpha) tracer

²³³U (ICP-MS) or U²³² (alpha) tracer

 243 Am tracer Ce carrier (1mg/mL) 1.25M Ca(NO₃)₂ 3.2M (NH₄)₂HPO₄

Deionized Water 2M Al(NO₃)₃ HNO₃ (70%) HCI (37%)

NH₄OH HF (49%) or NaF NaNO₂ Denatured ethanol Sulfamic Acid Ascorbic Acid

Oxalic acid/Ammonium oxalate
Hydroxylamine Hydrochloride

Equipment

Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)

Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)

Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)

Yellow Outer Tips (Eichrom AR-1000-OT)

Resolve Filters in Funnel (Eichrom RF-DF25-25PP01)

600mL Glass beakers

Stainless Steel Planchets with adhesive tape

Alpha Spectrometry System

ICP-MS System

50mL and 250mL Centrifuge Tubes

Centrifuge

Heat Lamp

Hot Plate

Analytical Balance

Vacuum Pump

Figure 1. Sample Preparation

Urine Sample + Tracers

Add 1mL 1.25M Ca(NO₃)₂ and 3mL 3.2M (NH4)₂HPO₄. Mix well.

Heating to near boiling can improve tracer/analyte equilibration.

Adjust to pH 9 with NH₄OH. Cool to room temperature.

Centrifuge 3500 rpm, 5 min.

Decant Supernate.

Dissolve precipitate in 8mL 6M HNO₃ and 8mL 2M Al(NO₃)₃.

Cool samples to room temp.

Fix valence by adding: (mix between steps)


-0.5mL 1.5M sulfamic acid

-40uL 50mg/mL Fe carrier

-1.5mL 1M ascorbic acid (Wait 3 min)

-1mL 3.5M NaNO₂

Figure 2. Actinide Separation on TEVA - DGA



^{*} Adding 50uL of 30% H₂O₂ to the 6M HNO₃ tube rinse can help improve U recoveries and decontamination in the Np and Pu fractions.

References

- 1) Sherrod L. Maxwell, Vernon D. Jones, "Rapid determination of Actinides in urine by ICP-MS and alpha spectrometry: A hybrid approach," *Talanta*, 80(1), 143-150 (2009).
- 2) Sherrod L. Maxwell, Brian K. Culligan, Vernon D. Jones, Sheldon T. Nichols, Gary W. Noyes, Maureen A. Bernard, "Rapid Determination of ²³⁷Np and Plutonium Isotopes by ICP-MS and Alpha Spectrometry," *Health Physics*, 101(2), 180-186 (2011).

^{**} Adding a 1mL UTEVA cartridge between TEVA and DGA can help improve uranium decontamination.