

## 90Y Generator

AN-1615-10

**Summary of Method** A method for the preparation of  $^{90}$ Y ( $t_{1/2} = 64.1$  hours) from  $^{90}$ Sr ( $t_{1/2} = 28.6$  years) source material is presented. The method employs 2mL cartridges of Sr and DGA resins to obtain high purity  $^{90}$ Y in small volumes of eluate while preserving valuable  $^{90}$ Sr source material. The source material, containing  $^{90}$ Sr/ $^{90}$ Y, in 4M HNO<sub>3</sub>, is loaded onto stacked 2mL cartridges of Sr and DGA resins.  $^{90}$ Sr is retained on Sr Resin, while  $^{90}$ Y is retained on DGA. The  $^{90}$ Sr source is recovered from Sr Resin with a small volume of 0.1M HCI. Following a suitable ingrowth period, the  $^{90}$ Sr can be acidified to 4M HNO<sub>3</sub> and used to produce additional  $^{90}$ Y. The  $^{90}$ Sr is preserved nearly indefinitely and continuously purified from chemical and radiologic impurities run to run.  $^{90}$ Y is recovered from DGA resin with 0.1M HCI. For applications where  $^{90}$ Y must be recovered in minimal volumes, DGA, Branched may be used in place of DGA, Normal.

#### Reagents

Sr Resin Cartridges (Eichrom SR-R50-S)

DGA, Normal Resin Cartridges (Eichrom DN-R50-S) or DGA, Branched Resin Cartridges (Eichrom DB-R50-S) Liquid Scintillation Cocktail

90Sr Source

**Deionized Water** 

HCI

HNO<sub>3</sub>


### **Equipment**


Glass vials for storage of 90Sr source.

Glass or plastic vials/bottles for collection of  $^{90}\mathrm{Y}$  and waste.

5, 10 or 20mL plastic luer lock syringes Liquid Scintillation system for measurement of <sup>90</sup>Sr and <sup>90</sup>Y.\*

\*90Y may also be measured by Cerenkov counting without the addition of scintillation cocktail.





# 90Sr/90 Separation

Sr

D

G

- (1) Precondition stacked 2mL cartridges of Sr and DGA Resins with 10mL 4M HNO<sub>3</sub>.
- (2) Acidify <sup>90</sup>Sr eluate from previous separation with 5mL conc. HNO<sub>3</sub>. (If new <sup>90</sup>Sr source, dilute to 20mL with 4M HNO<sub>3</sub>.)\*
- (3) Load <sup>90</sup>Sr and <sup>90</sup>Y in 20mL 4M HNO<sub>3</sub>.
- (4) Rinse Sr/DGA with 5mL 4M HNO<sub>3</sub>.
- (5) Separate Sr and DGA cartridges.


(6) Rinse DGA with 10mL 8M HCI.


D

G

Sr


- (7) Strip <sup>90</sup>Y with 10mL 0.1M HCI. (For higher <sup>90</sup>Y recovery in smaller volume, strip in opposite direction or use DGA, Branched)
- (8) Strip 90Sr from Sr Resin cartridge with 15mL 0.5M HCI. (More complete recovery may be obtained by stripping in opposite direction of load.) Save 90Sr for future use.





[HCI], M

\*Adding 1mg of stable Sr to the 90Sr source can help improve 90Sr recovery from Sr Resin (do only once, not each time).



#### References

1) McAlister and Horwitz, "Chromatographic Generator Systems for the actinides and natural decay series elements," *Radiochimica Acta*, 99:1-9 (2011).