

## Metal Impurities in Uranium, Plutonium and Mixed Oxides

**Introduction:** Metal impurities in uranium, plutonium and mixed oxides are measured in nuclear power plant fuel, weapons and waste. ICP-AES, ICP-MS, TIMS and alpha spectroscopy are used to measure a wide variety of impurities. New methods have been developed to shorten the analysis time, decrease waste generation and improve precision of these analyses.

**The Approach:** To avoid the need to maintain gram per liter Pu and U standards for matrix matching, these new methods remove the uranium and/or plutonium from the sample prior to instrumental analysis. The approach is quite versatile, allowing U only to be removed, while Pu and other metals pass, or capturing both U and Pu together while allowing the other matrix impurities to pass and be analyzed.

For example, the specification for neptunium and plutonium in nuclear fuel is less than, 3.3 Bq/g U. Removal of the uranium from the sample matrix makes the use of alpha spectroscopy possible by eliminating uranium self adsorption and the source of an interfering spectra from  $^{234}\text{U}$  which overlaps with  $^{237}\text{Np}$  at 4.7 MeV.

Uranium and/or Plutonium are removed from the sample by the use of an extraction chromatographic resin, marketed under the name of UTEVA<sup>®</sup> Resin. This resin is comprised of diamyl amyolphosphonate sorbed on an inert polymeric support. A nitric acid uptake curve for UTEVA Resin is shown in Figure 1. Virtually no other metals sorb on UTEVA when loaded from nitric acid.




Figure 1. UTEVA Resin Acid Dependency

### References:

- C. Good, Third International Uranium Hexafluoride Conference, Paducah KY, (1995)
- S. L. Maxwell, III, Eichrom Western Users' Group Workshop, Albuquerque, MN (2000)
- M. J. Brisson, ACS Spring National Meeting, Orlando, FL (2002)

### Case History 1: Analysis of impurities in Pu/U.

A method has been developed and implemented at the Washington Savannah River site to characterize legacy U and Pu solids stored for over a quarter of a century. Analysis was required as part of a stabilization and repackaging program associated with long term storage. Over 1600 cans of material were analyzed containing various amounts of U and Pu-239. Figure 2 is a schematic of the U/Pu removal using UTEVA Resin. Uranium and plutonium samples were spiked with a variety of metals. Their ICP-MS recoveries are presented in Table 1.




Figure 2. U/Pu matrix removal using UTEVA Resin for trace metal analysis.

### Case History 2: Analysis of transuranic elements in UF<sub>6</sub>

A variation on the method in case 1 allows for the collection and measurement of plutonium with the trace metal fraction. This is accomplished by dissolving the uranium sample in both 2 M nitric acid and 0.1 M oxalic acid. After a 15 mL 0.1 M oxalic acid/ 2 M nitric acid rinse, the combined load and rinse are passed through one additional UTEVA column. This accomplishes the  $>10^8$  decontamination factor for uranium based on a 0.5 gram sample of 4.4% enriched U.

Table 1.

Element	% Recovery
Ag	106
Al	101
As	88
B	89
Ba	106
Be	90
Ca	NA
Cd	94
Ce	108
Cr	103
Cu	106
Fe	106
Ga	101
Hf	90
HG	77
K	102
La	108
Li	101
Mg	103
Mo	101
Na	98
Nb	98
Ni	103
P	154
Pb	99
Se	87
Si	132
Ta	84
V	104
W	113
Zn	91
Zr	63

# eichrom

Eichrom Technologies, Inc.

8205 S. Cass Ave.

Suite 106

Darien, Illinois

60561 USA

+1-630-963-0320 (tel)

[info@eichrom.com](mailto:info@eichrom.com)

[www.eichrom.com](http://www.eichrom.com)

Eichrom Europe

Campus de Ker Lann

Parc de Lormandière – Bât. C

rue Maryse Bastié

35170 Bruz – FRANCE

+33 (0)2 23 50 13 80 (tel)

[EichromEurope@eichrom.com](mailto:EichromEurope@eichrom.com)