

RADIUM-228 IN WATER

1. SCOPE

- 1.1. This is a method for separation and measurement of ^{228}Ra in water via its beta emitting ^{228}Ac daughter.
- 1.2. This method does not address all aspects of safety, quality control, calibration or instrument set-up. However, enough detail is given for a trained radiochemist to achieve accurate and precise results for the analysis of the analyte(s) from the appropriate matrix, when incorporating the appropriate agency or laboratory safety, quality and laboratory control standards.

2. SUMMARY OF METHOD

- 2.1. A barium sulfate precipitation technique is used to concentrate radium from water samples. ^{133}Ba is used to monitor chemical recovery and correct results to improve precision and accuracy. After ingrowth, ^{228}Ac , a daughter of ^{228}Ra , is separated by Ln Resin, prior to measurement by low background gas flow proportional counter.

3. SIGNIFICANCE OF USE

- 3.1. This is a rapid, reliable method for measurement of ^{228}Ra in water samples that is based on chemistry similar to EPA procedure RA-05.

4. INTERFERENCES

- 4.1. Potential beta emitters such as bismuth, yttrium and thorium would be retained on the Ln Resin, while eluting actinium from the resin.
- 4.2. Interferences from other radioactive rare earth elements are eliminated under the stripping conditions of 0.35M HNO_3 .

5. APPARATUS

- Analytical balance- 0.0001 g sensitivity
- Beakers, glass
- Bunsen burner
- Centrifuge, with rotor and carries for 15mL tubes

- Centrifuge tubes, 50mL
- Column rack, Eichrom Part: AC-103
- Extension funnels, 25mL, Eichrom Part: AC-120
- Fume hood
- Gamma spectrometry system (for determination of ^{133}Ba recovery)
- Hotplate
- Liquid scintillation vials
- Low background gas flow proportional counter with appropriate sample carriers and planchets
- Stir rods, glass
- Test tubes, Pyrex glass, 15mL
- Vortex mixer
- Watch glasses

6. REAGENTS

Note: Analytical grade or ACS grade reagents are recommended.

^{133}Ba tracer (~3000 dpm/mL)
Barium chloride dihydrate, $\text{BaCl}_2 \cdot 2\text{H}_2\text{O}$
Deionized water, all reagents are prepared with deionized water
Hydrogen peroxide (30%), concentrated H_2O_2
LN [®] resin- 2mL prepacked column, 100-150 μm , Eichrom Part LN-C50-A
Nitric acid (70%), concentrated HNO_3
Potassium carbonate, K_2CO_3
Sulfuric acid (96%), concentrated H_2SO_4

- 6.1. *Barium carrier (30 mg/ml)*- Dissolve 13.3g $\text{BaCl}_2 \cdot 2\text{H}_2\text{O}$ in 200mL deionized water. Dilute 250mL with water.
- 6.2. *Nitric acid solution (0.095 M)*- Add 5.9mL of concentrated HNO_3 to 900mL of water. Dilute to 1L with water. Check that pH of solution is ~1.
- 6.3. *Nitric acid solution (0.35 M)*- Add 21.9mL of concentrated HNO_3 to 800mL of water. Dilute to 1L with water.
- 6.4. *Potassium carbonate (50wt%)*- Dissolve 50g K_2CO_3 in 50mL of water.

7. PROCEDURE

7.1. Water Sample Preparation:

- 7.1.1. If required, filter the sample through a 0.45 micron filter.
- 7.1.2. Aliquot 500 to 1000mL of the sample (or enough to meet required detection limit) into an appropriate size beaker. Add Ba-133 tracer and 1mL of barium chloride carrier.
- 7.1.3. Add 10mL of concentrated H₂SO₄ to sample. Place each beaker on a hot plate.
- 7.1.4. Cover each beaker with a watch glass. Heat samples at about 90°C for 2 hours.
- 7.1.5. Allow samples to cool. A fine white precipitate is formed. Allow precipitate to settle until solution can be decanted or centrifuge.
- 7.1.6. Decant supernate and discard to waste.
- 7.1.7. Transfer the precipitate to a 15mL Pyrex glass centrifuge tube using water. Centrifuge for 10 minutes at 2000 rpm.
- 7.1.8. Decant supernate and discard to waste.
- 7.1.9. Wash the precipitate with 4-5 mL of water. Mix well and centrifuge for 5-10 minutes. Check the pH of the supernate. The pH should be about 6, if necessary decant supernate and repeat wash/centrifuge steps.

7.2. Conversion of Ba(Ra)SO₄ to Ba(Ra)CO₃:

- 7.2.1. Add 1 mL of K₂CO₃ solution and 2-3 mL of water, submerging the precipitate completely.
- 7.2.2. Heat the bottom of the Pyrex tube over a low flame while stirring the tube periodically. Continue heating, mixing until volume is reduced to ~1mL. (This step may also be performed using microwave heating or a water bath.)
- 7.2.3. Let the tube cool and repeats step 7.2.1 and 7.2.2 two times.
- 7.2.4. Add 15-20 mL of water. Mix well. Centrifuge and check the pH of the supernate. The pH of at this point should be approximately 12.
- 7.2.5. Discard the supernate and continue washing with water and centrifuging until a pH of 7 is achieved. Discard supernate.

Note: The precipitate now should be primarily BaCO_3 . The water washing of the BaCO_3 precipitate removes excess CO_3^{2-} and SO_4^{2-} ions. Any sulfate remaining in the sample can result in incomplete dissolution during the following steps and decreased Ba/Ra yields.

- 7.2.6. Once the pH of 7 is achieved add 5 mL of 0.095M HNO_3 to the tube and gently heat as necessary to dissolve the residue in the tube.

Note: If the entire residue does not dissolve in the load solution, then centrifuge the solution and check the supernatant for ^{133}Ba yield recovery. If acceptable (>80%) recovery is produced then proceed further if not then take the residue (undissolved) and repeat steps 7.2.1 through 7.2.5. Combine the dissolved solutions and check for ^{133}Ba recovery.

- 7.2.7. Transfer the dissolved solution into a 20mL liquid scintillation vial. Rinse the Pyrex tube with an additional 5 mL of 0.095M HNO_3 and add the tube rinse to the same vial.
- 7.2.8. Count the ^{133}Ba by gamma spectrometry. Record ^{133}Ba recovery.
- 7.2.9. Let the solution sit in the vial for at least 30 hours for ingrowth of ^{228}Ac from ^{228}Ra .

7.3. ^{228}Ac Separation Using Ln Resin:

Note: ^{228}Ac has a 6.13 hour half-life. The following steps should be performed quickly to achieve the lowest detection limits.

- 7.3.1. Place an LN Resin column in the column rack for each sample.
- 7.3.2. Remove the cap and bottom plug from each column, push the top frit down to the top of the resin bed, and allow each column to drain.
- 7.3.3. Add 5mL of 0.095M HNO_3 into each column to precondition resin. Allow each column to drain.
- 7.3.4. Transfer each solution from step 7.2.9. into the appropriate LN Resin column reservoir.
- 7.3.5. Allow the load solution to drain through column.
- 7.3.6. Add 5mL of 0.095M HNO_3 into the sample vial and transfer this rinse to the appropriate column reservoir. Allow each column to drain. Record the date and time at the start of this rinse for decay correction.

- 7.3.7. Rinse each column with 5mL of 0.095M HNO₃. Allow each column to drain. Discard the eluent (or save this radium fraction in case issues arise with the subsequent separation and measurement of ²²⁸Ra via ²²⁸Ac).
- 7.3.8. Repeat step 7.3.7.
- 7.3.9. Ensure that clean, labeled beakers or vials are below each column.
- 7.3.10. Add 10mL of 0.35M HNO₃ to elute actinium from each column.
- 7.3.11. Prepare sources for counting by rare earth micro precipitation using Eichrom Method SPA-01.

8. CALCULATIONS

Calculate ²²⁸Ra in pCi/L:

$$^{228}\text{Ra (pCi/L)} = \frac{A}{2.22 \times E \times V \times Y \times e^{-\lambda t_1}} \times \frac{\lambda t_2}{1 - e^{-\lambda t_2}}$$

where:

- A = net count rate, cpm
E = counting efficiency expressed as fraction
Y = ¹³³Ba (Ra) yield expressed as fraction
V = Sample volume (liters)
t₁ = decay time of ²²⁸Ac, from start of rinse until start of counting (minutes)
t₂ = counting time (minutes)
λ = decay constant of ²²⁸Ac (1.88*10⁻³ min⁻¹)

9. PERFORMANCE DATA


Reported by Burnett, et al. (1995)

Sample	Expected Result Ra-228 (pCi/L)	Experimental Result Ra-228 (pCi/L)
EMSL Ra in Water, 7/17/92	16.7 ± 3.3	17.6 ± 0.7
EMSL Standard Ra-228 Solution	59.3 ± 2.1	58.0 ± 1.3
Decontamination Experiment ^a	29.7 ± 1.0	31.7 ± 0.5
EMSL Performance Evaluation "A" 4/19/94 ^b	20.1 ± 5.0	21.6 ± 1.2

- a. The decontamination experiment consisted of using approximately 800 mL of EMSL Performance Evaluation "B" (April 20, 1993) which contained the following nuclides: ^{90}Sr (^{90}Y) = 397 dpm; ^{60}Co = 466 dpm; ^{134}Cs = 438 dpm; ^{137}Cs = 438 dpm; and ^{226}Ra 213 dpm. In addition the following were added: approximately 8000 dpm ^{133}Ba , 900 dpm ^{207}Bi , 500 dpm ^{210}Pb , and 65.9 dpm ^{228}Ra .
- b. Ra-226 was analyzed via radon emanation by collection of the sample load and rinse from the Ln Resin column. Their result of 20.1 +/- 1.1 pCi/L compared very well to the expected result of 20.0 +/- 3.0 pCi/L.

10. REFERENCES

- 1) Burnett, W.C., P.H. Cable, and Russ Moser, "Determination of Radium-228 in Natural Waters Using Extraction Chromatographic Resins," *Radioactivity & Radiochemistry*, Vol. 6, No. 3, pp. 36-43 (1995).
- 2) Maxwell, S.L., et al., "Rapid method for determination of 228Ra in water samples," *J. Radioanal. Nucl. Chem.* 295, 2181-2188, (2013).



- 1) Aliquot up to 1000mL of water into glass beaker.
- 2) Add Ba-133 tracer and 1mL of Ba carrier.
- 3) Add 10mL of concentrated H₂SO₄.
- 4) Cover with watch glass and heat to 90°C for 2 hours.
- 5) Allow samples to cool. Allow fine white BaSO₄ precipitate to settle or centrifuge samples to separate more quickly.
- 6) Decant supernate to waste.
- 7) Transfer precipitate to 15mL Pyrex glass centrifuge tube.
- 8) Wash precipitate with 4-5mL water. Centrifuge. Check pH of supernate. Decant supernate and discard as waste. Repeat until pH ~6.
- 9) Add 1mL 50% K₂CO₃ and 2-3mL of water, completely covering precipitate. Heat over flame until dry. Allow tube to cool. Repeat two more times.
- 10) Add 15-20mL of water. Mix well. Centrifuge. Check pH. Decant supernate to waste. Repeat until supernate pH ~7.
- 11) Add 5mL of 0.095M HNO₃. Gently heat to dissolve residue. If solids remain. Centrifuge. Decant supernate to 20mL LSC vial.
- 12) Check Ba-133 recovery. If less than 80%, repeat 9-11 on solids from step 11.
- 13) Allow at least 30 hours for Ac-228 ingrowth.

