

# SEPARATION OF AMERICIUM FROM RARE EARTHS

## 1. SCOPE

- 1.1. This is a method for separation of americium (and curium) from rare earths prior to alpha spectrometry source preparation.
- 1.2. This method does not address all aspects of safety, quality control, calibration or instrument set-up. However, enough detail is given for a trained radiochemist to achieve accurate and precise results for the analysis of the analyte(s) from the appropriate matrix, when incorporating the appropriate agency or laboratory safety, quality and laboratory control standards.

## 2. SUMMARY OF METHOD

- 2.1. Americium is separated from rare earths using Eichrom TEVA resin prior to measurement by alpha spectrometry. Am and Cm are retained on TEVA from a thiocyanate solution, while rare earths pass through the resin. Am and Cm are then recovered in a matrix that will allow alpha spectrometry source preparation by rare earth fluoride micro precipitation (Eichrom SPA01) or electrodeposition (Eichrom SPA02).

## 3. SIGNIFICANCE OF USE

- 3.1. This method is used to improve the alpha spectrometry peak resolution for americium (and curium) samples containing significant amounts (>100µg) of rare earths.

## 4. INTERFERENCES

- 4.1. Rare earths present in samples may not be separated from americium or curium using standard separation chemistries on TRU or DGA Resins. Light lanthanides, such as La and Ce, can be separated from Am/Cm on DGA resin by rinsing with 10-15 bed volumes of 3M HCl. However, heavier rare earths and yttrium and scandium will co-elute with Am/Cm. The presence of rare earths (>100µg) can lead to poor alpha spectrometry peak resolution through self-absorption in alpha spectrometry sources. This method will effectively separate Am and Cm from rare earths.

## 5. APPARATUS

- Analytical balance, 0.0001 g sensitivity
- Beakers, glass
- Centrifuge tubes, 50mL
- Column Rack (Eichrom Part: AC-103)
- Extension funnels, 25mL (Eichrom Part: AC-120)
- Fume hood
- Hotplate with stirrer
- Vortex mixer

## 6. REAGENTS

**Note: Analytical grade or ACS grade reagents are recommended.**

<i>Ammonium thiocyanate</i> , NH <sub>4</sub> SCN
<i>Deionized water</i> , all reagents are prepared with deionized water
<i>Formic acid</i> , HCOOH
<i>Hydrochloric acid (37%)</i> , concentrated HCl
<i>TEVA® resin</i> , 2mL prepacked column, 100-150 μm, Eichrom Part TE-C50-A

- 6.1. *Ammonium thiocyanate (4M)* - formic acid (0.1M) - Dissolve 304g NH<sub>4</sub>SCN in water, add 4.25mL formic acid and dilute to 1L with water. Prepare fresh prior to use.
- 6.2. *Ammonium thiocyanate (1.5M)* - formic acid (0.1M) - Dissolve 114g NH<sub>4</sub>SCN in water, add 4.25mL formic acid and dilute to 1L with water. Prepare fresh prior to use.
- 6.3. *Hydrochloric acid (1M)* - Add 83mL concentrated HCl to 750mL of water and dilute to 1L with water.

## 7. PROCEDURE

- 7.1. Evaporate each americium/curium sample to dryness in a glass beaker.
- 7.2. Dissolve the residue in 5mL 4M NH<sub>4</sub>SCN/0.1M formic acid
- 7.3. For each sample dissolved, place a TEVA column in the column rack. Place a waste reservoir below each column.

- 7.4. Remove the cap and bottom plug from each column, push the top frit down to the top of the resin bed, and allow each column to drain.
- 7.5. Add 5mL of 4M NH<sub>4</sub>SCN/0.1M formic acid into each TEVA column to condition the resin. Allow the solution to drain.
- 7.6. Transfer each solution from step 7.2. to the appropriate TEVA column. Allow the load solution to drain through column.
- 7.7. Add 3mL of 4M NH<sub>4</sub>SCN/0.1M formic acid to rinse each sample beaker. Add the rinse solution to the appropriate TEVA column. Allow the rinse to drain through each column.
- 7.8. Add 10mL of 1.5M NH<sub>4</sub>SCN/0.1M formic acid to rinse each TEVA column. Allow the rinse solution to drain through each column. Discard the solution collected to this point as waste.

**Note: Ammonium thiocyanate is incompatible with KClO<sub>3</sub> and Pb(NO<sub>3</sub>)<sub>2</sub>. Do not mix the waste solution with either of these chemicals.**

- 7.9. Ensure that clean, labeled centrifuge tubes are below each column.
- 7.10. Add 20mL 1M HCl to elute americium and curium.
- 7.11. Prepare samples for measurement of americium and curium by electrodeposition (Eichrom SPA02) or rare earth fluoride micro precipitation (Eichrom SPA01).

## 8. REFERENCES

- 1) Horwitz, E.P., et al. "Separation and Preconcentration of Actinides from Acidic Media by Extraction Chromatography," *Analytica Chimica Acta*, 281, 361-372 (1993).
- 2) Chiarizia, R., et al. "Am(III) and Eu(III) Extraction by Aliquat-336 and Benzyl Substituted Quaternary Ammonium Salts from Nitrate and Thiocyanate Solutions," *Solvent Extraction and Ion Exchange*, 13(4), 614-645 (1995).

