



**The Capacity Factor (k'):
How it is measured and
what it does not tell you.**

**E. Philip Horwitz
Senior Consulting Scientist
Eichrom Technologies, Inc.
Darien, IL 60561**

Acid dependency of k' for various ions at 23-25°C.
Sr Resin



Dry Weight Distribution



Weigh Resin (g)

Add Volume (mL) of
Solution with A_0

Mix

Withdraw Volume
of Solution to Measure A_s

$$D_w ? \frac{A_0 - A_s}{w(g)} \bigg/ \frac{A_s}{v(\text{mL})}$$

How is k' calculated from D_w ?

One needs to know:

- 1. Density of the extractant or combination of extractants used to prepare a given resin.**
- 2. Weight of extractant sorbed on the inert support.**
- 3. Density of the loaded resin.**
- 4. Packing density (weight of resin per unit volume of column).**
- 5. Bed volume of column.**

Relationship Between k' and D_w

$$D_w = \frac{A_o - A_s}{A_s} \cdot \frac{\text{mL}}{\text{g}} \quad (1)$$

$$\frac{D_w}{D} = \text{mL of stationary phase per gram of resin} \quad (2)$$

$$k' = D \cdot \frac{v_s}{v_m} \quad (3)$$

$$k' = D_w / \text{Conversion Factor} \quad (4)$$

- D_w = weight distribution ratio
 D = volume distribution ratio
 k' = number of F.C.V. to peak maximum
 $A_o - A_s$ = activity sorbed on g grams of resin
 A_s = activity in mL of solution
 v_s and v_m = volume of stationary and mobile phase, respectively, in a packed column

What k' does not tell you.

Breakthrough!

**Breakthrough is a function
of:**

- 1. k'**
- 2. Particle size**
- 3. Flow rate**
- 4. Temperature**

Comparison of Elution Curves for Sr^{2+} for Two Particle sizes of Sr Resin
Eluent 3.2 M HNO_3 , 23-24°C

