

Rapid Determination of Pu, Np, and U in 1-8L Seawater Samples

Summary of Method Plutonium, Neptunium, and Uranium are separated and concentrated from up to 8L samples of seawater with a hydrous titanium oxide precipitation, enhanced with 5mg of lanthanum and 125mg of ferric iron. A second precipitation with lanthanum fluoride removes additional matrix ions, and Uranium and Pu+Np are separated from potentially interfering radionuclides in the sample using stacked 2mL cartridges of Eichrom TEVA and TRU Resins. Isotopic U and Pu+Np are measured by alpha spectrometry following cerium fluoride microprecipitation onto Eichrom Resolve® Filters. Chemical yields are determined by recovery of ^{232}U and ^{242}Pu (or ^{236}Pu if measuring ^{237}Np) tracers. Recoveries of ^{232}U average $95 \pm 6\%$, while ^{236}Pu average $90 \pm 9\%$. Measured values of ^{238}U , ^{239}Pu , and ^{237}Np typically agree to within 10% of reference value. A single operator can process batches of 12 samples through alpha source preparation in 6-8 hours. Alpha spectrometry count times will vary depending on desired detection limit and data quality objectives.

Reagents

TEVA Resin, 2mL Cartridges (Eichrom TE-R50-S)
 TRU Resin, 2mL Cartridges (Eichrom TR-R50-S)
 Nitric Acid (70%)
 Hydrochloric Acid (37%)
 Hydrofluoric Acid (49%) or Sodium Fluoride
 Ammonium Hydroxide (listed as 28% NH_3 or 56% NH_4OH)
 Iron Carrier (50mg/mL Fe, as ferric nitrate)
 Lanthanum and Cerium Carriers (1mg/mL)
 ^{232}U and ^{242}Pu (or ^{236}Pu if meas. ^{237}Np) tracers
 Oxalic acid/Ammonium Oxalate
 Deionized Water H_2O_2 (30%)
 10% (w:w) TiCl_3 2M $\text{Al}(\text{NO}_3)_3$
 Boric acid Sulfamic Acid
 NaNO_2 Ascorbic Acid
 Denatured Ethanol 1.25M $\text{Ca}(\text{NO}_3)_2$

Equipment

Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)
 Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)
 Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)
 Yellow Outer Tips (Eichrom AR-1000-OT)
 Resolve Filters in Funnel (Eichrom RF-DF25-25PP01)
 50mL and 250-500mL Centrifuge Tubes
 Centrifuge
 Stainless Steel Planchets with adhesive tape
 Alpha Spectrometry System
 Analytical Balance
 Vacuum Pump
 Heat Lamp

Figure 1. Sample Preparation








Figure 2. TEVA-TRU Separation and Alpha Source Preparation

<p>(1) Adjust valence states of actinides by adding the following reagents in the order listed (mix between steps):</p> <ul style="list-style-type: none"> -0.2mL 1.5M Sulfamic acid -0.01mL 50mg/mL Fe carrier -1.5mL 1M Ascorbic acid -1mL 3.5M NaNO₂ <p>(2) Precondition stacked 2mL TEVA + TRU cartridges with 5mL 3M HNO₃.</p> <p>(3) Load sample solution at ~1mL/min.</p> <p>(4) Rinse sample tube with 5mL 3M HNO₃. Add tube rinse to cartridges.</p> <p>(5) Rinse cartridges with 10mL 3M HNO₃.*</p> <p>(6) Separate TEVA and TRU cartridges.</p>	<p>(10) Rinse TRU with 20mL 4M HCl-0.2M HF.</p> <p>(11) Rinse TRU with 12mL 10M HNO₃.</p> <p>(12) Strip U from TRU with 15mL 0.1M ammonium bioxalate.</p> <p>(13) Add 0.5mL 10% TiCl₃ to U samples and 0.5mL 30% H₂O₂ to Pu/Np samples.</p> <p>(14) Add 50-100ug Ce carrier to each sample. Mix well. Add 1mL 49% HF. Mix well. Wait 15-20 minutes.</p> <p>(15) Set up Resolve® Filter Funnel on vacuum box.</p> <p>(16) Wet filter with 3mL 80% ethanol followed by 3mL DI water.</p> <p>(17) Filter sample.</p> <p>(18) Rinse sample tube with 5mL DI water and add to filter.</p> <p>(19) Rinse filter funnel with 3mL DI water.</p>	<p>(20) Rinse filter funnel with 1-2mL 100% ethanol.</p> <p>(21) Draw vacuum until filter is dry.</p> <p>(22) Remove filter from funnel assembly and mount filter on stainless steel planchet with adhesive tape.</p> <p>(23) Dry filter under heat lamp for 3-5 minutes.</p> <p>(24) Measure actinides by alpha spectrometry.</p>
--	---	--

*Adding 50uL of 30% H₂O₂ to tube rinse can improve Uranium recoveries and decontamination in Pu(Np) fractions.

Method Performance Pu, Np and U from Seawater

Analyte	Volume, L	Replicates	Tracer	% Tracer		Analyte(mBq/L)		% Bias
				Recovery	Reference	Measured		
²³⁹ Pu	2	5	²³⁶ Pu	91 ± 9	33.8	32.6 ± 1.4	-3.6	
²³⁹ Pu	4	1	²³⁶ Pu	86	16.9	16.2	-4.1	
²³⁹ Pu	8	2	²³⁶ Pu	87 ± 3	27.8	27.6 ± 0.5	-0.7	
²³⁷ Np	2	5	²³⁶ Pu	91 ± 9	17.4	17.7 ± 1.5	1.7	
²³⁷ Np	4	1	²³⁶ Pu	86	8.7	7.2	-17	
²³⁷ Np	8	2	²³⁶ Pu	87 ± 3	4.4	4.2 ± 0.4	-4.5	
²³⁸ U	2	5	²³² U	99 ± 2	51.8	49.3 ± 1.5	-4.8	
²³⁸ U	4	1	²³² U	86	25.9	25.0	-3.6	
²³⁸ U	8	2	²³² U	92 ± 5	96.3	94 ± 3	-2.4	

16 hour count times

References

1) Sherrod L. Maxwell, Brian K. Culligan, Jay B. Hutchinson, Robin C. Utsey, Daniel R. McAlister, "Rapid determination of actinides in seawater samples," *J. Radioanal. Nucl. Chem.*, 300(3), 1175-1189 (2014).