



Ari Lehmusvuori, PhD

29.4.2021

## Liquid Scintillation Counters



Hidex  
600 SL



Hidex  
300 SL



Hidex  
Triathler

## Microplate Reader



Hidex Sense  
Beta

## Gamma Counter



Hidex Automatic  
Gamma Counter

## Radiowater generator



## Sample preparation instruments



Hidex  
600 OX Oxidizer



Hidex  
Q-ARE 100

HIDEX

# Sample preparation




600 OX Oxidizer




Q-ARE 100plus

# Sample preparation




# Sample




soil




water



Nuclear power plant  
decommissioning




plant



tissue

 = radioisotope


# Sample




soil




water



Nuclear power plant  
decommissioning




plant



tissue



After sample prep



● = radioisotope

Homogenous and clear liquid sample  
for Liquid Scintillation Counting

# Why sample preparation?

- > Solid sample transformation into homogenous liquid solution for LSC measurement
- > Removal of impurities to decrease quenching and luminescence background
- > Concentration of the target radionuclide
- > Isolation of the target radionuclide from mixture of several different radionuclides

# Hidex Q-ARE

QUICK  
AUTOMATED  
RADIONUCLIDE  
EXTRACTION





# Extraction chromatography

Extraction chromatography (EXC) combines the selectivity of liquid-liquid extraction with ease-of-use of the solid phase extraction (SPE) chromatography




Figure 1. Extraction chromatography resin structure. The stationary phase that contains liquid extractant compound specific to the target radionuclide is impregnated on the inert support.

# Extraction chromatography

Extraction chromatography (EXC) combines the selectivity of liquid-liquid extraction with ease-of-use of the solid phase extraction (SPE) chromatography





Figure 1. Extraction chromatography resin structure. The stationary phase that contains liquid extractant compound specific to the target radionuclide is impregnated on the inert support.

# Extraction chromatography

Steps:

## 1. Conditioning

- Strong acid for activation of the resin

## 2. Sample loading


- Sample is applied and the target radionuclides are caught by the resin

## 3. Washing

- Impurities are washed away by acid

## 4. Elution

- Radionuclide is eluted and collected



# Extraction chromatography

Steps:

## 1. Conditioning

- Strong acid for activation of the resin

## 2. Sample loading

- Sample is applied and the target radionuclides are caught by the resin

## 3. Washing

- Impurities are washed away by acid


## 4. Elution

- Radionuclide is eluted and collected


Acid




1.




2.



3.



4.



● = target radionuclide

# Extraction chromatography

Steps:

## 1. Conditioning

- Strong acid for activation of the resin

## 2. Sample loading


- Sample is applied and the target radionuclides are caught by the resin

## 3. Washing

- Impurities are washed away by acid

## 4. Elution

- Radionuclide is eluted and collected



# Extraction chromatography

## Steps:

### 1. Conditioning

- Strong acid for activation of the resin

### 2. Sample loading


- Sample is applied and the target radionuclides are caught by the resin

### 3. Washing

- Impurities are washed away by acid

### 4. Elution

- Radionuclide is eluted and collected



# Extraction chromatography

Steps:

## 1. Conditioning

- Strong acid for activation of the resin

## 2. Sample loading


- Sample is applied and the target radionuclides are caught by the resin

## 3. Washing

- Impurities are washed away by acid

## 4. Elution

- Radionuclide is eluted and collected



# Extraction chromatography

Steps:

## 1. Conditioning

- Strong acid for activation of the resin

## 2. Sample loading


- Sample is applied and the target radionuclides are caught by the resin

## 3. Washing

- Impurities are washed away by acid

## 4. Elution

- Radionuclide is eluted and collected





# Traditional way. Manual Extraction chromatography



# Traditional way. Manual Extraction chromatography



# Traditional way. Manual Extraction chromatography




# Traditional way. Manual Extraction chromatography



# Automated Extraction chromatography

Q-ARE  
100plus



# Q-ARE 100plus

12 Reagent solutions

8 Samples

Elution fraction  
collection bottles  
(total 20)

1 pump for each  
sample

8 extraction  
chromatography  
columns



# Q-ARE 100plus

12 Reagent solutions

8 Samples

Elution fraction  
collection bottles  
(total 20)

1 pump for each  
sample

8 extraction  
chromatography  
columns



# Q-ARE 100plus

12 Reagent solutions

8 Samples

Elution fraction  
collection bottles  
(total 20)

1 pump for each  
sample

8 extraction  
chromatography  
columns





# Q-ARE 100plus

12 Reagent solutions  
6 + 6

8 Samples

1 pump for each  
sample

Elution fraction  
collection bottles  
(total 20)

8 extraction  
chromatography  
columns



# Q-ARE 100plus

12 Reagent solutions

8 Samples

Elution fraction  
collection bottles  
(total 20)

1 pump for each  
sample

8 extraction  
chromatography  
columns



# Q-ARE 100plus

Acid resistant and fume hood free design


Acid resistant pumps, valves, tubing and fittings (PEEK, PTFE)




Acrylic doors covering the pumps and columns

# Hidex Q-ARE liquid handling


Basic view



# Hidex Q-ARE 100plus liquid handling



# Automated Tandem chromatography



# Hidex Q-ARE 100plus software

Home > Single Column Mode

· Column Selection: **Left** | All | Right  
 · Number of Sequence: 1  
 · Process Name: Conditioning  
 · Reagent: 1) D.W.  
 · Mounted Position: E1  
 · Flow Rate: 0.5 ml / min  
 · Total Volume: 0.5 ml  
 · Lower Volume of Reagent: 0 ml

Left Column

1

Q = 0.1\*R+0.0

2

Q = 0.2\*R+0.0

3

Q = 0.1\*R+0.0

4

Q = 0.2\*R+0.0

SEQ No.	Process Name	Reagent	Mounted Position	Flow Rate (ml/min)	Total Volume (ml)	Lower Volume of Reagent (ml)
1	Conditioning	1) D.W	E1	0.5	0.5	0
2	Conditioning	1) D.W	E1	0.5	0.5	0
3	Conditioning	1) D.W	E1	0.5	0.5	0
4	Conditioning	1) D.W	E1	0.5	0.5	0
5	Conditioning	1) D.W	E1	0.5	0.5	0
6	Conditioning	1) D.W	E1	0.5	0.5	0
7	Conditioning	1) D.W	E1	0.5	0.5	0
8	Conditioning	1) D.W	E1	0.5	0.5	0
9	Conditioning	1) D.W	E1	0.5	0.5	0
10	Conditioning	1) D.W	E1	0.5	0.5	0

Right Column

5

Q = 0.1\*R+0.0

6

Q = 0.1\*R+0.0

7

Q = 0.2\*R+0.0

8

Q = 0.2\*R+0.0

SEQ No.	Process Name	Reagent	Mounted Position	Flow Rate (ml/min)	Total Volume (ml)	Lower Volume of Reagent (ml)
1	Conditioning	7) D.W	E1	0.5	0.5	0
2	Conditioning	7) D.W	E1	0.5	0.5	0
3	Conditioning	7) D.W	E1	0.5	0.5	0
4	Conditioning	7) D.W	E1	0.5	0.5	0
5	Conditioning	7) D.W	E1	0.5	0.5	0
6	Conditioning	7) D.W	E1	0.5	0.5	0
7	Conditioning	7) D.W	E1	0.5	0.5	0
8	Conditioning	7) D.W	E1	0.5	0.5	0
9	Conditioning	7) D.W	E1	0.5	0.5	0
10	Conditioning	7) D.W	E1	0.5	0.5	0

# Columns & Resins

- Column size 1 – 20 ml
- Prepacked and self-packed columns
- Eichrom resins and cartridges

Product Name	Applications		
Actinide Resin	Group actinide separations/gross alpha measurements	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
Cl Resin	Cl, I	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
Cs Resins	Cs	<a href="#">Technical Info</a>	
Cu Resin	Cu	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
DGA Resins	Actinides, Lanthanides, Y, Ra	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
Ion Exchange Resins	Analytical grade cation and anion exchange resins	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
Ln Resins	Lanthanides, Ra-228	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
MnO <sub>2</sub> Resin	Ra	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
Nickel Resin	Ni	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
Pb Resin	Pb	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
Pre-filter Resin	Organics removal	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
RE Resin	Th, U, Np, Pu, Am, Cm, rare earth elements	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
Sr Resin	Sr, Pb	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
TEVA® Resin	Tc, Th, Np, Pu, Am/lanthanides	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
Tritium Column	<sup>3</sup> H	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
TRU Resin	Fe, Th, Pa, U, Np, Pu, Am, Cm	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>
UTEVA® Resin	Th, U, Np, Pu	<a href="#">Technical Info</a>	<a href="#">Part Numbers</a>



# Hidex Q-ARE

## Q-ARE 100plus

- 8 samples simultaneous processing
- 8 pumps and 8 columns



## Q-ARE 50

- 4 samples simultaneous processing
- 4 pumps and 4 columns



# Automated sample preparation

- ✓ Reduced worker exposure to radioactivity and hazardous chemicals > Improved user safety
- ✓ Improved reliability
- ✓ Higher throughput
- ✓ Lower overall costs
  - higher throughput with less man-hours
- ✓ Shorter analysis time



For more information,  
please contact us

[ari.lehmusvuori@hidex.com](mailto:ari.lehmusvuori@hidex.com)

Lemminkäisenkatu 62  
FI-20520 Turku

[info@hidex.com](mailto:info@hidex.com)  
[www.hidex.com](http://www.hidex.com)

Tel. +358 10 843 5570