

eichrom[®]

A BRAND OF
EICHROM TECHNOLOGIES



Fundamentals of metal ion separations: Nuclear Medicine Examples

Daniel McAlister, Ph.D.

30 April 2021

Outline

- $^{90}\text{Sr}/^{90}\text{Y}$ and ^{82}Sr (PPT vs IX vs EXC, extractant bleed)
- $^{225}\text{Ac}/^{213}\text{Bi}$ (Standard COW vs MSIG)
- ^{225}Ac production (Targeting minor component)
- $^{177}\text{Lu}/^{161}\text{Tb}$ (Low separation factors)
- ^{68}Ga (Inorganic adsorbants and bonded silicas)
- ^{227}Ac from stainless steel encapsulated neutron sources
 - Unique selectivity and synergistic enhancement of extraction

Eichrom Website – Nuclear Medicine



PRODUCTS **TECHNICAL INFO** **RESOURCES**

+1 (630) 963-0320

CONTACT



Radiochemistry

Learn More >



Nuclear Medicine

Learn More >

Eichrom Website – Nuclear Medicine

eichrom

A BRAND OF
EICHROM TECHNOLOGIES

PRODUCTS

TECHNICAL INFO

RESOURCES

+1 (630) 963-0320

CONTACT

For more than 25 years, Eichrom scientists and partners have developed chromatographic separations for a wide range of nuclear medicine radionuclides. Click the links below for more information on key radionuclides.

- [Actinium-225](#)
- [Bismuth-213](#)
- [Fluorine-18](#)
- [Gallium-68](#)
- [Germanium-68](#)
- [Indium-111](#)
- [Lead-203/212](#)
- [Lutetium-177](#)
- [Radium-223](#)
- [Scandium-44/47](#)
- [Strontium-82](#)
- [Terbium-161](#)
- [Thallium-201 \(Lead-201\)](#)
- [Thorium-227](#)
- [Yttrium-86](#)
- [Yttrium-90](#)
- [Zirconium-89](#)



$^{90}\text{Sr}/^{90}\text{Y}$ and ^{82}Sr

- Precipitation vs IX vs EXC
- Multicolumn separations
- Extractant Bleed

⁹⁰Sr Separations

	Fuming HNO ₃ ppt	Cation Exchange	EXC (Sr Resin)
Method	Evaporate sample. Series of ppt: Fuming HNO ₃ , Fe(OH) ₃ , Ba(CrO ₄).	Load from dilute acid. Rinse with acid gradients to remove impurities. Recover Sr in high acid concentration.	Concentrate sample using IX or ppt. Dissolve in 8M HNO ₃ . Rinse with HNO ₃ – oxalic acid. Recover Sr in 0.05M HNO ₃ .
Limitations	Labor Intensive. Dangerous chemicals.	Low separation factors. Large columns and elution volumes.	Lower capacity than IX. Extractant bleed (1-octanol).
Advantages	???	Less labor. Less hazardous.	High separation factors. Small columns.
Waste	Highly acidic, large volumes. Mixed waste.	Large volumes of acidic aqueous waste. Spent IX resin.	Small volumes of acidic aqueous waste. Spent EXC resin.


Fuming Nitric Acid (Separation of Sr from Ca, Ba and Fission Products)

- 1) Add 10 mg Sr, Ba, Ca carriers.
- 2) Concentrate water samples by evaporation or carbonate ppt.
- 3) Dissolve sample in enough fuming nitric acid (90-95%) to make 80% HNO₃.
- 4) Mix, heat 5 min. Centrifuge. Decant Supernate
 - Note: **Explosions are likely to occur** if the supernate waste is mixed with other wastes that may contain organic compounds.
- 5) Repeat 2x
- 6) Fe(OH)₃ ppt to remove Ca (2-3x)
- 7) Ba(CrO₄) ppt to remove Ba (2-3x)
- 8) Sr-oxalate ppt to concentrate and mount final Sr source.

Simplified scheme needed for separation of ⁹⁰Sr/⁹⁰Y.

Coryell and Sugarman, Radiochemical Studies: The Fission Products (McGraw-Hill Book Co., New York, 1951), N.N.E.S. Div. IV, Vol. 9, Book 3, Paper 236.

Cation Exchange



Small separation factors


($\alpha_{Y/Sr} \sim 20$)

-Large columns

-Large elution volumes



IX very useful for concentrating metal ions from large volumes and convert from low to high acidity.

Remove common matrix ions such as Na⁺ and K⁺.



^{90}Y (Sr Resin, RE Resin)

Dietz, M.L., Horwitz, E.P., 2000. Applications of extraction chromatography in the development of radionuclide generator systems for nuclear medicine, *Ind. Eng. Chem. Res.* 39, 3181-3188.



Sr-90 retained on Sr Resin

Sr-90 recovered in 0.05M HNO₃.

Y-90 retained on RE Resin.

Y-90 recovered in dilute HCl.

Relatively high extractant bleed (octanol and TBP).

HNO₃ and Fe carryover into HCl strip of Y-90.⁹

Options for extractant bleed

- Bonded ion exchangers (Use after EXC)
 - Stable in highly acidic conditions.
- Bonded silica ion exchangers (Use after EXC)
 - Stable from pH ~ 2-10.
- Polymeric scavengers (Use after EXC)
 - Stable in highly acidic conditions
- Bonded silica scavengers (C18, use after EXC)
 - Stable from pH ~2-10
- Alumina / Inorganic ion exchangers (Use after EXC)
 - Can leach metal ions

Relative extractant bleed

High





Low

• Sr Resin	1-octanol	
• Pb Resin	isodecanol	
• TRU/RE	TBP	
• Ac Resin	Dipex	higher at high pH
• LN	HDEHP	higher at high pH
• LN2	HEH[HEP]	higher at high pH
• LN3	H[DTMPP]	higher at high pH
• TEVA	Aliquat-336	higher at high pH
• UTEVA	DAAP	
• DGA, Normal	TODGA	
• DGA, Branched	TEHDGA	

^{90}Y (Sr Resin, DGA Resin)

D.R. McAlister, E.P. Horwitz, "Extraction of Selected Metal Ions by Mixtures of Diglycolamides and Crown Ethers," Solv. Extr. Ion Exch., accepted (2020), <http://dx.doi.org/10.1080/07366299.2020.1831249>.



Sr-90 retained on Sr Resin

Sr-90 recovered in 0.05M HNO_3 .

Y-90 retained on DGA Resin.


Y-90 recovered in dilute HCl.

Lower extractant bleed (octanol).

2-3M HCl rinse prior to recovery of Y-90 eliminates HNO_3 carryover and Sr impurity.

Higher Y retention = smaller columns

^{90}Y (DGA-Sr Resin), TODGA replaces 1-octanol diluent



Sr-90 and Y-90 retained on single resin.

Sr-90 recovered in 2M HCl or 0.05M HNO₃.

Y-90 recovered in dilute HCl.


Lowest extractant bleed (no octanol).

2-3M HCl rinse prior to recovery of Y-90 eliminates HNO₃ carryover and Sr impurity.


Higher Y retention = smaller columns

D.R. McAlister, D. Silvestri, E. Rush, E.P. Horwitz, "Extraction of Selected Metal Ions by Mixtures of Diglycolamides and Crown Ethers," *Solv. Extr. Ion Exch.*, 39(2), 184-203 (2021).

Sr-82 Flowsheet (50 grams Rb target)



Acid dependency of k' for various ions at 23-25°C.



Dunin, A.V., Nerozin, N.A., Togaeva, N.R., Khamyanov, S.V., Shapovalov, V.V. 2014. Extraction of Sr-82, Raw Material for Radiopharmaceutical Production. Pharmaceutical Chemistry Journal, 48(6), 395-397.

Sr-82 Flowsheet (50 grams Rb target)

TABLE 1. Technical Characteristics of Strontium-82

Parameter	Norm
Specific activity of Sr-82, TBq/g (Ci/g), at least	0.9(25)
Activity concentration of Sr-82, GBq/mL (mCi/mL), at least	2.2(60)
Total chemical impurities, $\mu\text{g}/\text{cm}^3$, less than	20
Radionuclidic impurities, mCi/mCi ^{82}Sr , less than	
Rb-83	0.0015
Rb-84	0.0001
Rb-86	0.0015
Sr-85	0.0001
HCl or HNO ₃ concentration	0.05 – 0.5


Dunin, A.V., Nerozin, N.A., Togaeva, N.R., Khamyanov, S.V., Shapovalov, V.V. 2014. Extraction of Sr-82, Raw Material for Radiopharmaceutical Production. Pharmaceutical Chemistry Journal, 48(6), 395-397.



$^{225}\text{Ac}/^{213}\text{Bi}$


- Standard Cow Generator
- Multicolumn Selectivity Inversion Generator

$^{225}\text{Ac}^{3+} / ^{213}\text{Bi}^{3+}$



<u>Nuclide</u>	<u>Half Life</u>	<u>Decay</u>	<u>Production</u>
^{225}Ac	10 d	α (5.0 – 5.8 MeV)	Decay ^{229}Th Proton Spallation ^{232}Th
^{213}Bi	45.6 m	α (5.6 – 5.9 MeV), 2.2% β^- (1423 keV), 97.8%	Decay ^{225}Ac

Bi-213 Generators (Standard Cow)



Scheme of Ac-225/Bi-213 generator with Ac-225 distribution profile

Elution with 0.1M HCl/NaI


30-40 mCi Ac-225:

- $75 \pm 10\%$ Yield of Bi-213

- $< 0.0002\%$ Ac-225 impurity

Application to clinical levels of ^{225}Ac could require multiple small generators or different technique.

Multi-Column Selectivity Inversion Generators (MSIG)



Store ^{225}Ac in solution

Improved radiation resistance.

Disposable columns

Likely requires automation to reduce operator dose.

D. R. McAlister and E. P. Horwitz, "Automated two column generator systems for medical radionuclides," *Applied Radiation and Isotopes*, 67, 1985-1991, (2009).

Bi-213 Generators (Summary Table)

Generator	Resins	Yield	Ac impurity	Matrix
Cow	MP50	76%	<2 E-5 %	0.6 mL 0.1M HCl/HI
MSIG	UTEVA 50Wx8 Prefilter	87%	<1 E-7%	2.0 mL 0.5M NaOAc 0.75M NaCl, pH 4.0

A. N. Vasiliev, V. A. Zobnin, Yu. S. Pavlov & V. M. Chudakov (2020) Radiation Stability of Sorbents in Medical $^{225}\text{Ac}/^{213}\text{Bi}$ Generators, Solvent Extraction and Ion Exchange, DOI: [10.1080/07366299.2020.1846892](https://doi.org/10.1080/07366299.2020.1846892)




$^{232/229}\text{Th}/^{225}\text{Ra}/^{225}\text{Ac}$ (^{225}Ac Production Methods)


- Solvent Extraction vs Ion Exchange vs EXC
- Importance of Selectivity
- Large Target Separations

$^{225}\text{Ac}^{3+}$ Sources

<u>Nuclide</u>	<u>Half Life</u>	<u>Decay</u>	<u>Production</u>
^{233}U	1.592 E5 y	α (4.5 – 4.8 MeV)	Thermal Breeder Reactors: $^1_0\text{n} + ^{232}\text{Th} \rightarrow ^{233}\text{Th} \rightarrow ^{233}\text{Pa} \rightarrow ^{233}\text{U}$
^{229}Th	7932 y	α (4.5 – 5.1 MeV)	Decay ^{233}U
^{225}Ac	10 d	α (5.0 – 5.8 MeV)	Decay ^{229}Th Proton Spallation ^{232}Th
^{225}Ra	14.9 d	β^- (356 keV)	Decay ^{229}Th Proton Spallation ^{232}Th
^{213}Bi	45.6 m	α (5.6 – 5.9 MeV), 2.2% β^- (1423 keV), 97.8%	Decay ^{225}Ac
^{227}Ac	21.77 y	α (4.4 – 5.0 MeV), 1.38% β^- (44.8 keV), 98.62%	Decay ^{235}U Proton Spallation ^{232}Th



Spallation



High energy protons strip neutrons and fragments from thorium forming lighter nuclides.

Fragments can also combine with thorium to form heavier nuclides.

Ac-225 Separation Schemes

$^{232}\text{Th}(p,x)^{225}\text{Ac}$ (^{227}Ac impurity)

$^{232}\text{Th}(p,x)^{225}\text{Ra}$ (β^-) ^{225}Ac

Primary Separation:
(Th Removal)

SX (DA[AP]) – HNO_3

Anion Exchange (MP1 or 1x8) - HNO_3

Cation Exchange (50Wx8) - H_2SO_4 or Citrate

Secondary Columns: UTEVA/DGA, cation exchange



- Remove remaining Th and spallation byproducts

- Separate Ac/Ra



- Remove key byproducts/impurities: ^{230}U , $^{231/233}\text{Pa}$, Ca, Fe, ^{140}La ,
 $^{227}\text{Ac}???$

SX with DA[AP]

J. Harvey, J. A. Nolen, T. Kroc, I. Gomes, E. P. Horwitz, D. R. McAlister, "Production of Ac-225 via high energy proton induced spallation of Th-232," Proceedings of Application of high energy proton accelerators, Fermilab, Chicago, IL, October 19-21, eds. Rajendran Raja and Shekhar Mishra, pp. 321-326 (2010).



Th Selective Separations



Extracting Th (10-50g) from HNO₃ requires very large columns (1-2 L of 1x8 or MP-1).

High acid concentration and volume may require evaporation prior to next steps.


Harvey, J.H., Nolen, J., Vandergrift, G., Kroc, T., Gomes, I., McAlister D.R., Horwitz, E.P. 2011. Production of Actinium-225 via High Energy Proton Induced Spallation on Thorium-232. Final Technical Report DE-SC0003602. <https://www.osti.gov/scitech/servlets/purl/1032445/>

Ac Selective Separations

A.H. Bond, E.P. Horwitz and D.R. McAlister, "A Multicolumn Selectivity Inversion Generator for the Production of High Purity Actinium-225 for Use in Therapeutic Nuclear Medicine," August 8, 2006, United States patent number 7,087,206.

D.R. McAlister, E.P. Horwitz, "Selective Separation of Radium and Actinium from Bulk Thorium Target Material on Strong Acid Cation Exchange Resin from Sulfate Media," *Applied Radiation and Isotopes*, 140, 18-23 (2018).

Mastren, T., Radchenko, V., Owens, A., Copping, R., Boll, R., Griswold, J.R., Mirzadeh, S., Wyant, L.E., Brugh, M., Engle, J.W., Nortier, F.M., Birnbaum, E.R., John, K.D., Fassbender, M.E. 2017. Simultaneous Separation of Actinium and Radium Isotopes from a Proton Irradiated Thorium Matrix. *Nature Scientific Reports*, 7, 8216. doi:10.1038/s41598-017-08506-9



Reverses selectivity

- Rejects anionic Th
- Retains Ac/Ra


Allows for smaller columns.

Ac/Ra recovered in 5M HNO₃.

No evaporate required prior to next steps.

Dissolving Th metal in H₂SO₄?

Dissolution of Th in H₂SO₄/HF



H₂SO₄

HF

H₂O

Heat

H₂O


Low Solubility??!!

“Selective Separation of Radium and Actinium from Bulk Thorium Target Material,” D.R. McAlister, E.P. Horwitz, R. Perron, D. Gendron, P. Causey, J.T. Harvey, 11th International Symposium on Targeted Alpha Therapy, Ottawa, Ontario, Canada, April 1-4, 2019.

Solubility curves of the hydrates of thorium sulphate.

Th-sulfate solubility limit
~0.15 moles/L at 20C.


Hexa-, Octa-, Nona-,
hydrates more soluble at
high temperature.



Dissolution of Th in H₂SO₄/HF




H₂SO₄
HF
H₂O




Heat
H₂O
Low Solubility??!!

Cool. Mix.



0.6M H₂SO₄
0.03 M HF
0.1 M Th
pH 0.8 – 1.0



Th, Ac, Ra separation from H₂SO₄


(1) Precondition 5 BV
0.4M (NH₄)₂SO₄, pH 2.0

(2) Load Sample

(3) Rinse 5 BV 0.4M
(NH₄)₂SO₄, pH 2.0

(4) Rinse 5 BV 0.05M HNO₃

(5) Strip 5 BV 5.0M HNO₃



(1) (4) Waste

(2) (3) Th

(5) Ac, Ra

75 – 100 mL cation exchange column

1-2 L of feed for (25-50 grams Th target)


Recover Ac/Ra in 400-500 mL 5M HNO₃

No evaporation required for further processing.

D.R. McAlister, E.P. Horwitz, “Sulfate based system for the separation of Actinium and Radium from irradiated Thorium Target,” *Applied Radiation and Isotopes*, 140, 18-23 (2018).

Ac Polishing Steps


3-6M HNO₃




- UTEVA (Phosphonate)
 - removes remaining Th, Pa, U
- DGA
 - retains Ac³⁺, rare earths, Ca
- Ra²⁺ passes both columns (Fe, Ba, Al, many fission products)

Ac Separation from Calcium and La


k' on DGA, Normal Resin vs HNO_3



k' on DGA, Normal Resin vs HCl




Ac Polishing Steps



- UTEVA and DGA (and other reagents) can have small traces of Ca impurity precondition separately
- Ac/Ca will coelute from DGA, unless special care is taken to separate them.
- Conditioning UTEVA/DGA separately can eliminate much of the Ca impurity.

^{225}Ac Separation (8M HNO_3)



50-100 μm DGA, Normal Resin


Column

4.2 cm length

0.7 cm diameter

- Complete removal of Fe, Ba(Ra), Pb, Sr, Bi
- >95% removal of Ca
- La co-elutes with Ac

^{225}Ac Elution (breakthrough in calcium removal)



50-100 μm DGA, Normal Resin

Column

4.2 cm length

0.7 cm diameter

Cartridge

2.7 cm length

0.9 cm diameter


90% Ac

95% Ac


>99% Ac

0.1M HCl

^{225}Ac w/ Ca Impurity




^{225}Ac Purity (Ca Removal)




Separation of Ac, La, Ce, Pr and Nd on DGA Resin

Slurry Packed 25-53 μm DGA Resin, Operating Temperature 50(1) $^{\circ}\text{C}$




ppm/mL vs. Bed Volumes of Eluate

Slurry Packed 25-53 μm LN Resin, Preconditioned with 0.50 M (Na,H)OAc, 50(1) $^{\circ}\text{C}$




Ac Separation from La

k' on DGA, Normal Resin vs HNO_3



k' on DGA, Normal Resin vs HCl



Mastren, T., Radchenko, V., Owens, A., Copping, R., Boll, R., Griswold, J.R., Mirzadeh, S., Wyant, L.E., Brugh, M., Engle, J.W., Nortier, F.M., Birnbaum, E.R, John, K.D., Fassbender, M.E. 2017. Simultaneous Separation of Actinium and Radium Isotopes from a Proton Irradiated Thorium Matrix. Nature Scientific Reports, 7, 8216. doi:10.1038/s41598-017-08506-9



^{177}Lu and ^{161}Tb

- Rare Earth Separations (low α)
- Effect of Column Loading on Peak Shapes
- Milligrams Target to No-carrier added

177Lu/161Tb

<u>Nuclide</u>	<u>Half Life</u>	β_{\max}	β_{avg}	<u>Photons</u>	<u>Production</u>
¹⁷⁷ Lu	6.7 days	497 keV	130 keV	208 keV (10.4%) 113 keV (6.2%)	¹⁷⁶ Yb(n, γ) ¹⁷⁷ Yb(β^-) ¹⁷⁷ Lu
¹⁶¹ Tb	6.9 days	593 keV	154 keV	45 keV (18%) 49 keV (17%) 75 keV (10.2%)	¹⁶⁰ Gd(n, γ) ¹⁶¹ Gd(β^-) ¹⁶¹ Tb



E. P. Horwitz, D. R. McAlister, A. H. Bond, R. E. Barrans, J. M. Williamson, "A Process for the Separation of ¹⁷⁷Lu from Neutron Irradiated ¹⁷⁶Yb Targets," *Applied Radiation and Isotopes*, 63, 23-36 (2005).

E. P. Horwitz, D. R. McAlister, M. L. Dietz, "Extraction chromatography versus solvent extraction: How similar are they?" *Sep. Sci. and Technol.*, 41(10), 2163-2182 (2006).


D. R. McAlister, E. P. Horwitz, "The Characterization and Novel Applications of Extraction Chromatographic Materials Containing Bis(2-ethyl-1-hexyl)phosphoric Acid, 2-ethyl-1-hexylphosphonic acid, mono 2-ethyl-1-hexyl ester and 2,4,4-trimethyl-1-pentylphosphonic acid," *Solv. Extr. Ion Exch.*, 25(6), 757-769 (2007).

A. Dash, M. Raghavan, A. Pillai, F.F. Knapp Jr., "Production of Lu-177 for targeted therapy: Available options," *Nucl Med Mol Imaging.*, 49, 85-107 (2012).


Acidic Phosphorus Extractants (LN Series)



Phosphoric
Acid



Phosphonic
Acid



Phosphinic
Acid



Basicity of P=O






Acidity of P-OH



Separation Factors






Separation Factors




Separation of Yb and Lu on LN2 Resin

~2% column loading capacity


Slurry Packed 25-53 μm LN2 Resin, Operating Temperature 50(1) $^{\circ}\text{C}$



Lu/Yb Separation on LN2 Resin, 50°C, 50 mg Yb





~20% column loading capacity



Scale up to 200 mL column (2%)?

Lu-177 LN2 vs DGA




DGA compliments LN2

Concentrate Yb/Lu between LN2 columns.

Reduce acidity to avoid evaporations.

Removes common impurities (Ca, Al, Fe, Na, K).


Lu-177 (20 mL LN2 / 2-5 mL DGA)



Lu/Yb Separation on HEH[EHP] Resin (LN2) vs Yb mass

Yb mass (mg) Parameter	Yb				Lu			
	50 mg Yb	25 mg Yb	5 mg Yb	0.5 mg Yb	50 mg Yb	25 mg Yb	5 mg Yb	0.5 mg Yb
Peak	3.00	3.34	5.10	6.38	7.33	8.85	10.13	10.62
Gaussian Width	2.02	1.67	1.30	1.11	6.63	2.99	1.77	1.81
Resolution (Yb/Lu)	1.00	2.36	3.28	2.90				
Number of Plates	35	64	246	529	20	140	524	551
Plate Height (cm)	0.5951	0.3281	0.0853	0.0397	1.0738	0.1498	0.0401	0.0381
% Column Loading	20%	10%	2%	0.2%	20%	10%	2%	0.2%

Tb Isotopes (Tb-161)



Larger separation factor for
Tb/Ga (4.7)

Better separation on smaller
column (5 mL) with higher
loading (75-80%)

>98% Tb

<0.06% Gd

D.F. = 1700 for single column



Purification of ^{68}Ga from SnO_2 generator.
Purification of ^{68}Ga labeled compounds.
(Bonded Silicas)

^{68}Ga Generator

Ge-68 adsorbed on a hydrous stannic oxide column.

Elute Ga-68 ($t_{1/2} = 68$ min) with 0.1M HCl.

Ga-68 needs to be purified to ensure good labeling and purity.

- Ge-68
- Stable metals Sn, Fe

Recover in a matrix suitable for labeling reaction.

Bonded functional groups to ensure low organic bleed.



Tworowska, I., Ranganathan, D., THamake, S., Delpassand, E., Mojtahedi, A., Schultz, M., Zhernosekov, K., Marx, S., (2016). Radiosynthesis of clinical doses of ^{68}Ga -DOTATATE (GalioMedixTM) and validation of organic-matrix-based $^{68}\text{Ge}/^{68}\text{Ga}$ generators, *Nuclear Medicine and Biology*, 43, 19-26.

^{68}Ga Generator


Collect Ga-68 from generator on SCX-silica cartridge in 0.1M HCl.

Recover Ga-68 from SCX with 5M NaCl/0.1M HCl.


- low acid to allow labeling rxn at pH ~4.5 with addition of acetate or ascorbate buffer

Label Ga-68 to small molecule.

Scavenge free Ga-68 from pH 4-5 buffer with WCX-silica.



Strong Cation Exchange




Weak Cation Exchange



Separation of ^{227}Ac from Ac(Be) neutron sources in stainless steel capsule.


Influence of FeCl_3 on the uptake of Am(III) from HCl



FeCl_3 increases the uptake of Am(III) more than expected for the additional Cl^-



Influence of Metal Ion on the uptake of Am(III) from 1M HCl





Other metal ions that form anionic chloride complexes show similar effect.

Bi, Tl, Fe, In, Ga.

Al shows only a modest increase due to additional Cl^- from AlCl_3 .


Equilibria for the synergistic extraction of Am(III) from HCl + FeCl₃



M.A. Antonio, D.R. McAlister, E.P. Horwitz, “Europium(III) Diglycolamide Complex: Insights into the Coordination Chemistry of Lanthanides in Solvent Extraction,” *Dalton Transactions*, 44(2), 515-521 (2015).


k' Ac(III) on DGA Resin, normal vs. HCl

50-100 mm, 1 h contact time, 22(1)°C




^{227}Ac from Stainless Steel

2 mL cartridge DGA, 2 mL cartridge prefilter, 2 mL/min, 22(1) $^{\circ}\text{C}$



^{227}Ac from Stainless Steel



^{227}Ac from Stainless Steel



1st Separation (DGA + Prefilter, HCl)

Metal Impurities in Ac-227 from 316 Stainless Steel

Fraction	Al (mg)	Cr (mg)	Fe (mg)	Mg (mg)	Mn (mg)	Mo (mg)	Ni (mg)
Load ¹	0.12	225	1100	<0.0005	31	153	170
Rinse ²	0.07	4	46	<0.0003	0.38	148	3.2
Strip ³	0.0006	0.006	0.2	<0.0002	0.0004	0.13	0.0016
DF	N/A	24000	2800	N/A	40000	190	70000

¹2 grams of 316 Stainless Steel, dissolved in HCl/H₂O₂, adjusted to 50 mL 3 M HCl

²25 mL of 3 M HNO₃

³20 mL of 0.1 M HCl

Polishing Step (DGA, HNO₃)

Metal Impurities in Ac-227 from 316 Stainless Steel

Fraction	Al (mg)	Cr (mg)	Fe (mg)	Mg (mg)	Mn (mg)	Mo (mg)	Ni (mg)
Load ¹	ND	0.0055	0.085	ND	ND	0.0065	ND
Oxalate Rinse ²	ND	0.003	0.043	0.016	0.0002	0.003	ND
HNO ₃ Rinse ³	ND	ND	ND	ND	ND	ND	ND
Strip ⁴	ND	ND	ND	ND	ND	ND	ND
Total DF	N/A	>10 ⁷	>10 ⁷	N/A	>10 ⁶	>10 ⁶	>10 ⁶

¹Product from 1st separation acidified to 40 mL 3 M HNO₃

²25 mL of 3 M HNO₃ + 0.05 M oxalic acid

³25 mL of 3 M HNO₃

⁴20 mL of 0.1 M HCl

Conclusions

- The unique selectivity of EXC resins can simplify purification of radionuclides for nuclear medicine applications.
- Organic leaching must be considered when using EXC resins for nuclear medicine applications.
- Combinations of multiple EXC and IX columns are often needed to meet the high purity requirements for nuclear medicine.
- Ion exchange or solvent extraction can compliment EXC when high capacity or low organic impurity levels are required.

References

- 1) “Separation of Selected Nuclear Medicine Isotopes Using Extraction Chromatography,” D.R. McAlister, E.P. Horwitz, 64th Radiobioassay and Radiochemical Measurements Conference, Santa Fe, NM, Oct 27- Nov 1, 2019. http://www.eichrom.com/wp-content/uploads/2019/11/7_McAlister_RRMC-2019_Nuclear-Medicine-Separations.pdf
- 2) “Selective Separation of Radium and Actinium from Bulk Thorium Target Material,” D.R. McAlister, E.P. Horwitz, R. Perron, D. Gendron, P. Causey, J.T. Harvey, 11th International Symposium on Targeted Alpha Therapy, Ottawa, Ontario, Canada, April 1-4, 2019.
- 3) “Separation and Counting Options for 89/90Sr,” D.R. McAlister, 63rd Radiobioassay and Radiochemical Measurements Conference, Portland, ME, May 20-25, 2018.
- 4) D. R. McAlister and E. P. Horwitz, “Automated two column generator systems for medical radionuclides,” Applied Radiation and Isotopes, 67, 1985-1991, (2009).

Thank you

Questions???