

# Validation of Extraction Chromatography Resin Material Identity Using Fourier Transform Infrared Spectroscopy (FTIR)

Madeleine A. Eddy, Jose Muñoz, and Daniel R. McAlister  
Eichrom Technologies, LLC Lisle, IL USA


## Abstract

Eichrom Technologies manufactures extraction chromatography resins for use in a range of applications, including nuclear medicine. These resins are used to purify radionuclides from target material and other contaminants prior to radiolabeling. Eichrom conducts resin-specific quality control designed to confirm the material identity and performance specifications for each product. However, the final resin materials are nearly identical to each other visually, and customer interest was expressed for a simple, non-destructive identity test as an additional check for all end-users. Herein, the use of Fourier Transform Infrared Spectroscopy (FTIR) to characterize Eichrom's extraction chromatography resins is presented as an option for identity validation of the materials. This simple technique can provide confirmation and distinguish between most Eichrom resins or classes of resins that use similar or identical extractants. Further work to improve the quality of the spectral library is ongoing, with hopes to improve the differentiation between SR/PB, RE/TRU, and DGA-N/DGA-B/DOODA resins.

## Eichrom EXC Resins


### Standard Formulation

- AmberChrom CG71M solid support
- 40% wt. liquid extractant or solid extractant + diluent
  - Physiosorbed to support




### Resins for Actinides

- TEVA: +4 actinides
- UTEVA: +4 actinides and uranium (UO<sub>2</sub><sup>2+</sup>)
- AC: high efficiency gross actinide extraction




### LN Series

- LN, LN2, LN3: similar chemistry
- Acidic organophosphorus extractants
- High SF for adjacent LN's
- Extracts from dilute acid




### DGA Family

- DGA-Normal, DGA-Branched, DOODA
- All C8-amides with ethylene linkers
- Extracts 3+ and 4+ actinides/lanthanides from acidic media




### Crown Ether Resins

- Selective for Sr and Pb
- Requires alcohol diluent to fully dissolve crown ether
  - SR = octanol
  - PB = isodecanol



### Resins for f-Elements

- TRU: transuranic actinides
- RE: same extractant + diluent as TRU, in different ratios
- RE shows improved uptake of +3 f-elements compared to TRU



## Motivation / Background

- Extraction chromatography resins are typically odorless white powders
  - Hard to distinguish different resins based on physical appearance
  - In-house resin QC employs radiochemical separations
- Nuclear medicine customers may require validation of resin ID on receipt of material
- A simple, non-destructive method for confirming the EXC material is the correct formulation is needed
- Recently acquired a CARY 630 FTIR instrument with ATR cell

## ATR-FTIR Spectroscopy

- FTIR spectra of resins represent the chemical features of the extractant, diluent (if present), and solid support
- Peak positions correspond to functional groups


Functional Group	Absorption (cm <sup>-1</sup> )	Intensity	Functional Group	Absorption (cm <sup>-1</sup> )	Intensity
Alkane	C-H 2850-2960	Medium	Carbonyl compound	C=O 1670-1780	Strong
Alkene	=C-H 3020-3100	Medium	Aldehyde	1730	Strong
Alkyne	≡C-H 3300	Strong	Ketone	1715	Strong
Alkyl halide	C-Cl 600-800	Strong	Ester	1735	Strong
	C-Br 500-600	Strong	Amide	1690	Strong
	C-I 3400-3650	Strong, broad	Carboxylic acid	1710	Strong
Alcohol	O-H 1050-1150	Strong	Carboxylic acid	O-H 2500-3100	Strong, broad
Arene	C-H 3030	Weak	Nitrile	C≡N 2210-2260	Medium
Aromatic ring	1660-2000	Weak	Nitro	NO <sub>2</sub> 1540	Strong, broad
	1450-1600	Medium	Amine	N-H 3300-3500	Medium
				C-N 1030-1230	Medium

## Identity Testing Methods

- Create an Eichrom library with extractants, resins, diluents, and solid supports
- Create a method to match-fit unknown samples using the Eichrom library
- Prepare unknown samples containing Eichrom resins or misc. resins not in Eichrom Library
- Recruit QC chemist (JAM) to measure the unknown samples and use match-fits from the Eichrom library to determine the correct resin ID

## Cary 630 FTIR

- Portable (3.5 kg)
- Agilent MicroLab software
- < 2 cm<sup>-1</sup> resolution
- Diamond ATR attachment



## Match-fitting



## Results

Sample	Resin	Resin Fit	First Fit	Second Fit	Third Fit	JAM Resin ID
C	Amberchrom	0.92292	Amberchrom	TEVA	ZR	Amberchrom
F	TEVA	0.87270	TEVA	WBEC	CL	TEVA
M	UTEVA	0.88755	UTEVA	RE	TRU	UTEVA
K	AC	0.90213	UTEVA	LN	LN2	AC
O	LN	0.94367	LN	TRU	UTEVA	LN
B	LN2	0.93417	LN2	UTEVA	RE	LN2
V	LN3	0.90316	LN3	TEVA	CL	LN3
R	DGA-N	0.85212	DOODA	DGA-N	TEVA	DGA-N
A	DGA-B	0.84132	DGA-N	DOODA	DGA-B	DGA-B
G	DOODA	0.85303	DOODA	DGA-N	DGA-B	DOODA
D	SR	0.88180	PB	SR	TEVA	SR
U	PB	0.89244	PB	SR	TEVA	PB
J	TRU	0.87032	TRU	RE	UTEVA	TRU
H	RE	0.85859	TRU	RE	UTEVA	RE

## Conclusions

- ATR-FTIR can determine the identity of most extraction chromatography resins
- Spectra are a combination of signatures from the amberchrom support, organic extractant, and diluent (if included) so similar resins have similar spectra
  - Some difficulty distinguishing between DGA resin family
  - Some difficulty distinguishing between SR and PB resin
  - Some difficulty distinguishing between TRU and RE resin
- In instances where the software incorrectly matched the resin, the user compared spectral fits manually and was able to correctly assign the resin ID
- This basic test case suggests that this technique could be used by nuclear medicine customers as a form of identity testing to confirm resin ID on receipt of materials, though some work to improve the fitting parameters may be required.

## References

McMurtry, John E. Organic Chemistry: A Tenth Edition - OpenStax adaptation 1. OpenStax, 2023.

"What is FTIR Spectroscopy? Principles Overview." Agilent Technologies, www.agilent.com/en/support/molecular-spectroscopy/ftir-spectroscopy/ftir-spectroscopy-basics-faqs.