

## Certificate of Analysis


No. CA12074

Product: Eichrom's Diphonix Resin  
Lot #: 033111DPM

Eichrom's Diphonix Resin, Lot 033111DPM, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specifications</u>	<u>Result</u>
% Solids	25-35 %	34.8
Phosphorous	$\geq 1.3$ eq/kg	1.33
Acid Capacity	$\geq 5.3$ eq/kg	5.35
Iron Loading Capacity	TBD	15.1 g/L



April 23, 2012

Jill Bryant  
Quality System Coordinator