



# Certificate of Analysis

No. CA 15107

Product: Eichrom's Anion Exchange Resin

Lot #: 051915A4C

Eichrom's Anion Exchange resin, Lot # 051915A4C has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White to light Yellow	Pass
Extractables	< 1 mg (residue) / g (resin)	n.d.
pH	≥ 4	Pass
Wet Volume Capacity	≥ 1.0 meq/mL	1.2

Jill Bryant  
Quality System Coordinator

June 12, 2015