



Certificate of Analysis

No. CA 14114

Product: Eichrom's Cation Exchange Resin

Lot #: 063014C8C

Eichrom's Cation Exchange resin, Lot # 063014C8C has been manufactured according to Eichrom standard operating procedures. This resin has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	Tan to orange-brown	Pass
Extractables	< 1 mg(residue) / g(resin)	0.55 mg/g
pH	≥ 4.0	pass

Jill Bryant
Quality System Coordinator

July 10, 2014