



Certificate of Analysis

No. CA 14166

Product: Eichrom's Anion Exchange Resin

Lot #: 082014A8M

Eichrom's Anion Exchange resin, Lot # 082014A8M has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White to light Yellow	Pass
Extractables	< 1 mg (residue) / g (resin)	0
pH	≥ 4	Pass
Wet Volume Capacity	≥ 1.2 meq/mL	1.4



Jill Bryant
Quality System Coordinator

September 24, 2014