

Certificate of Analysis


No. CA 17200

Product: Eichrom's Cation Exchange Resin; 50Wx8, 100-200 mesh, H⁺ Form
Date of Manufacture: 15-Sep-2017
Lot #: 091517C8F

Eichrom's Cation Exchange resin, Lot 091517C8F, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

Test Value	Specification	Result
Extractables	< 1 mg(residue) / g(resin)	0.30 mg/g
pH	≥ 4.0	pass
Wet Volume Capacity	min. 1.7 meq/mL	1.9 meq/mL



Jill Bryant
Quality System Coordinator

December 21, 2017