



Certificate of Analysis

No. CA 13152


Product: Eichrom's Cation Exchange Resin

Lot #: 091813C8F

Eichrom's Cation Exchange resin, Lot # 091813C8F has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	Tan to orange-brown	Pass
Extractables	< 1 mg(residue) / g(resin)	0.53 mg/g
pH	≥ 4.0	pass



Jill Bryant
Quality System Coordinator

September 24, 2013