

# Certificate of Analysis

No. CA 13191

Product: Eichrom's Cation Exchange Resin

Lot #: 100913C8M

Eichrom's Cation Exchange resin, Lot # 100913C8M has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	Tan to orange-brown	Pass
Extractables	< 1 mg(residue) / g(resin)	0.03 mg/g
pH	≥ 4.0	pass
Wet Vol. Capacity	Minimum 1.7 meq/mL	2.1

  
Jill Bryant  
Quality System Coordinator

November 11, 2013