

Certificate of Analysis


No. CA 17058

Product: Eichrom's Anion Exchange 2cc Resin in 7mL Cartridge
Part: A47ML-R50-M-CL-2CC-G
Lot: A4MR17A

Eichrom's Anion Exchange 2cc resin in 7mL cartridge, lot A4MR17A, has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, 110916A4M, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White to light Yellow	Pass
Extractables	< 1 mg (residue) / g (resin)	n.d.
pH	≥ 4	Pass
Wet Volume Capacity	≥ 1.0 meq/mL	1.0



Jill Bryant
Quality System Coordinator

April 11, 2017