



**eichrom**  
TECHNOLOGIES

A GCI COMPANY

# Certificate of Analysis

No. CA 14190

Product: Eichrom's Anion Exchange Resin cartridges

Lot #: A8MR14E

Eichrom's Anion Exchange resin cartridges, Lot # A8MR14E has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

| <u>Test Value</u> | <u>Specification</u>       | <u>Result</u> |
|-------------------|----------------------------|---------------|
| Appearance        | White to light Yellow      | Pass          |
| Extractables      | < 1 mg(residue) / g(resin) | 0.02 mg/g     |
| pH                | ≥ 4                        | Pass          |

Jill Bryant  
Quality System Coordinator

October 16, 2014