

## Certificate of Analysis

No. CA06154


Product: Eichrom's Actinide Resin Columns

Lot #: AC6A

Eichrom's Actinide Resin Columns Lot # AC6A has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.5 M HCl acid	$\geq 300$	426



Jill Bryant  
Quality Control Chemist

Prepared: December 20, 2006