

## Certificate of Analysis

No. CA09076


Product: Eichrom's Actinide Resin

Lot #: ACA060209

Eichrom's Actinide Resin Lot # ACA060209 has been manufactured according to Eichrom standard operating procedures. This resin lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.5 M HCl acid	$\geq 300$	419

  
\_\_\_\_\_  
Jill Bryant  
Quality System Coordinator

June 17, 2009