

Certificate of Analysis


No. CA11138

Product: Eichrom's Actinide Resin
Lot #: ACA061311

Eichrom's Actinide Resin Lot # ACA061311 has been manufactured according to Eichrom standard operating procedures. This resin lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.5 M HCl acid	≥ 300	468



Jill Bryant
Quality System Coordinator

June 16, 2011