

Certificate of Analysis

No. CA 15177

Product: Eichrom's Cation Exchange Resin Cartridges
Lot #: C8MR15E

Eichrom's Cation Exchange resin cartridges, Lot # C8MR15E has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

Test Value	Specification	Result
Extractables	< 1 mg(residue) / g(resin)	n.d. mg/g
pH	≥ 4.0	pass
Wet Volume Capacity	min. 1.7 meq/mL	2.0 meq/mL (Per vendor cert)

Jill Bryant
Quality System Coordinator



September 22, 2015