

Certificate of Analysis

No. CA08258


Product: Eichrom's DGA Resin-Branched

Lot #: DBS03255

Eichrom's DGA Resin-Branched Lot # DBS03255 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. The resin meets all Eichrom quality specifications.

ANALYSIS RESULTS

Test Value	Specification	Result
Appearance	White beads	Pass
Eu Activity Balance	95 – 105 %	97 %
Eu Breakthrough	< 2 % in Load	0.4 %
Eu Recovery	> 95 % in 3 BV	95.1 %



Jill Bryant
Quality Control Chemist

December 17, 2008