

**Certificate of Analysis**  
No. CA12110

Product: Eichrom's DGA-Normal Resin  
Lot #: DBS052112

Eichrom's DGA-Normal Resin Lot # DBS052112 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

**ANALYSIS RESULTS**

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
% <sup>241</sup> Am stripped in 0.01M HNO <sub>3</sub>		
Least acceptable mean (4 rep)	90.0	98.0
Max. standard deviation	2.5	0.3
% <sup>241</sup> Am stripped in 0.5M HNO <sub>3</sub>		
Maximum acceptable mean (4 rep)	2.0	0.0
Max. standard deviation	1.0	0.3
% <sup>233</sup> U stripped in 0.5M HNO <sub>3</sub>		
Least acceptable mean (4 rep)	95.0	98.3
Max. standard deviation	2.5	0.1
% <sup>233</sup> U stripped in 0.01M HNO <sub>3</sub>		
Maximum acceptable mean (4 rep)	2.0	0.59
Max. standard deviation	1.0	0.1

  
Sarah Teichma  
Quality Supervisor

6-20-12  
June 20, 2012