



Certificate of Analysis  
No. CA15231

Product: Eichrom's DGA-Branched Resin Cartridges  
Lot #: DBSR15C

Eichrom's DGA-Branched Resin cartridge Lot # DBSR15C has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
% <sup>241</sup> Am stripped in 0.01M HNO <sub>3</sub>		
Least acceptable mean (4 rep)	90.0	98.2
Max. standard deviation	2.5	0.3
% <sup>241</sup> Am stripped in 0.5M HNO <sub>3</sub>		
Maximum acceptable mean (4 rep)	2.0	0.0
Max. standard deviation	1.0	0.3
% <sup>233</sup> U stripped in 0.5M HNO <sub>3</sub>		
Least acceptable mean (4 rep)	95.0	97.8
Max. standard deviation	2.5	0.2
% <sup>233</sup> U stripped in 0.01M HNO <sub>3</sub>		
Maximum acceptable mean (4 rep)	2.0	1.0
Max. standard deviation	1.0	0.1

Jill Bryant  
Quality System Coordinator

November 25, 2015

