


Certificate of Analysis
No. CA08105

Product: Eichrom's DGA Resin-Branched Cartridges
Lot #: DBSR5A

Eichrom's DGA Resin-Branched Cartridges Lot # DBSR5A has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

Test Value	Specification	Result
Appearance	White beads	Pass
Eu Activity Balance	95-105%	97.0 %
Eu Breakthrough	< 2% in Load	0.4 %
Eu Recovery	> 95% in 3 BV	95.1 %



Jill Bryant

June 3, 2008

Jill Bryant
Quality System Coordinator