

Certificate of Analysis


No. CA07042

Product: Eichrom's DGA Resin-Branched Cartridges
Lot #: DBSR7A

Eichrom's DGA Resin-Branched Cartridges Lot # DBSR7A has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

| Test Value | Specification | Result |
|---------------------|---------------|--------|
| Appearance | White beads | Pass |
| Eu Activity Balance | 95-105% | 97.0 % |
| Eu Breakthrough | < 2% in Load | 0.4 % |
| Eu Recovery | > 95% in 3 BV | 95.1 % |



Jill Bryant

Jill Bryant
Quality Control Chemist

Prepared: March 21, 2007