


**Certificate of Analysis**  
No. CA07148

Product: Eichrom's DGA Resin-Normal  
Lot #: DNS09257

Eichrom's DGA Resin-Normal Lot # DNS09257 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

**ANALYSIS RESULTS**

Test Value	Specification	Result
Appearance	White beads	Pass
Eu Breakthrough	< 2 % in Load	< 0.4 %
Eu Recovery	> 95 % in 3 BV	96.9 %



Jill Bryant  
Quality Control Chemist

October 10, 2007