

**Certificate of Analysis**  
No. CA06190

Product: Eichrom's DGA Resin-Normal Cartridges  
Lots: DNSR6A, DNSR6B, DNSR6C


Eichrom's DGA Resin-Normal Cartridges have been manufactured according to Eichrom's standard operating procedures. The lot of resin used to manufacture these lots has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

**ANALYSIS RESULTS**

| <u>Test Value</u> | <u>Specification</u> | <u>Result</u> |
|-------------------|----------------------|---------------|
| Appearance        | White beads          | Pass          |
| Eu Breakthrough   | < 2 % in Load        | 0.4 %         |
| Eu Recovery       | > 95 % in 3 BV       | 96.2 %        |

**cpm/mL vs. Bed Volumes of Eluate**

**50-100  $\mu$ m DGA Resin, normal, Preconditioned with 4.0 M  $\text{HNO}_3$ , 22(1)  $^\circ\text{C}$**



*Jill Bryant*  
Jill Bryant  
Quality Control Chemist

Prepared: January 5, 2007