


Certificate of Analysis
No. CA07071

Product: Eichrom's DGA Resin-Normal cartridges
Lot #: DNSR7A

Eichrom's DGA Resin-Normal cartridge Lot # DNSR7A has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

Test Value	Specification	Result
Appearance	White beads	Pass
Eu Breakthrough	< 2 % in Load	< 0.1 %
Eu Recovery	> 95 % in 3 BV	99.4 %



Jill Bryant
Quality Control Chemist

Prepared: May 25, 2007