

Certificate of Analysis

No. CA08024

Product: Eichrom's DGA Resin-Normal cartridges

Lot #: DNSR8A


Eichrom's DGA Resin-Normal cartridge Lot # DNSR8A has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. The resin used to manufacture this lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

Test Value	Specification	Result
Appearance	White beads	Pass
Eu Breakthrough	< 2 % in Load	< 0.4 %
Eu Recovery	> 95 % in 3 BV	96.9 %

cpm/mL vs. Bed Volumes of Eluate

50-100 μm DGA Resin, normal, Preconditioned with 4.0 M HNO_3 , 22(1) $^\circ\text{C}$



Jill Bryant
Quality Control Chemist

February 11, 2008