

Certificate of Analysis

No. CA08244

Product: Eichrom's DGA Resin-Normal cartridges


Lot #: DNSR8F

Eichrom's DGA Resin-Normal cartridge Lot # DNSR8F has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom's quality assurance policy and procedures and meets all Eichrom quality specifications.

ANALYSIS RESULTS

Test Value	Specification	Result
Appearance	White beads	Pass
Eu Breakthrough	< 2 % in Load	0.0 %
Eu Recovery	> 95 %	98.8 %

DGA Normal Elution Curve



Jill Bryant

December 2, 2008

Jill Bryant
Quality System Coordinator