


Certificate of Analysis
No. CA 17135Product: Eichrom's Ln3 Resin
Lot #: L3S081117

Eichrom's Ln3 Resin Lot # L3S081117 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
%La in load/rinse	95-105	99.9
%Be in load/rinse	95-105	99.9
%U in load/rinse	< 3	1.9
%La in HEDPA strip	< 1	< 0.1
%Be in HEDPA strip	< 1	< 0.1
%U in HEDPA strip	95-105	98.1


Sarah McAlister
Quality Manager
August 17, 2017