



Certificate of Analysis

No. CA13230

Product: Eichrom's LN Resin columns
Lot #: LN13A

Eichrom's LN Resin column Lot # LN13A has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.1 N nitric acid	>700	2178.2


Jill Bryant
Quality System Coordinator

January 10, 2014