

Certificate of Analysis


No. CA09031

Product: Eichrom's LN Resin
Lot #: LNA031209

Eichrom's LN Resin Lot # LNA031209 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.1 N nitric acid	>700	1858



March 27, 2009

Jill Bryant
Quality System Coordinator