

Certificate of Analysis

No. CA12102

Product: Eichrom's LN Resin
Lot #: LNA060512

Eichrom's LN Resin Lot # LNA060512 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.1 N nitric acid	>700	935



Sarah Fejchma
QC Supervisor

6-7-12
June 7, 2012