

Certificate of Analysis


No. CA10128

Product: Eichrom's LN Resin
Lot #: LNA070810

Eichrom's LN Resin Lot # LNA070810 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.1 N nitric acid	>700	1634



July 13, 2010

Jill Bryant
Quality System Coordinator