

Certificate of Analysis

No. CA09155


Product: Eichrom's LN Resin

Lot #: LNA102009

Eichrom's LN Resin Lot # LNA102009 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.1 N nitric acid	>700	1856



Jill Bryant
Quality System Coordinator

October 27, 2009